

Light Sources & Laser Safety

How to Attach and Install

## Laser Shield Curtain Guide

Application Systems

Optics & Optical Coatings

Opto-Mechanics

Bases

Manual Stages

Actuators & Adjusters

Motoeized Stages

Light Sources &

(1) Wipe the glass clean.(2) Spray water onto the glass surface.(3) Apply the product to the wet glass.

(1)

laser light.

(4) Push out water and air from underneath the product by moving a rubber spatula on top of the product from the center to its edges.

depending on the solvent type. If stained, wash with water containing a neutral detergent or wipe with alcohol.

Also wipe with alcohol when curtains become cloudy over time with oily exudations (plasticizer).

These products are intended for protection or shielding from accidental exposure to scattered

• Effective as safety measures for expected and unexpected visitors, since laser injuries can occur instantaneously.

density (OD) of laser light to be absorbed are inscribed on these products.

Apply to windows, inset windows or partitions in laser controlled areas in laboratories and factories.

Absorb indirect scattering light of laser light to protect the eye. The type, wavelength and optical

Replace curtains periodically because the optical density may deteriorate depending on the usage or storage environment (direct sunlight, high temperature and high humidity) or due to scratches.

These products are made of flexible PVC, which may be vulnerable to degradation by organic solvents, acids, and alkalies,

Installation method of YLC-1 laser shield curtain




Guide


Lasers


Detectors


Laser Safety Equipments

Light Sources









## Attention

- > Do not use with incompatible lasers or wavelengths. (Even if laser names are the same, their wavelengths might be different.)
- Do not use products that are damaged or after they have received large laser energy.
- Never subject laser (shield) curtains to direct laser beam exposure. Direct exposure may damage the curtains.
- ▶ These are not protective equipments that completely absorb laser light. (Refer to the absorption characteristic graph.)
- Do not directly look into the laser beam through laser (shield) curtains.

