
FC-114 FC-414 FC-514

# Instruction manual

- FS Series Stage -



Target FV01.032 ~ Target EN01.06





# ■ Notes regarding these materials

- These materials are intended as a reference to assist our customers in the use of the SIGMAKOKI CO., LTD. Product best suited to the customer's application; they do not convey any license under any intellectual property rights, or any other rights, belonging to SIGMAKOKI CO., LTD. or a third party.
- SIGMAKOKI CO., LTD. assumes no responsibility for any damage, or infringement of any third-party's rights, originating in the use of any product data, diagram, charts, programs, or algorithms contained in these materials.
- All information contained in these materials, including product data, diagrams, charts, programs and algorithms represents information on products at the time of publication of these materials, and are subject to change by SIGMAKOKI CO, LTD, without notice due to product improvements or other reasons.
- When using any or all of the information contained in these materials, including product data, diagrams, charts, programs, and algorithm, please be sure to evaluate all information and products. SIGMAKOKI CO., LTD. assumes no responsibility for any damage, liability or other loss resulting from the information contained herein.
- SIGMAKOKI CO., LTD. products are not designed or manufactured for use in equipment or system that
  is used under circumstances in which human life is potentially at stake, SIGMAKOKI CO., LTD. products
  cannot be used for any specific purposes, such as apparatus or systems for transportation, vehicular,
  medical, aerospace, nuclear, or undersea repeater use.
- The prior written approval of SIGMAKOKI CO., LTD. is necessary to reprint or reproduce in whole or in part these materials.
- If these products or technologies are subject to the Japanese export control restrictions, they must be exported under a license Japanese government and cannot be imported into a country other than the approved destination. Any diversion or reexport contrary to the export control laws and regulations of Japan and/ or the country of destination is prohibited.

# Contact

 Tokyo Head Office
 TEL +81-3-5638-6551

 Osaka Branch
 TEL +81-6-6307-4835

 Kyushu Sales Office
 TEL +81-92-481-4300

#### Precautions for use

Please read this instruction manual carefully before use to ensure correct use. Keep this manual in a safe place and use as a reference when operating this equipment.

- This equipment is for controlling the associated feedback stage in closed loop and open loop. Do not use
  it for any other purpose.
- · This instrument is a precision instrument. Please handle with care.
- · Protect from impact or excessive force from being applied to the equipment.
- · Maintain the installation environment described in this instruction manual.



# ■ Safety Please be sure to observe

Explains what you must observe to prevent harm to people and others, and damage to property.

· The degree of harm or damage caused by improper use is classified and explained in the following display.

<b>!</b> DANGER	This symbol indicates that a risk of death or serious injury is imminent.
<b>WARNING</b>	This symbol indicates that there is a possibility of death or serious injury.
<b>A</b> CAUTION	This symbol indicates that there is a possibility of physical damage to the equipment.

 $\boldsymbol{\cdot}$  Important notices are categorized and explained by the following symbols.

0	Indicates required content that must be executed.
0	Indicates prohibited content that must not be executed.
<u> </u>	Indicates a Caution that you should be aware of. Also displayed on the equipment.

· General usage requirements.

	<u></u> <u></u> WARNING
0	Use within the power supply voltage range specified in this instruction manual. Input voltage outside the specified range may cause smoke or fire.
0	Do not use a damaged power cable as this may result in electric shock, short circuit, or fire.
0	Make sure that the power cable plug is fully inserted. Incomplete insertion may cause smoke or fire due to electric shock or overheating.
0	Use the fuse specified in this instruction manual. Using anything other than the rated value may cause smoke or fire.
0	Make sure to ground the protective earth terminal to earth ground to prevent electric shock.
	Do not use in an explosive atmosphere. Never use it in a place where flammable or explosive gas or vapor is present because this may cause an explosion or fire.
$\Diamond$	Do not place in a corrosive atmosphere as this can cause corrosion and poor contact of the electrical connections and can cause malfunction or failure, possibly resulting in a fire.
	Do not use in dusty places. Dust accumulation on the power cable plug or insulation can contribute to damage due to multiple causes, possibly resulting in smoke or fire.
	Do not use the product where it will be exposed to water as this can cause an electric shock or malfunction.
	Do not open the cover. It may cause an electric shock or malfunction.
0	Do not plug or unplug the power cable with wet hands due to a risk of electric shock.
	Do not connect stages other than those described in this instruction manual as this may cause an electric shock or malfunction, risk of a runaway stage, smoke, or fire.
0	Do not use the product if there are any signs of damage or malfunction due to risk of electric shock, smoke, or fire. In case of damage or malfunctions, immediately unplug the power cable from the outlet.



# ■ The symbols on the product



Indicates a protective ground terminal.



Indicates a frame terminal.

# ■ Precautions for moving and shipping

The following are precautions for moving and transporting the equipment. Please observe these precautions,

- · Turn off the MAIN POWER switch, remove all connection cables, and then move the instrument,
- · Do not move the equipment by stacking equipment,
- · When shipping the equipment, please use the original packing materials.
- Be sure to attach this instruction manual when moving or s the product,

#### Precautions for install

The following are precautions for installing the equipment. Please be sure to observe.

- Avoid places where the temperature and humidity are high, where the product is exposed to direct sunlight, or where the temperature changes rapidly.
- Do not install it in a place where it will be exposed to water.
- · Install it on a sturdy and horizontal base.
- There is an intake slit on the right side and an exhaust slit for the cooling fan on the left side. To prevent the internal temperature from rising, provide a space of 50 mm or more on both sides.
- · When using it in a closed space such as a rack, make sure that the ambient temperature of the equipment
- · does not exceed 40 °C.
- · Do not place anything over 20 kg on this equipment.
- · Be sure to ground the protective earth terminal to earth.
- · Connect the stage ground to the frame terminal of the instrument.
- Do not use the product in a place where there is a strong magnetic field or electric field, or where there is a lot of distortion or noise in the input power supply waveform.
- · Allow for space to access the MAIN POWER switch.
- Do not plug the power cable into an outlet that makes it difficult to insert or remove the plug, or place anything in front of the plug that makes it difficult to connect or disconnect.

# Precautions for peripheral equipment connection

The following are precautions for connecting this equipment and peripheral components.

- $\cdot$  Do not connect or disconnect any connectors while the power button is turned on.
- Ony connect feedback stages that conform to specifications our feedback stage equipped with slowdown sensor or origin sensor.
- · Do not connect any active power source to the EMG connector. Passive, mechanical switch only.
- Be sure to connect models that are compatible with the minimum resolution for the feedback stage and feedback stage controller.

# Cleaning

To clean use a soft cloth with a mild detergent diluted with water and wipe gently,

MARNING When cleaning, be sure to turn off the power with the POWER button, turn off the MAIN POWER switch, and unplug the power cable from the outlet.

Use only neutral detergent diluted with water. Discoloration and roughness of the painted surface, fading of printed characters, and cloudiness of the acrylic board may occur.



# Contents

	Notes regarding these materials	<b>-</b> - i
	Contact:	
	Precautions for use	
	Safety Please be sure to observe	
	The symbols on the product	- <b>-</b> iii
	Precautions for moving and shipping	
	Precautions for install	<b>–</b> iii
	Precautions for peripheral equipment connection	<b>– -</b> iii
	Cleaning	- <b>-</b> iii
1.	Overview	<b></b> ·1
2.	Package Contents	<del>-</del> 1
3.	Option:	<del>-</del> 1
	Part names and functions	
	4-1.Front panel·	<b>-</b> 2
	4-3,Right side panel	<b>1</b> 0
	4-4.Left side panel·	<b>-</b> ·10 <b>-</b> ·11
5.	Operations:	<del>.</del> .12
	5-1.List of modes:	<b>-</b> ∙12
	5-2.Common contents	- ·12 - ·13
	5-4,LOCAL Mode:	<b>-</b> ·13
	5-6 TFACH mode:	<b>-</b> ·15
	5-7.MOVEMENT TEST mode:	••16 ••16
	5-9.INFORMATION mode	<b>-</b> ·17
	Commands:	
	6-1. List of commands:	<b>-</b> ∙21
	6-2. Command description format	<b>-</b> ∙25
	6-4. Command description	
	Parameters:	
	7-1. List of parameters ————————————————————————————————————	<b>-</b> 63
	7-3. Procedure for transition to GENERAL parameters———————————————————————————————————	66
	7-5 How to read parameter descriptions:	<b>-</b> .71
	7-6. Parameter description	
	Teaching function	
	8-1. About channels:	<b>-</b> ∙90 - •90



8-3. Teaching operation	92
9.Home return·	
9-1. Mechanical origin return 9-2. Electric home return	97 98
10.About Movement Test·	99
10-1.Movement Test·	99
11.About status·	
11-1.Operation related status	
12.Specification	
12-1.Basic performance	107
13.Dimensions·	118
13-1.Front panel·	<b></b> 118
14. Trouble shooting·	120
15.Update history·	·122
16,Index:	123



1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
	_	S	4	$\cup$	O	- 1	0	9	10	1 1	1 /	10	14	10	10

# 1. Overview

This equipment is for controlling the feedback stage in closed loop and open loop. The communication interface is GP-IB, USB, or Ethernet, and remote operation is possible using ASCII communication commands. Using the included teaching function, this equipment can automatically execute programs written to the internal memory. The teaching function can be operated from the front panel, general-purpose I / Q, all communication interfaces, and the jog controller.

# 2. Package Contents

Listed below are the items that are included in the shipping container. If anything is missing or damaged, contact us or your distributor.

Name	Model or specification	Quantity	Remarks
Feedback stage controller	FC-114 FC-414 FC-514	1	Any one
Power cable	Plug: KP-218 Connector: KS-16A Code: VCTF3 × 0.75mm2	1	2.3m, Tracking resistance
Protective ground wire, frame connection wire	AWG18	2	3m, Green/Yellow
Fuse	250V, 2.5A, Time lag	4	2 pre-installed at the time of shipment
Instruction manual	-	1	This book

#### · Power cable

For power supplies other than AC100 V, use an approved power supply cable which meets the safety standards where the controller will be used. If you have any questions, please contact the distributor. For the specifications of the connectable power cable and the procedure for inserting the plug of the power cable into the power outlet, refer to the following

Connector	Code	Plug	Length
IEC-60320-C13	Type SJT, No16 AWG Min.	NEMA6-15P	2.3m or less
Rated: 7A, 250V	3-Conductors	Tandem blade	
UL, CSA Approved	(Single phase: 2-current carrying & ground)	Rated: 7A, 250V	
	UL, CSA Approved	UL, CSA Approved	

Compatible commercial power supply
AC100V ~ 240V, 50/60Hz

# 3. Option

See below for connectable options. Change the parameter "Option type" according to the option to connect. For details, see our catalog or website.

### · Connectable options

Model	Parameter "Option type" setting values
JC-01, JC-01-3, JC-01-4	ТуреО
MD-400	Type1


⚠WARNING Do not connect other than the above options.



# 4. Part names and functions

#### 4-1.Front panel

Description of the names and functions on the front panel.



# (1) Display unit



The display contents are menus such as coordinate values, status, and F1-F8 button operations. For details, See "5. Display and operation".

 $(2)F1/ \Rightarrow \text{ button}$   $(3)F2/ \Rightarrow \text{ button}$   $(4)F3/ \Rightarrow \text{ button}$   $(5)F4/ \Rightarrow \text{ button}$   $(5)F4/ \Rightarrow \text{ button}$ 

F1 to F4 execute the functions displayed on the right side of the display unit. The arrows are used in parameter mode. For F1 to F4, See "5. Display and operation" for more information. For arrows, See "5. Display and operation" and "7. Parameters".

#### (6)F5 / Enter / Page button



F5 performs the function displayed on the right side of the display. Enter is used in parameter mode, and Page is used in INFORMATION mode. For F5, See "5. Display and operation" and "7. Parameters".

# (7)F6/Cancel / Return button



F6 performs the function displayed on the right side of the display. Cancel is used in parameter mode, and Return is used in INFORMATION mode. For F6 Return, See "5. Display and operation" and "7. Parameters".

(8)F7 button



(9)F8 button

F7 and F8 execute the functions displayed on the right side of the display unit. General is used to display parameter mode, and Setup is used to display setup mode. For F7, F8 and Setup, See "5. Display and operation". For General, See "7. About parameters".



# (10) Restart button

Restart

"Restart" will restart the controller. When the button is pressed and released, the startup operation is the same as when the power is turned on. The state after restart is shown below.

ltem	Contents					
Parameters	arameters will not be reset to default values, but if they are being changed, ne changed content will not be reflected. See "7. About parameters".					
Mode	After Restart, the mode will be the mode set by the parameter "Ini Mode".					
Coordinate value	Both axes will be set to zero.					
Stage	If a stage is moving, it stops. If it is stopped, it remains stopped. At startup, the motor coils will be powered off momentarily, then powered back up.					
Command	Any pending command will be discarded.  If the command "A:" or "M:" is executed without setting the command "D:" or "ACC:" after the restart, the speed will be set to parameters "Max Speed" and "Acc Time". See "7. About parameters".					
Teaching registration contents	It will not be deleted, but if the command "T_OFF:" is not executed after editing, the edited content will not be reflected and will be discarded.					
Communication	Communication will be disconnected. After the restart is complete, it will be in the initial state, so it will be the mode set by the parameter "Ini mode".  Depending on the control software used, communication may resume after the equipment has been restarted or after the mode has been changed, resulting in unintended stage operation.					
Link during USB communication	Will be disconnected. To relink, the software must close the COM port and then reconnect the cable. For USB, See "(2) USB".					
Jog controller	Reboot according to the mode.					

(11) POWER button



(12) On lamp

# (13)Off lamp

Operates the equipment power. Press to turn the power on when the MAIN POWER switch on the rear panel is ON. Press to turn off the power when the power is on or in sleep mode. However, it cannot be used when the "On lamp", "Off lamp", and "Sleep lamp" are off. The operation method is shown below.

Power on	Powe	er off	Usage prohibited
Sleep	Sleep Sleep		Sleep
On Off	On On Off		On U
<b>↓</b>	•	Ļ	
Press for over 0,2 seconds	Press for over	er 1 seconds	
Sleep		Sleep	
On ● Off	On Off	Ф	

If the power does not turn on even if you press it for more than 0.2 seconds when turning on <a href="https://www.wishen.com/www.nc-numbers.com/ww

• Be sure to read "Installation" and "Peripheral equipment connection" before turning on the power.



 If you do not intend to use the unit for a long time, turn off the MAIN POWER switch on the rear panel.

Before turning off the MAIN POWER switch on the rear panel, turn off the power with this
key. If the MAIN POWER switch on the rear panel is turned OFF without turning off the power
with this button, the state of this button will remain ON. When the MAIN POWER switch is
turned on, the power to the instrument turns on without operating this button.



(14) LOCAL lamp

(15) USB lamp

(16) Ethernet lamp

(17) GP-IB lamp

The green indicators light up in the currently selected mode. In TEACH mode, both the "LOCAL lamp" and "REMOTE lamp" are lit.

LOCAL		REMOTE		TEACH			
LOCAL	USB		GP-IB	USB	Ethernet	GP-IB	
Mode Local ●	Mode Local	Mode Local	Mode Local	Mode Local	Mode Local ●	Mode Local •	
USB Ethernet GP-IB	USB  Ethernet GP-IB	USB Ethernet ● GP-IB	USB Ethernet GP-IB	USB  Ethernet GP-IB	USB Ethernet ● GP-IB	USB Ethernet  GP-IB	
- LSTN	- LSTN	- LSTN	- LSTN	- LSTN	- LSTN	- LSTN	
TALK SRQ	SRQ SRQ	SRQ	SRQ SRQ	SRQ	SRQ	SRQ SRQ	

(18) LSTN lamp

(19) TALK lamp

(20) SRQ lamp

When communication is performed and the communication interface is set to GP-IB, the LSTN lamp, TALK lamp, and SRQ lamp will be lit depending on the status. The LSTN lamp lights when this equipment receives a command, and the TALK lamp Lights when sent to. The SRQ lamp lights when a service request is executed from this equipment.

LSTN	TALK	SRQ		
Mode	Mode	Mode		
Local	Local	Local		
USB	USB	USB		
Ethernet	Ethernet	Ethernet		
GP-IB ●	GP-IB ●	GP-IB ●		
LSTN • TALK SRQ	LSTN TALK SRQ	LSTN TALK SRQ		

# (21) Motor 1, 2 lamp

This shows the motor excitation status of the stage. The axis whose green lamp is lit is energized and the axis whose light is off is demagnetized.

First axis excitation On Second excitation On	First axis excitation On Second excitation Off	First axis excitation Off Second excitation On	First axis excitation Off Second excitation Off
Motor	Motor	Motor	Motor
1 •	1 •	1	1
2 •	2	2 •	2

#### (22) Close 1, 2 lamp

Indicates the stage control loop status. An axis whose green lamp is lit is in a closed loop state, and an axis that is not lit is in an open loop state.

First axis Close	First axis Close	First axis Open	First axis Open		
Second axis Close	Second axis Open	Second axis Close	Second axis Open		
Close 1 2	Close	Close	Close		
	●1	1	1		
	2	• 2	2		



1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
	_			_	_			_			. —	. —			

# (23) SLEEP lamp

Sleep

Lights when the equipment is in sleep mode. At this time, the On lamp goes off.

### (24) Sleep function

In sleep mode, the display turns off and the motor is demagnetized. Modes, coordinate values, statuses, and commands retain the state prior to sleep execution. When returning from the sleep state, if the stage table has moved to a different position from before sleep execution, it will move to the coordinate value before sleep execution by feedback control. However, if an error other than a limit error, overflow error, or emergency stop occurs in the sleep state, the coordinates before the sleep execution cannot be restored. In addition, when the instrument is restarted due to a power failure or instantaneous voltage drop, all retained modes, coordinate values, statuses, and commands are discarded. For error, See "11. Status". For Restart, See "(10) Restart button". Becomes a little warm in sleep mode, but this is not a malfunction.

# (25) JC connector



A connector for connecting options. Change the parameter "Option type" according to the option to connect.

#### Connectable options

涇

	Model	
JC-01, JC-01-3, JC-01-4		
MD-400		


MARNING Do not connect other than the above options.

When removing the jog controller from the instrument, shut off the power supply to the jog controller.




#### 4-2.Rear panel

Names and functions on the rear panel,



#### (1) AC connector



Connector for power cable. When the power supply is AC100 V, use the included power cable. For power supplies other than AC100 V, use an approved power supply cable which meets the safety standards where the controller will be used. If you have any questions, please contact the distributor. For the specifications of the connectable power cable and the procedure for inserting the plug of the power cable into the power outlet, refer to the following

Connector	Code	Plug	Length
IEC-60320-C13	Type SJT, No16 AWG Min.	NEMA6-15P	2.3m or less
Rated: 7A, 250V	3-Conductors	Tandem blade	
UL, CSA Approved	(Single phase: 2-current carrying & ground)	Rated: 7A, 250V	
	UL, CSA Approved	UL, CSA Approved	

Step1 Make sure the power supply is within the AC 100 V to 240 V range.

Step2 Check that the MAIN POWER switch is OFF.
Step3 Connect the power cable to the AC inlet.
Step4 Insert the power cable plug into an outlet.


 $ilde{\mathbb{M}}$ WARNING Connect to an overvoltage category II power supply.



- Use only the power cable that comes with the product or that meets local safety standards.
- Do not use the power cable supplied with this equipment as the power cable for other equipment.
- · Please follow the procedure for plugging the power cable into the outlet.



# (2) Fuse holder




The fuse holder is under the AC connector. Refer to the following for fuse specifications and replacement procedure.

Voltage	Current	Characteristic	Size	
AC250V	2.5A	Time lag	$\phi$ 5mm $ imes$ 20mm	

Step1 Check that the MAIN POWER switch is OFF.
Step2 Unplug the power cable from the outlet.
Step3 Disconnect the power cable from the AC inlet,

Step4 Hook the tip of a thin flathead screwdriver or similar tool into the groove in part A below, and apply force in the direction of the arrow to remove the holder.



Step5 Replace both fuses with the above specifications.

Step6 Push the holder with the replaced fuse firmly into the AC inlet.

Step7 Connect the power cable to the AC inlet.

Step8 Connect the power cable according to the AC inlet insertion procedure.

- · Do not use the included fuse for other equipment.
- When replacing, always replace both fuses

#### (3) Protective earth terminal



Electrical ground terminal. Connect the round crimp terminal of the included protective conductor wire to this terminal. When the power cable includes a ground connection, this terminal does not need to be used.

When using a power supply cable without a ground terminal, a separate ground connection MARNING must be used due to this being a class I apparatus. There is a risk of electric shock if not grounded.



- · Use the protective grounding cable supplied with the instrument for grounding.
- Do not use the included the protective grounding cable for other equipment.





# (4) Enclosure terminal





Connect this terminal to a surface plate or metal mount that holds the stage to provide a common ground between the controller and the stage. Connect the round crimp terminal side to this terminal.

MARNING Do not use as a protective earth terminal.

MCAUTION If not connected, the system may not function correctly.

 $\triangle$ 

- Use the attached frame connection line for connection.
- $\cdot$  Do not use the frame connection cable supplied with this equipment for other equipment,

# (5) MAIN POWER switch





OFF

Turns the main power on and off. This switch alone does not turn on the power to the instrument. The POWER button on the front panel must also be pressed to turn on the instrument.


- Be sure to read "Installation" and "Connection" before turning on the power with this switch.
- · If not using the product for a extended time, set this switch to OFF.



Turn off the power with the POWER button on the front panel before turning off this switch.
If this switch is turned off without turning off the power with the POWER button, the state of
POWER button will remain ON. When the MAIN POWER switch is turned on, the power to the
instrument turns on without operating this button.

#### (6) Motor cable connectors

1 MOTOR 2




Connections for the stage motor cables. Connect the axis one stage to the left connector and the axis two stage to the right connector.

A CALITION

Be sure to correctly connect the motor cable and scale cable to the first and second axes with the power off. If the connections are incorrect, unintended operation will occur.



# (7) Scale cable connectors



Connections for the stage scale cables. Connect the axis one stage to the left connector and the axis two stage to the right connector..

 $ilde{\Lambda}$ WARNING Please connect only approved feedback stages that meet our stage specifications.

ACAUTION Be sure to correctly connect the motor cable and scale cable to the first and second axes with the power off. If the connections are incorrect, unintended operation will occur.

# (8) GP-IB connector



Connection for the GP-IB cable, Use this when you want to operate the instrument via the GP-IB interface, For details, see  $^{''}(1)$  GP-IB $^{''}$ .

# (9) USB connector (B type)



Connection for the USB cable. Use this when you want to operate the instrument via the USB interface. For details, see  $^{''}(2)$  USB $^{''}$ .

#### (10) Ethernet connector (LAN)



Connection for the Ethernet cable. Use this when you want to operate the instrument via the Ethernet interface, For details, see "(3) Ethernet".

(11)I/O connector



The controller has a general-purpose input and output connector which includes digital input and output, teaching operation input and status output, scale division signal output, in-position signal output, and alarm signal output. Refer to "(4) General purpose I/O" for details.

# (12) EMG connector



Connection for an emergency stop switch. Use this when you want to stop the stage operation in an emergency, Refer to "(5) Emergency stop" for details,

MARNING To enable the emergency stop switch, change the parameter "EMG Connecter" to Enable and perform a test operation to confirm that it active.


↑ CAUTION Do not connect any device that supplies power. Passive switch only.



-1	2	2	4	5	6	7	0	0	10	11	10	12	11	15	16
l l	_	0	4		O	- 1	0	9		1 1			14		

### 4-3. Right side panel

Names and functions on the right side panel.



# (1) Ventilation slots


# (2) Dial switch for motor driver setting

Dial switch for motor driver setting, For details, refer to "(6) Motor driver setting switch".

CAUTION Depending on your environment, you may be asked to change the settings. Do not change any other settings. If it is changed incorrectly, unintended stage operation may occur.

# 4-4.Left side panel

Names and functions on the left side panel.



# (1) Ventilation slots

CAUTION Do not block the intake slit. See "Installation" for details.

# (2) Cooling fan



1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16

# 4-5. Audio Indicator

An audible beep indicates certain operations. The beep can be turned off with the parameter. Note that the beep may be interrupted at times; this is not a malfunction.

Length of beep	Operation
About 0.1 seconds	Button operation
About 0.1 seconds	Completion of command "RESET:"
About 0.1 seconds twice	Start of command "RESET:"
About 0.3 seconds (two-step sound)	Transition to sleep mode and return
About 0.4 seconds	POWER button operation
About 0.4 secorius	Restart button operation
About 1 seconds	Entering stage stroke limits
About 2 seconds	Error occurrence
Continuous tone	Internal power supply output short circuit or overload The volume depends on the degree of overload of the internal power supply. It also sounds in sleep mode. See "14. Troubleshooting" for details.



# 5. Operations


# 5-1.List of modes

The operational modes are as follows.

Mode	Contents
SETUP	Switch modes.
LOCAL	Local operation using the front panel buttons and jog controller.
REMOTE	Remote operation using a communication interface.
TEACH	Automatic execution with commands registered in the internal memory is possible.
MOVEMENT TEST	Execute a test operation to verify stage function.
PARAMETER	Edit operational Parameters.
INFORMATION	Displays controller information.

# 5-2.Common contents

The display contents are common to all modes.




No	ltem	Contents												
1	Axis name	The upper disp The display cor												two.
		<b>)</b>						Motion completed successfully, Position is Stable. *1						
2	Desitioning status	<b>&gt;</b>	READY All operations are accepted.			Motio *1	n di	d no	ot cc	mpl	ete :	succ	ess	iully.
_	Positioning status	>	BUSY Operation is denied.			Motic comp				s an	d ha	as n	ot b	een
		no display	BUSY	Ope deni	eration is ied.	Comi occur			pera	atio	n or	eri	or	has
		The upper display value containe. See the property that is a second to the property of the pro	an be either	the	read value									
13	Coordinate value	Model	Minim	Example (unit :mm)										
	TO SO TO SIN TOLON	FC-114	100nm					0		0	0	0	1	m
		FC-414	FC-414 50nm				0		0	0	0	0	5	m
		FC-514	10nm				0		0	0	0	0	1	m
		The upper disp displayed is set are set to millim	by the para											
		Unit	Cor	Example (Model: FC-114)										
4	Unit	n	Nanom	eter	(nm)	1	2	3	4	5	6	0	0	n
		u	Microm	eter	(um)		1	2	3	4	5		6	u
		m	Millime	ter (ı	mm)		1	2		3	4	5	6	m
	ii.	٥	Degre				_1_	2		3	4	5	6	•
			Minimum digit is					1	2	3	4	5	6	
5	Status	See "11. Status" for more information.												
	Jog controller	J	"J" in the blo	ck	Power is su	pplying	to t	he Jo	og cc	ontro	ller.			
6	connection		Block only		The power									off.
	confirmation mark		o display		The Jog co	ntroller	con	nect	or is	not d	conn	ecte	d.	
7	Menu	This is a function	on assigned f	rom	his is a function assigned from F1 to F7.									

<sup>\*1</sup> See parameter "INPos Range" for in-position range. \*2 BEC button and command "BEC:" can be used. \*3 BEC button and command "BEC:" cannot be used.



#### 5-3,SETUP Mode

SETUP mode is used to choose the item of each mode, sleep and the Jog Controller's power supply. To switch to this mode, press the F7 and F8 buttons simultaneously when the stage is stopped in LOCAL, REMOTE, or TEACH mode.



No	Item	Contents
1	Mode	Block only
2	Menu	This is a function assigned from F1 to F7.


Button	Display	Contents							
F1	MOD	Select a mode, Each time you press it, it switches,							
ГІ		LOCAL $\rightarrow$ REMOTE $\rightarrow$ TEACH $\rightarrow$ Back to top							
F2	PRM	Move to PARAMETER mode.							
F3	SLP	Transit to sleep state. For details, see "(24) Sleep function". To return from sleep mode, press any of the F1 to F8 buttons,							
F4	INF	Move to INFORMATION mode,							
F5	JON * 1	Supply power to the Jog controller.							
F6	JOF * 1, 2	Power supply to the jog controller is stopped.							
F7	RTN	Exit the SETUP mode,							
F8	-	None							

<sup>\* 1</sup> If the jog controller connector is not connected, it will not be displayed. If it is connected, either will be displayed.d.

st 2 Can be used only when positioning is completed or an error occurs and the stage is stopped.

#### 5-4.LOCAL Mode

In LOCAL mode, the stages can be operated from the front panel of the instrument. In this mode, menus 1 to 3 are used to operate the stages from the front panel. Use the F7 or F8 button to switch menus. For the procedure to switch to this mode, refer to "5-3. SETUP mode".



No	ltem	Contents
11		Represents the quantity of movement pulses output to the stage with a single click of the CW or CCW button. Use "PLS" to select the pulse amount. Only menu 1 is displayed.
2	Mode	L (LOCAL) is displayed.
3	Menu	Displays the functions assigned to F1 through F7.



Menu	Button	Target axis	Display	Contents
	F1	1	CCW	Moves the stage in the CCW direction (Motor side). Press and hold for continuous operation, and press once (single click) for pulse operation. See LS" for pulse amount selection. Flashes when at the limit of travel.
	F2	2		The control status during continuous operation differs depending on the setting of the parameters "Jog Cont" and "Stage Cont Type" . *1
	F3	1	3, 2, 1	Sets the speed of the stage when operated by the CW or CCW button. The operation speed changes each time the button is pressed. The speed depends on the setting of parameters "Jog Speed 3", "Jog Speed 2" and og Speed 1". 1 at startup. Axis display depends on the parameter "AXIS
	F4	2		Name"setting. In this case, the first axis side is X and the second axis side is Y.  [Speed stage1 => Speed stage2 => Speed stage3 => Back to Speed stage1]
				Moves the stage in the CW direction (Anti motor side). Press and hold for
1	F5	1	CW	continuous operation, and press once (single click) for pulse operation. See LS" for pulse amount selection. Flashes when at the limit of travel.
	F6	2		The control status during continuous operation differs depending on the setting of the parameters "Jog Cont" and "Stage Cont Type" . *1
	F7	-	PLS, PST	Sets the number of movement pulses output to the stage when the CW or CCW button is clicked once. When the display is PLS, the value changes each time it is pressed. It is "1" at startup. The speed depends on the parameter "Jog Speed 1". When the stage is moving, the display changes to blinking PST, and the function of this button switches to stop movement. When movement is completed or stopped by this button, the display returns to PLS. [1pulse => 10pulse => 100pulse => 100pulse => Back to top]
	F8	_	MEU	Switch to menu 2.
	F7+F8	-	-	If both button are pressed at the same time switches to the SETUP mode. For details on the SETUP mode, see "5-3. SETUP mode".
	F1	1	ZEO	Set the coordinate value of the target axis to zero. *1
	F2	2		
	F3	_	MTR	Set the motor power to on or off. It changes each time it is pressed. For state, check the Motor lamp on the front panel.  [Motor coils excited - Power on : Exc]  Moter 1 Exc => Exc => Not => Not => Back to top  Moter 2 Exc Not => Not => Not
2	F4	-	CLS	Sets the positioning control method (Closed loop or Open loop). It changes each time it is pressed. For state, check the Close lamp on the front panel.  Close 1 Close => Close => Open => Open => Back to top  Close 2 Close Open
	F5	1	DEC	Clear Error (when positioning status display is $\tilde{S}$ ) or Cancel ESTOP (when emergency stop is executed). When limit error, overflow error or TEACHING
	F6	2	BEC	command error, errors can be canceled without affecting the coordinate values. To cancel the emergency stop, press either button.
	F7	-	RTN	Return to Menu 1.
	F8	-	MEU	Switch to menu 3.
	F7+F8	_	-	If both button are pressed at the same time switches to the SETUP mode. For details on the SETUP mode, see "5-3, SETUP mode".
	F1	1	ORG	Return the stage to the mechanical origin . Even if the software limit is set, the software limit is ignored. The return method can be selected by the parameter "ORG Mode Sel" . Refer to "9. Return to origin" for the return
	F2	2		method. For details on software limits, see "3. Options", "+ Soft LMT Pos" and "-Soft LMT Pos". *1
3	F3 F4	1, 2	STP	Stops movement of machine origin return and the electric origin return.
	F5 F6	1 2	EOG	Return the stage to electrical origin, Refer to "9. Return to origin" for the return method, $*1$
	F7	-	RTN	Return to Menu 1.
	F8	-	MEU	Return to Menu 1.
	F7+F8	_	-	If both button are pressed at the same time switches to the SETUP mode, For details on the SETUP mode, see "5-3. SETUP mode".

<sup>\* 1</sup> While one axis is operating, the other axis can be operated.

Before demagnetizing the motor or operating the stage, make sure that there is no influence on the surroundings.


ı															
	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15

# 5-5.REMOTE mode

Remote mode allows operation of the stages by remote control from an external computer. For the procedure to switch to this mode, refer to " ① SETUP mode". For parameters, refer to "7. Parameters". For each communication interface, refer to "12. Specifications".

 $\triangle$ 

Connect only the communication cable appropriate to the communication interface set in the parameters.



No	Item	Contents
1	Mode	R (REMOTE) is displayed. When shifting to the parameter setting mode, P (PARAMETER) is displayed. For details, refer to "6. Commands".

Button	Display	Contents
F1-F8	-	None
F7+F8		If both button are pressed at the same time switches to the SETUP mode. For details on the SETUP mode, see "5-3, SETUP mode".

# 5-6.TEACH mode

TEACH mode allows the operation program registered in the internal memory of the equipment to be executed. For the procedure to switch to this mode, refer to " ① SETUP mode". For the operation, refer to "8. TEACHING function".




1	2	3	4	5	6	7	8	9	10	11	12	13	14	15

# 5-7,MOVEMENT TEST mode

Movement Test mode moves the stage with one of three prepared movement patterns. This mode can be switched by double clicking the F5 button after entering the TEACH mode.

Restart is necessary to exit this mode, For TEACH mode, refer to "5-6. TEACH mode". For details of operation, refer to "10. Movement test".



No	Item	Contents
1	Menu	It is the function of the F1 to F7 buttons.

Button	Target axis	Display	Contents
F1	1, 2	SMT	Perform step movement, Pitch is ten percent of the stage's stroke,
F2	1, 2	LMT	Perform reciprocating motion between limit sensors.
F3	1	3, 2, 1	Select the operation speed of the test operation. The operation speed
F4	2		changes each time the button is pressed. The speed depends on the setting of parameters "Jog Speed 3", "Jog Speed 2" and "Jog Speed 1". 1 at startup.
F5	1, 2	MMT	Perform step movement, Pitch is 1mm.
F6	1, 2	STP	Stops test operation.
F7	-	RST	Reset the equipment to exit this mode.
F8	-	-	None

Before operating the stage, make sure that there is no influence on the surroundings.

# 5-8.PARAMETER mode

PARAMETER mode is used to check or change the parameters. Refer to "5-3. SETUP mode" for how to enter this mode, For the operation, see "7. Parameters".




# 5-9.INFORMATION mode

INFORMATION mode displays information regarding this equipment on the display panel. For the procedure to switch to this mode, refer to  $^{\prime\prime}5$ -3. SETUP mode $^{\prime\prime}$ .

# (1) Equipment information

Displays information such as the controller model name.



No	ltem		Contents					
1	Page	Page nu	age number					
2	Category	Model	Indicates Equipment information.					
		Vendor	Vendor name					
		Model		FC-114				
			Model name	FC-414				
3	Information			FC-514				
		S/N	Serial Number					
		E/N	Equipment Number					
		F/V	Firmware Version					


Button	Display	Contents
F1-F4	-	None
F5 (Page)	-	Switch to next page.
F6 (Return)	-	Returns to the mode before the transition to INFORMATION mode.
F7-F8	-	None



1	2	3	4	5	6	7	8	9	10	11	12	13	14	15

# (2) GP-IB interface information


If the parameter "I / F Sel" is set to GP-IB, the GP-IB settings are displayed.



No	Item		Contents						
1	Page	Page nu	ge number						
2	Category	IF	It means the communication information						
		IF	The configured communication interface.						
		Addr	GP-IB address setting value						
3	Information	Del	GP-IB delimiter						
		EOI	GP-IB end of Identify						
		SRQ	GP-IB service request						

Button	Display	Contents
F1-F4	-	None
F5 (Page)	-	Return to first page.
F6 (Return)	-	Returns to the mode before the transition to INFORMATION mode.
F7-F8	-	None

(3) USB interface information If the parameter "I / F Sel" is set to USB, the USB settings are displayed.




No	ltem		Contents					
1	Page	Page nu	ge number					
2	Category	Ŀ	It means the communication information					
3	Information	E	The configured communication interface.					
٥	Information	Del	USB delimiter					

Button	Display	Contents
F1-F4	-	None
F5 (Page)	-	Return to first page.
F6 (Return)	-	Returns to the mode before the transition to INFORMATION mode.
F7-F8	-	None

1	2	3	4	5	6	7	8	9	10	11	12	13	14	15

### (4) Ethernet interface information

If the parameter "I / F Sel" is set to Ethernet, the Ethernet settings are displayed.



No	Item		Contents						
1	Page	Page nu	age number						
2	Category	IF	It means the communication information						
		F	The configured communication interface.						
		Del	Ethernet delimiter						
		IP	Ethernet IP address						
3	Information	DGW	Ethernet default gateway						
		SNM	Ethernet subnet mask						
		MAC	Ethernet MAC address						
		PORT	Ethernet Port Number (Fixed)						

Button	Display	Contents
F1-F4	-	None
F5 (Page)	_	Moves to the next page for the second page and to the first page for the third page.
F6 (Return)	-	Returns to the mode before the transition to INFORMATION mode.
F7-F8	-	None



1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16

# 6. Commands

Commands are sent and received from the communication interface in ASCII code format,

 $\triangle$ 

- Before control, check that the communication interface is operating normally.
- Please familiarize yourself with the functions before using the commands.
- · Set the command transmission interval to 10msec or more.
- · Please check clearance around of the stages before operation.

# 6-1. List of commands

The modes that allow these commands to be used are LOCAL, REMOTE, and TEACH modes. Commands cannot be used in other modes.

# (1) Available modes

Describes the modes in which the command can be used.

Symbol	Contents					
L	It indicates that it can be used in the LOCAL mode.					
R It indicates that it can be used in the REMOTE mode.						
R	Indicates that it can be used after switching from REMOTE mode to PARAMETER mode.					
T	It indicates that it can be used in the TEACH mode.					
Т	Indicates that it can only be registered in the internal memory as a teaching program,					

# (2) Available positioning states

Describes the positioning states in which the command can be used.

Symbol	Contents					
R It indicates that it can be used in the positioning completed state.						
B Indicates that it can be used during stage move.						
_	This command is not related to the positioning status.					

# (3)List

Item	No,	Command name	Description	٨	/lod	е	State	Page
Control	01	RESET:	Reset controller	L	R	Т	R B	26
	02	*RST	Reset controller	L	R	Т	R B	26
	03	RESET_RET:	Reset controller ("END" is returned after completion)	L	R	Т	R B	26
	04	SLEEP:	Turn on or off the sleep function	L	R	Т	R	26
	05	SLEEP_RET:	Turn on or off the sleep function (Status is returned after completion)	L	R	Т	R	26
	06	SLEEP?	Get the sleep status	L	R	Т	R	27
	07	MODE:	Change mode	L	R	Т	R	27
	80	MODE?	Get current mode	L	R	Т	R B	27
	09	F:	Change the positioning control method (Closed or Open loop)	R		Т	R	28
	10	FR:	Get positioning control method (Close or Open loop)		R		R	28
	11	C:	Change the excitation state of the motor	R		Т	R	28
		CR:	Get the excitation state of the motor		R		R	29
	13	BEC:	Cancel busy and error conditions	L	R	Т	R B	29
	14	CMDR:	Get the last sent command (Exclude this command)	L	R	Т	R B	29
	15	ECHO:	Change command echo back state		R		R	30
		ECHOR:	Get command echo back status		R		R	30
Information	17	*IDN?	Get the equipment information	L	R	Т	R	30
		VENDOR:	Get the vendor name		R		R	30
		MODEL:	Get the model name		R		R	31
		SN:	Get the serial number		R		R	31
		EN:	Get the equipment number		R		R	31
		FV:	Get the firmware varsion		R		R	31
		RESO:	Get the resolution		R		R	32
		LIMR:	Get the stage stroke		R		R	32
		AN:	Get the axis name		R		R	33
	26	UNT:	Get the unit		R		R	33



Item	No,	Command name	Description	Mod	de	State	Page
Motion status	27	Q:	Get the status (Position coordinate values, Error, Motion, and Positioning)	R	Т	R B	34
Information	28	SRQ:	Get the status (Error, Motion, and Positioning)	R	Т	R B	37
ii ii oi i ii atioi i		P:	Get the position coordinate values	R	÷	R B	39
		ER:	Get the error status	R	Ť	R B	40
	31	STS:	Get the motion status	R	÷	R B	42
	32	!:	Get the positioning status	R	Ť	R B	43
Origin		H:	Return to mechanical origin	R	T	R	44
	34	Z:	Return to electrical origin	R	Ť	R	44
		R:	Set the position coordinate values to zero	R	Ť	R	45
		LIMG:	Perform the stage stroke detection movement	R		R	45
Motion		L:	Operation stop and emergency stop	R	_	R B	45
		ACC:	Set the acceleration and deceleration time	R	Т	R	46
	39	ACCR:	Get the acceleration and deceleration time	R		R B	46
	40	D:	Set the movement speed of the stage	R	Т	R	47
	41	DR:	Get the movement speed	R		R B	48
	42	A:	Set the position to move in absolute motion	R	Т	R	49
		M:	Set the amount of movement to move in relative motion	R	Ė	R	50
		G	Start moving (The setting values is lost after execution)	R	Ť	R	51
	45	GN:	Start moving (The setting values is not lost,)	R	†	R	51
	46	GC:	Delete the setting values of the "A:" and "M:" commands	R		R	52
	47	GR:	Get the setting values of the "A:" and "M:" commands	R		R B	52
	48	JG:	It moves according to the set number of pulses	R	Т	R	53
	49	JY:	Start moving without specifying a target point	R	' '	R B	53
Position		PIT_DEL:	Delete location information registered for the specified number	R		R	54
registration		PIT_SET:	<del>-</del>	R		R	54
1061311 411011		PIT_GET:	Register current coordinate value to specified number	R			54
	52	PII_GET	Get the coordinate value registered in the specified number	n		R	54
	53	PITG:	Start Moving to the coordinate value registered in the specified number	R		R	54
General	54	<b> </b> :	Get general-purpose input status	R	Т	R B	55
purpose I/O	55	O:	Set general-purpose Output status	R	Т	R B	55
Teaching	56	T_ON:	Enter edit mode	Т		R	56
	57	T_OFF:	Save contents and exit	Т		R	56
	58	T_DEL:	Delete content	Т		R	56
	59	T_SET:	Set content	Т		R	56
	60	T_GET:	Get content	Т		R	56
	61	TC:	Select a channel	Т		R	56
	62	TCR:	Get current channel	Т		R B	57
	63	TQ:	Get teaching status	Т		R B	57
	64	TG:	Teaching start	Т		R	57
	65	TP:	Pause	Т		R B	57
	66	TO:	Execute line by line	Т		R	57
	67	TL:	Stop execution	Т		R B	58
	68	TR:	Get channel subscription status	Т		R B	58
	69	TFR:	Get the loop count	Т		R B	58
	70	TM:	Set the Teaching monitor function	Т		R	59
	71	TMR:	Get the Teaching monitor function setting	Т		R	60
	72	TNR:	Get current line number	Т		R B	60
		TACR:	Get current execution command	Т		R B	60
Teaching		FS:	Loop setting	T		_	61
registration	75	FE:	Set loop end	Т		-	61
only	76	END:	Set the execution end line	Т		-	61
	77	T:	Set the wait time	Т		-	61
	78	GIS:	Wait until the specified general-purpose input state is	Т		-	61
Paramatar		PRM_ON:	reached Switch to parameter setting mode	R		R	62
Parameter settings		PRIVI_ON: PRM_OFF:	Switch to parameter setting mode				
Soun igo		PRIVI_OFF: PRM_SET:	Exit parameter setting mode	R R		R B	62
		PRIVI_SET:	Set parameters Get parameters			R B	62
	02	I DIVI_GET	Get parameters	R	l	LU R	62




1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16

### 6-2. Command description format

Explains how to read the command description.

#### (1) Send command format

Describes the send command format. Sending means sending a command from the computer to the equipment. There is no reply from the equipment to the computer after receiving these commands. The command format is not case sensitive.



No,	Item	Subcommand	Description
		LOCAL	Switch to LOCAL mode
1	Control	REMOTE	Switch to REMOTE mode
		TEACH	Switch to TEACH mode
<u> </u>			

√6

No	Ite	em	Contents	Description
1	Command N	lumber	-	The command number.
2	Command N	lame	-	The name of the command.
3	Feature Des	cription	-	Describes the feature of the command.
4	Send Forma	t	-	This is the format when sending commands.
	Conditions	Mode	REMOTE REMOTE TEACH TEACH	It indicates that black characters are usable, and gray characters are disabled. The outline character of REMOTE can be used after sending the command "PRM_ON:".  The outline character of TEACH can be used only for registration in the internal memory as a Teaching program.
5	that can be sent	Control state	READY * 1 BUSY * 2	It indicates that black characters are usable, and gray characters are disabled.
	00110		Send	Indicates that the command is a send-only command.
		Command	Send/Reply	Indicates that the command has a reply after sending the command.
		type	Teaching	It can be used only when registering in the internal memory as a
			registration only	Teaching program.
			No,	Format block number.
6	Details		Subcommand	The name of the subcommand.
			Description	Describes the subcommand.


<sup>\* 1</sup> All operations are accepted.



<sup>\*</sup> 2 Operation related to operation is denied.

# (2) Send / Reply command format

Explains the format of the send /reply command. With the send/reply command, there is a reply from the equipment after sending. Sending means sending a command from the computer to the equipment. Reply means that there is a reply from the device to the computer. The command format is not case sensitive..



**EDIT** 

**TEST** 

**SETUP** 

PRM

**INFO** 

#### \* Example 1 LOCAL mode

Send	Replay
MODE ?	LOCAL N

\* Example 2 PARAMETER mode (From "PRM\_ON;" command)

**TEACH** 

Send	Replay
MODE ?	REMOTE,PRMCMD

#### 16 ECHOR:


LOCAL REMOTE TEACH READY **BUSY** 

3

MOVEMENT TEST mode

SETUP mode

PARAMETER mode (From SETUP mode) **INFORMATION** mode



No,	Item	Reply	Description
4	Ctatus	0	Not Echo back
'	Status	1	Echo back
		1	2

No	ltem	Contents	Description
1	Reply format	-	Reply format
		No,	Block number for reply format.
2	Details	Reply	The reply content.
		Description	A description of the reply.
3	Reply example	-	This is a reply example.



It can be confirmed by the reply contents by the command "PRM\_GET: A15"

1	2	3	4	7	6	7	8	9	10	11	12	13	14	15	16

# 6-3. Echo back

When echo back is turned on, the command sent immediately before is echoed back. If the command sent from the equipment is a command error, "CMD ER" is returned.

Note: Please be sure to clear the receive buffer of the personal computer before sending the next command.

# (1) Command without a reply

Send format	MODE: LOCAL
No,	1
Reply format	MODE: LOCAL

No,	Description
1	Reply send command

# (2) Command with a reply

Send format MODE?

No,			1			2			3		
Reply format	Μ	0	D	Ε	?		L	0	С	Α	L

No,	Description
1	Reply send command
2	Space
3	Reply contents of "MODE?" Command reply

# (3) When a command error occurs due to a send command


Send format \*\*\*

No, 1
Reply format C M D E R

\* There is a space between "CMD" and "ER".

No,	Description
1	Reply command error







1 2 3 4 5 6 7 8 9 10 11 12 13 14 15	5 16	15	14	13	12	11	10	9	8	7	6	5	4	3	2	1

#### 06 SLEEP?

#### LOCAL REMOTE TEACH READY Send/Reply

Description Check the sleep mode status. For details on sleep, see "(24) Sleep function".

Send format SLEEP?

No, Reply format 0

No,	ltem	Reply	Description		
1	Status	Ctatus 0		Sleep mode OFF	
		Status 1	Sleep mode ON		

07 LOCAL REMOTE TEACH READY BUSY MODE: Send

Description Change the mode of the equipment.

No, Send format

MODE: LOCAL REMOTE TEACH

No,	ltem	Subcommand	Description
		LOCAL	Switch to LOCAL mode
1	Control	REMOTE	Switch to REMOTE mode
		TEACH	Switch to TEACH mode

#### 08 MODE? LOCAL REMOTE TEACH READY BUSY Send/Reply

Description Check the current mode of this equipment.

Send format M O D E ?

No. Re

Ο,			_
leply format	***		
		,	****
			-

No,1	No,2	Description			
	None	LOCAL mode			
1.0041	SETUP	SETUP mode			
LOCAL	PRM	PARAMETER mode (From SETUP mode)			
	INFO	INFORMATION mode			
	None	REMOTE mode			
	PRMCMD	PARAMETER mode (From "PRM_ON:" command)			
REMOTE	SETUP	SETUP mode			
	PRM	PARAMETER mode (From SETUP mode)			
	INFO	INFORMATION mode			
	JOGCMD	TEACH mode (JOG and command operation) $*$ 1			
	Ю	TEACH mode (I/O port operation) * 1			
	EDIT	TEACH edit mode (From "T_ON:" command)			
TEACH	TEST	MOVEMENT TEST mode			
	SETUP	SETUP mode			
	PRM	PARAMETER mode (From SETUP mode)			
	INFO	INFORMATION mode			
	LOCAL	REMOTE SETUP PRM INFO None PRMCMD SETUP PRM INFO JOGCMD IO EDIT TEACH TEST SETUP PRM			

<sup>\*</sup> 1 It can be confirmed by the reply contents by the command "PRM\_GET: A15" .

#### \* Example 1 LOCAL mode

Send	Reply
MODE ?	LOCAL

# \* Example 2 PARAMETER mode (From "PRM\_ON:" command)

Send	Reply
MODE ?	REMOTE,PRMCMD



09 F: LOCAL REMOTE TEACH READY BUSY Send

Description Change the feedback stage control method.

No,	Item	Subcommand	Description	
		1	First axis	
1	Axis	2	Second axis	
			W	W
0	Caustural	0	Open loop	
	Control	1	Closed loop	

10 FR: LOCAL REMOTE TEACH READY BUSY Send/Reply

Description Get the feedback stage control method.

No, 1
Send format FR: 1
2
W

No,	ltem	Subcommand	Description	Reply format block No,
	Axis	None	Depends on the parameter "AXIS Sel" * 1	None
1		1	First axis	1
		2	Second axis	1
		W	Both axis	1: First axis 2: Second axis

<sup>\* 1</sup> It can be checked by "PRM\_GET: GO1" command.

No, 1 O ,

No,	Item	Reply	Description
1,2	Chartain	0	Open loop
	Status	Status 1	Closed loop

11 C: LOCAL REMOTE TEACH READY BUSY Send

Description Change the current motor excitation state.

No, Send format C: 1 C 2 1 W

No,	Item	Subcommand	Description
	Axis	1	First axis
1		2	Second axis
		W	Both axis
0	Control	Countrial	Non-excitation
_		1	Excitation



12 CR: REMOTE TEACH READY Send/Reply

Description

Get the current motor excitation status.

No,

Send format



No,	ltem	Subcommand	Description	Reply format block No,		
		None	Depends on the parameter $^{\prime\prime}$ AXIS Sel $^{\prime\prime}$ * 1	None		
1	Axis	1	First axis	1		
		2	Second axis	1		
		W	Both axis	1 : First axis 2 : Second axis		

<sup>\* 1</sup> It can be checked by "PRM\_GET: G01" command.

No, Reply format

1		2
0	,	0
1		1

	No,	ltem	Reply	Description
	1,2	Status 0	0	Non-excitation
			Excitation	

13 BEC:

#### LOCAL REMOTE TEACH READY BUSY Send


Description

Forcibly completes positioning during operation. Cancel in case of error or emergency

SLOP.	
	Limit error
Resolvable errors * 2	Overflow error
	Teaching command error
Releasable state	Emergency

<sup>\* 2</sup> Command errors are excluded because they are cleared when a normal command is received.

No, Send format



No,	ltem	Subcommand	Description
	Axis    None   Depends on the parameter "AXI	Depends on the parameter "AXIS Sel" * 3	
1		First axis	
		Second axis	
		W	Both axis

<sup>\* 3</sup> It can be checked by "PRM\_GET: GO1" command.

#### <del>14</del> CMDR:

LOCAL REMOTE TEACH READY BUSY

Description

Gets the last normal command sent other than this command. If a command has not been sent before issuing this command, "\*" is returned.

Send format

CMDR:

No,

Reply format

\*\*\*\*

No,	ltem	Reply example	Description
1	Command	BEC:	Last command sent



2	2 3	4 5	6 7 8	9 10	11	12   13	3 14	15   16
15	ECHO:		LOCA	L REMOTE	TEACH	READY	BUSY	Send
	Description No, Send format		back after sending	command. For	details, se	e "6-3. Ech	oback".	
	No,	ltem	Subcommand			Description		
	1	Control	0			cho back Ol		
	·		1 1		E	Echo back O	N	<u>-</u>
16	ECHOR:		LOCA	L REMOTE	TEACH	READY	BUSY	Send/Rep
	Description Send forma		back setting status					
	No, Reply forma	t 0 1						
	No,	Item	Reply			Description		
	1	Status	0			cho back Ol		
			1 1			Echo back O	'IN	
(2) lr	nfomation co	ommands						
17	*IDN?		LOCA	L REMOTE	TEACH	READY	BUSY	Send/Rep
	Description Send format		ment information.					
	No,	1	2 3	4	5			
	Reply format	**** ,	**** , ****	, **** ,	****			
	No,	Item	Reply example			Description		
	1		SIGMAKOKI		\	Vendor nam	е	
	2		FC-114			Model name		
	3	Information	00000			Serial Number		
	4		00.00			uipment Nur		
	5		00.000		Fir	mware Vers	sion	
18	VENDOR:		LOCA	L REMOTE	TEACH	READY	BUSY	Send/Rep
	Description Send format		endor name of this	equipment.				
	No, Reply forma	1 ****						
	No,	ltem	Reply example			Description		

No,	ltem	Reply example	Description
1	Information	SIGMAKOKI	Vendor name



1   1	2 3 4 5	6 7 8	9 10	11 12	13   14	15   16
19	MODEL:	LOCA	L REMOTE	TEACH REAL	DY BUSY	Send/Reply
		Model name of this e	quipment.			
	No, 1 Reply format *****					
	No, Item	Reply		Descrip	otion	
	1 Information	FC-114 FC-414 FC-514		Model r	name	
20	SN:	LOCA	L REMOTE	TEACH REAL	DY BUSY	Send/Reply
	Description Get the Send format S N	Serial Number name	of this equipme	nt.		
	No, 1 Reply format ****					
	N.I. I.					
	No, Item	Reply example		Descrip		
	No, Item  1 Information	Reply example 0000		Descrip Serial Nu		
21			L REMOTE	Serial Nu	umber	Send/Reply
21	1 Information EN:	0000		Serial Nu	umber	Send/Reply
21	1 Information  EN:  Description Get the	LOCA		Serial Nu	umber	Send/Reply
21	1 Information  EN:  Description Get the Send format EN  No, 1  Reply format *****	OOOO  LOCA  Equipment Number of		Serial Nu TEACH REAL  nt.  Descrip	DY BUSY	Send/Reply
21	1 Information  EN:  Description Get the Send format EN  No, 1  Reply format *****	0000 LOCA Equipment Number o		Serial Nu	DY BUSY	Send/Reply
21	1 Information  EN:  Description Get the Send format EN  No, 1  Reply format *****	OOOO  LOCA  Equipment Number of	of this equipmen	Serial Nu TEACH REAL  nt.  Descrip	DY BUSY  potion  Number	Send/Reply  Send/Reply
	1 Information  EN:  Description Get the Send format EN  No, 1 *****  No, Item 1 Information  FV:  Description Get the	OOOO  LOCA Equipment Number of	of this equipmen	Serial Nu TEACH REAL  nt.  Description	DY BUSY  potion  Number	
	1 Information  EN:  Description Get the Send format EN  No, 1 *****  No, Item 1 Information  FV:  Description Get the	Reply example 00.00	of this equipmen	Serial Nu TEACH REAL  nt.  Description	DY BUSY  potion  Number	
	1 Information  EN:  Description Get the I E N  No, 1 *****  No, Item 1 Information  FV:  Description Get the I F V  No, 1 *****	Reply example 00.00	of this equipmen	Serial Nu TEACH REAL  nt.  Description	option  Number  DY BUSY  BUSY  BUSY	



## 23 RESO:

LOCAL REMOTE TEACH READY BUSY Send/Reply

Description

Get the Minimum resolution of this equipment.

No,

Send format

					1
R	Ε	S	0	:	/
					1
					2
					W

No,	Item	Subcommand	Description	Reply format block No,
		None	Depends on the parameter "AXIS Sel" * 1	None
1	1 Axis	1	First axis	1
		2	Second axis	1
		W	Both axis	1: First axis 2: Second axis

<sup>\* 1</sup> It can be checked by "PRM\_GET: G01" command.

No, Reply format

1		2
****	,	****

No,	Item	Reply	Description
		100	A discission and a discission
1,2	Status	50	Minimum resolution (unit: nm)
		10	(Griff: Filli)

24 LIMR:


LOCAL REMOTE TEACH READY BUSY Send/Reply

Description

Get the stage stroke, It is necessary to execute the "LIMG:" command in advance, If the "LIMG:" command has not been executed, " $\star$ " is returned.

No,

Send format



No,	ltem	Subcommand	Description	Reply format block No,
		None	Depends on the parameter $^{\prime\prime}$ AXIS Sel $^{\prime\prime}$ * 2	None
1	Axis	1	First axis	1
		2	Second axis	1
		W	Both axis	1 : First axis 2 : Second axis

<sup>\* 2</sup> It can be checked by "PRM\_GET: G01" command.

No, Reply format

1		2
****	,	****

	No,	ltem	Reply example	Description
Ī	1,2	Stage information	200000 * 3	Stroke

<sup>\* 3</sup> The minimum digit is the minimum resolution digit. For FC-114, it is 20.0000mm.



1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16

# 25 AN:

LOCAL REMOTE TEACH READY BUSY Send/Reply

Description

No,
Send format

Get the axis name,

1

1

2

No,	ltem	Subcommand	Description	Reply format block No,
		None	Depends on the parameter "AXIS Sel" * 1	None
1	Axis	1	First axis	1
		2	Second axis	1
		W	Both axis	1 : First axis 2 : Second axis

<sup>\* 1</sup> It can be checked by "PRM\_GET: GO1" command.

No, 1 2
Reply format \* , \*

No,	Reply example	Description
1	×	First axis name
2	Υ	Second axis name

26	UNT:	LOCAL	REMOTE	TEACH	READY	BUSY	Send/Reply

Description

No,
Send format

Get the unit.

1

U N T :

1

2

W

No,	ltem	Subcommand	Description	Reply format block No,
		None	Depends on the parameter "AXIS Sel" * 2	None
1	Axis	1	First axis	1
		2	Second axis	1
		W	Both axis	1 : First axis 2 : Second axis

<sup>\* 2</sup> It can be checked by "PRM\_GET: GO1" command.

No, 1 2 Reply format \* , \*

No,	ltem	Reply	Description						
		N	Nanometer						
		U	Micrometer						
1,2	Unit	M	Millimeter						
		D	Degree						
		Р	No unit (minimum digit is minimum resolution digit)						





# (3) Motion status Information commands

# 27 Q: LOCAL REMOTE TEACH READY BUSY Send/Reply


Description Get the coordinate values, the error, motion, and positioning status.

There are two types of reply format, and it choose by block number 1 of the send format. Reply format 1 is conventional.

Reply format 2 can choose the type and unit of coordinate value.

Also, it is possible to check the positioning status of each axis and all errors that occurre.


No, Send format



\* 1 If there is "No," to be set to "None", pad the specified No, additional command to the left. However, when No,1 is "None", No,2 and No,3 cannot be selected.

No,	ltem	Subcommand	Description					
1	Poply format	None	Format1					
* 1	Reply format	S	Format2					
		None	Reply in the unit set in parameter "UNIT Sel" $*2$					
		Ν	Nanometer					
2 * 1	Unit	U	Micrometer					
* 1		M	Millimeter					
		D	Degree					
		Р	No unit (minimum digit is minimum resolution digit)					
	Casualisasta valva	None	Reply with the value set in the parameter "Count Sei" $*3$					
3	Coordinate value type	E	Encoder values					
. 1	i, pe	С	Command values					

- \* 2 Can be confirmed with the reply by Command "PRM\_GET: AO3" or "PRM\_GET: AO4".
- \* 3 Can be confirmed with the reply by Command "PRM\_GET: G15" or "PRM\_GET: G16".



No,	Item	Reply	Description					
1 (First axis)	Ciaro	+	Plus					
3 (Second axis)	Sign		Minus					
2 (First axis) 4 (Second axis)	Coordinate value	00000001 * 4	Coordinate value					

<sup>\* 4</sup> The response content is an example. In the case of FC-114, it represents 100nm.



No,	Item	Reply	Description
		K	Normal (No error)
		1	Command error
		2	Scale error
		3	Limit stop
		4	Over speed error
		5	Overflow error
5	Status (Error)	6	Emergency stop
		7	Interpolator error
		8	Limit error
		9	System error
		Α	Slowdown sensor area
		В	Software Limit stop
		Т	TEACHING command error
		K	Normal stop
		М	During command move
		F	Out of the in-position range (After positioning is completed) $*1$
		G	During fine adjustment * 2, * 3
		Р	During electrical origin return
		0	During mechanical origin return
		С	CW side limit stop
6 (First axis)	Status (Motion)	W	CCW side limit stop
7 (Second axis)	Status (Motion)	В	CW side software limit stop
		V	CCW side software limit stop
		Α	CW side slowdown sensor area
		U	CCW side slowdown sensor area
		E	Error occurred
		Н	Motor is transitioning to excitation
		I	Motor is transitioning to non-excitation
		D	Disabled axis (Not set by parameter "AXIS Sel")
8	System reservation	0000	System reservation
9	Status	R	All axes have been positioned and no errors have occurred. * 1
9	(Positioning)	В	Positioning of all axes is incomplete or an error has occurred. $*2$

- \* 1 All operations are accepted. (This state is READY)
- \* 2 Refusing operation related to the move of the stage. (This state is BUSY)
- \* 3 It can be READY state by "BEC:" command.

No,	1	2	3		4	5	6		7		8		9		10		11		12	13
Reply format2	Ν	+	****	,	Ν	+	****	,	****	,	****	,	K	,	K	,	**	,	R	, R
	U	_	·		U	-							М		М				В	В
	M		•		М								F		F					
	D				D	İ							G		G					
	Р				Ρ	İ							Р		Р					
													0		0					
													С		С					
													W		W					
													В		В					
													<u>\</u>		$\overline{\vee}$					
													A		Ă					
													Û		Û					
													E		E					
													_		_					
													Н		Н					
													ı		Ш					
													D		D					

No,	ltem	Reply	Description		
		Z	Nanometer		
1 (First axis) 2 (Second axis)		U Micrometer			
	Unit	М	Millimeter		
2 (Second axis)		D	Degree		
		Р	No unit (minimum digit is minimum resolution digit)		
2 (First axis)	Sign	+	Plus		
5 (Second axis)	JELI	_	Minus		



No,	Item	Reply		Description							
3 (First axis) 6 (Second axis)	Coordinate value	****		Depends on the instruction	n unit. * 4						
			Bit	Reply: 1	Reply: 0						
			1 (MSB)	Normal	Error occurred						
			2	Command error	No occurred						
			3	Scale error	No occurred						
			4	Limit stop	Other						
			5	Over speed error	No occurred						
7 (First axis)	Status	1,0	6	Overflow error	No occurred						
8 (Second axis)	(Error)	1,0	7	Emergency stop	No occurred						
			8	Interpolator error	No occurred						
			9	Limit error	No occurred						
			10	System error	No occurred						
			11	Slowdown sensor area	Other						
			12	Software limit stop	Other						
			13 (LSB)	TEACHING command error	No occurred						
		K		Normal stop							
		М	During command move								
		F	Out of the in-position range (After positioning is completed)								
		G		During fine adjustment,	* 2, * 3						
		Р		During electrical origin	return						
		0		During mechanical origin							
	Status	С	CW side limit stop								
9 (First axis)		W	CCW side limit stop								
10 (Second axis)	(Motion)	В		CW side software limit	<u> </u>						
				CCW side software lim	<u>`</u>						
		Α	CW side slowdown sensor area								
		U		CCW side slowdown sen	sor area						
		E		Error occurred							
		Н		Motor is transitioning to e	excitation						
		I		Motor is transitioning to nor							
		D	Disabled axis (Not set by parameter "AXIS Sel")								
11	System reservation	0000		System reservatio							
12 (First axis)	Status	R		ave been positioned and no en							
13 (Second axis)	(Positioning)	В	Positioning of all axes is incomplete or an error has occur								

- \* 1 All operations are accepted. (This state is READY)
- \* 2 Refusing operation related to the move of the stage. (This state is BUSY)
- \* 3 It can be READY state by "BEC:" command.

### \* 4 Example 1 : When the coordinate value is 12.3456mm. (use FC-114)

Unit	Reply coordinate value
nm	12345600
um	12345.6
mm	12,3456
None	123456

# \* 4 Example 2 : When the coordinate value is 1.23455 $^{\circ}\,$ , (use FC-414)

Unit	Reply coordinate value								
Degree	1,23455								
None	123455								

Example: First axis is stopped at -12.34567mm, and the unit is nm. The second axes are operating at a position of 0.12345 mm, and the unit is  $\mu$  m. (use FC-514)

Send	Reply							
Q:	-1234567,+0012345,K,K,M,0000,B							
Q:S	N-1234567,U+123.45,1000000000000,100000000000,K,M,0000,R,B							
Q:SM	M-12.34567,M+0.12345,1000000000000,100000000000,K,M,0000,R,B							
Q:SE	N-1234567,U+123.45,1000000000000,100000000000,K,M,0000,R,B							
Q:SUC	U-12345,67,U+123.45,1000000000000,100000000000,K,M,0000,R,B							

<sup>\* 5</sup> The number of digits of the coordinate value varies depending on the unit.



## 28 SRQ:

## LOCAL REMOTE TEACH READY BUSY Send/Reply

Description

Get the status, Select Reply format in block No. 1 of Send format. Reply format1 can acquire the conventional method, and Reply format2 can acquire the positioning status and error status for each axis.

No, Send format



No,	ltem	Subcommand	Description					
4	Danily formers	None	Format1					
1 1	Reply format	S	Format2					

No,	1		2		3		4		5
Reply format1	1 2 3 4 5 6 7 8 9 A B	,	X M F G P O C S B > A D H T - D	,		,	**	,	RB

No,	Item	Reply	Description			
		K	Normal (No error)			
		1	Command error			
		2	Scale error			
		3	Limit stop			
		4	Over speed error			
	Charta an	5	Overflow error			
1	Status (Error)	6	Emergency stop			
	(LITOI)	7	Interpolator error			
		8	Limit error			
		9	System error			
		Α	Slowdown sensor area			
		В	Software Limit stop			
		Т	TEACHING command error			
		K	Normal stop			
		М	During command move			
		F	Out of the in-position range (After positioning is completed)			
		G	During fine adjustment, * 2, * 3			
		Р	During electrical origin return			
		0	During mechanical origin return			
		С	CW side limit stop			
2 (First axis)	Status	W	CCW side limit stop			
3 (Second axis)	(Motion)	В	CW side software limit stop			
		V	CCW side software limit stop			
		Α	CW side slowdown sensor area			
		U	CCW side slowdown sensor area			
		E	Error occurred			
		Н	Motor is transitioning to excitation			
		1	Motor is transitioning to non-excitation			
		D	Disabled axis (Not set by parameter "AXIS Sel")			
4	System reservation	0000	System reservation			

<sup>\* 1</sup> All operations are accepted. (This state is READY)



<sup>\*</sup> 2 Refusing operation related to the move of the stage. (This state is BUSY)

1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16

No,	ltem	Reply	Description
נ	Status	R	All axes have been positioned and no errors have occurred. * 1
5	(Positioning)	В	Positioning of all axes is incomplete or an error has occurred. * 2

No,	1		2		3		4		5		6		7
Reply format2	****	,	****	,	K	,	Κ	,	**	,	R	,	R
					М		Μ				В		В
					F		F					_	
					G		G						
					Р		Р						
					0		0						
					С		С						
					W		W						
					В		В						
					$\vee$		V						
					Α		Α						
					U		U						
					Е		Ε						
					П		Н						
					П		Τ						
					D		D						
								ı					

No,	Item	Reply		Description					
			Bit	Reply: 1	Reply: 0				
			1 (MSB)	Normal	Error occurred				
			2	Command error	No occurred				
			3	Scale error	No occurred				
			4	Limit stop	Other				
			5	·					
1 (First axis)	Status	1,0	6	Overflow error	No occurred				
2 (Second axis)	(Error)	1,0	7	Emergency stop	No occurred				
			8	Interpolator error	No occurred				
			9	Limit error	No occurred				
			10	System error	No occurred				
			11	Slowdown sensor area	Other				
			12	Software limit stop	Other				
			13 (LSB)	TEACHING command error	No occurred				
		K		Normal stop					
		М		During command move					
		F	Out of the in-position range (After positioning is completed) *						
		G		During fine adjustment, $*2, *3$					
		Р		During electrical origin	return				
		0		During mechanical origin					
		С		CW side limit stop	)				
3 (First axis)	Status	W		CCW side limit sto	þ				
4 (Second axis)	(Motion)	В		CW side software limit	t stop				
		V	CCW side software limit stop						
		Α		CW side slowdown sens	sor area				
		U		CCW side slowdown sen	sor area				
		E	Error occurred						
		Н	Motor is transitioning to excitation						
		1		Motor is transitioning to non-excitation					
		D	Disa	abled axis (Not set by param	eter "AXIS Sel")				
5	System reservation	0000		System reservatio					
6 (First axis)	Status	R	All axes have been positioned and no errors have occurred. * 1						
7 (Second axis)	(Positioning)	B		of all axes is incomplete or an	error has occurred, * 2				

- \* 1 All operations are accepted. (This state is READY)
- \* 2 Refusing operation related to the move of the stage. (This state is BUSY)
- \* 3 It can be READY state by "BEC;" command.

Example: When axis 1 is stopped and axis 2 is operating.

LAMITIPIC: WHICH MAIS	Example: Which axis it is stopped and axis 2 is operating.							
Send	Reply							
SRQ:	K,K,M,0000,B							
SRQ:S	100000000000,10000000000,K,M,0000,R,B							



### 29 F


# OCAL REMOTE TEACH READY BUSY Send/Reply

Description

Get the coordinate values. There are two types of reply format, and it choose by block number 2 of the send format. Reply format 1 is conventional. Reply format 2 can choose the type and unit of coordinate value.

No,

Send format



\* 1 If there is "No," to be set to "None", pad the specified No, additional command to the left. However, when No,2 is "None", No,3 and No,4 cannot be selected.

No,	ltem	Subcommand	Description	Reply format block No,
		None	Depends on the parameter "AXIS Sel" * 2	None
		1	First axis	1,2 (Reply format1)
		2	Second axis	1,2,3 (Reply format2)
1 ※1	Axis	W	Both axis	1,2: First axis 3,4: Second axis (Reply format1) 1,2,3: First axis 4,5,6: Second axis (Reply format2)
2	Reply	None	Format1	
<u> </u>	format	S	Format2	_
		None	Reply in the unit set in parameter "UNIT Sel" $*3$	
		N	Nanometer	
3	Unit	U	Micrometer	_
<b>※</b> 1	Oriit	M	Millimeter	
		D	Degree	
		Р	No unit (minimum digit is minimum resolution digit)	
4	Coordinate	None	Reply with the value set in the parameter "Count Sel" * 4	_
<b>※</b> 1	value type	E	Encoder values	
		С	Command values	

- \* 2 It can be checked by "PRM\_GET: G01" command.
- $\ast$  3 Can be confirmed with the reply by Command "PRM\_GET: A03" or "PRM\_GET: A04" .
- $\ast$  4 Can be confirmed with the reply by Command "PRM\_GET: G15" or "PRM\_GET: G16" .


No, Reply format1

1	2		3	4
+	****	,	+	****
_			_	

No,	Item	Reply	Description
1,3	Cign	+	Plus
1,0	Sign		Minus
2,4	Coordinate value	00000001 * 5	Coordinate value

\* 5 Contents are examples. For FC-514, it represents 10nm.

No, Reply format2



No,	Item	Reply	Description
		Ν	Nanometer
		U	Micrometer
1,4	Unit	М	Millimeter
		D	Degree
		Р	No unit (minimum digit is minimum resolution digit)



1	2	3	4	7	6	7	8	9	10	11	12	13	14	16

No,	ltem	Reply	Description
2,5	Sign	+	Plus
2,5	OBLI		Minus
3,6	Coordinate value	****	Varies by unit * 1

### \* 1 Example 1: When the coordinate value is 12,3456mm. (use FC-114)

Unit	Reply coordinate value
nm	12345600
um	12345.6
mm	12,3456
None	123456

### \* 1 Example 2 : When the coordinate value is 1.23455 $^{\circ}\,$ . (use FC-414)

Unit	Reply coordinate value
Degree	1,23455
None	123455

Example: First axis is stopped at -12.34567mm, and the unit is nm. The second axes are operating at a position of 0.12345 mm, and the unit is  $\mu$  m. (use FC-514)

Send	Reply
P:	-1234567,+0012345
P:S * 2	N-1234567,U+123,45
P:1S	N-1234567
P:2S	U+123.45
P:WSN *2	N-1234567,N+12345

<sup>\* 2</sup> The number of digits of the coordinate value varies depending on the unit.

30 ER: LOCAL REMOTE TEACH READY BUSY Send/Reply

Description

Get the error status. There are two types of reply format, and it choose by block number 1 of the send format. Reply format 1 is conventional. Reply format 2 can check all errors that occurred.

No, Send format



\*3 If there is "No," to be set to "None", pad the specified No, additional command to the left. However, when No,1 is "None", No,2 cannot be selected.

No,	Item	Subcommand	Description	Reply format2 block No,
		None	Depends on the parameter "AXIS Se1" $*$ 4	None
1	Axis	1	First axis	1
* 3	AXIS	2	Second axis	•
		W	Both axis	1: First axis 2: Second axis
2	David of favorant	None	Format1	
*3	Reply format	S	Format2	_

<sup>\* 4</sup> It can be checked by "PRM\_GET: GO1" command.



1	2	3	4	7	6	7	8	9	10	11	12	13	14	15	16

No, Reply format1

No,	Item	Reply	Description
		K	Normal (No error)
		1	Command error
		2	Scale error
		3	Limit stop
		4	Over speed error
	0	5	Overflow error
1	Status (Error)	6	Emergency stop
	(LITOI)	7	Interpolator error
		8	Limit error
		9	System error
		Α	Slowdown sensor area
		В	Software Limit stop
		Т	TEACHING command error

No, 1 2
Reply format2 \*\*\*\*\* , \*\*\*\*\*

No,	Item	Reply	Description					
			Bit	Reply: 1	Reply: 0			
			1 (MSB)	Normal	Error occurred			
			2	Command error	No occurred			
			3	Scale error	No occurred			
			4	Limit stop	Other			
	Status (Error)		5	Over speed error	No occurred			
1,2		1,0	6	Overflow error	No occurred			
1,∠			1,0	1,0	7	Emergency stop	No occurred	
			8	Interpolator error	No occurred			
			9	Limit error	No occurred			
					10	System error	No occurred	
			12	Software limit stop	Other			
			13 (LSB)	TEACHING command error	No occurred			

Example: Both axes are normal

Example: Beth axes are normal								
Send	Reply							
ER:	К							
ER:S	100000000000,10000000000							



1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16

# 31 STS:

# LOCAL REMOTE TEACH READY BUSY Send/Reply

Description
No,
Send format

Get the motion status,

1
2
W

	No,	Item	Subcommand	Description	Reply format block No,
Ī			None	Depends on the parameter "AXIS Sel" * 1	None
	1	A:-	1	First axis	4
1		Axis	2	Second axis	l l
١			W	Both axis	1: First axis 2: Second axis

<sup>\* 1</sup> It can be checked by "PRM\_GET: G01" command.

No,	1	2
Reply format	Κ,	K
	M	М
	F	F
	G	G
		Р
	0	0
	С	С
	W	W
	В	В
	$\vee$	$\vee$
	Α	Α
	U	U
	E	E
	Н	Н
	D	D

No,	Item	Reply	Description	
			K	Normal stop
		M	During command move	
		F	Out of the in-position range (After positioning is completed) * 2	
		G	During fine adjustment, * 3, * 4	
		Р	During electrical origin return	
		0	During mechanical origin return	
		С	CW side limit stop	
1 (First axis)	Status	W	CCW side limit stop	
2 (Second axis)	(Motion)	В	CW side software limit stop	
		V	CCW side software limit stop	
		Α	CW side slowdown sensor area	
		U	CCW side slowdown sensor area	
		E	Error occurred	
		Н	Motor is transitioning to excitation	
		1	Motor is transitioning to non-excitation	
		D	Disabled axis (Not set by parameter "AXIS Sel")	

<sup>\* 2</sup> All operations are accepted. (This state is READY)



 $<sup>\</sup>ast$  3 Refusing operation related to the move of the stage. (This state is BUSY)

<sup>\* 4</sup> It can be READY state by "BEC:" command.

1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 | 13 | 14 | 15 | 16

## 32 !

## LOCAL REMOTE TEACH READY BUSY Send/Reply

Description

Get the positioning status. There are two types of reply format, and it choose by block number 1 of the send format. Reply format 1 is conventional. Reply format 2 can check the positioning status of each axis.

No,

Send format



\* 1 If there is "No," to be set to "None", pad the specified No, additional command to the left. However, when No,1 is "None", No,2 cannot be selected.

No,	Item	Subcommand	Description	Reply format block No,
		None	Depends on the parameter "AXIS Se1" $*$ 2	None
1	Anda	1	First axis	1
* 1	Axis	2	Second axis	l l
		W	Both axis	1: First axis 2: Second axis
2	Reply format	None	Format1	
*1	neply format	S	Format2	_

<sup>\* 2</sup> It can be checked by "PRM\_GET: GO1" command.

No,

Reply format1



I	No,	Item	Reply	Description		
Ī	1	Status	R	All axes have been positioned and no errors have occurred. * 3		
ĺ		(Positioning)	В	Positioning of all axes is incomplete or an error has occurred. * 4		

No, Reply format2



No,	Item	Reply	Description
4.0	Status	R	All axes have been positioned and no errors have occurred. * 3
1, ∠	(Positioning)	В	Positioning of all axes is incomplete or an error has occurred. * 4

<sup>\* 3</sup> All operations are accepted. (This state is READY)

Example: When axis 1 is stopped and axis 2 is operating.

Send	Reply
!:	В
!:S	R,B
!:1S	R
!:WS	R,B



<sup>\*</sup> 4 Refusing operation related to the move of the stage. (This state is BUSY)

## (4) Commands related to the origin

33 H: LOCAL REMOTE TEACH READY BUSY Send

Description

Performs machine home position return. For details, see "9. Home Return" . When the axis to be executed is non-excitation, the command error occurs, If  $\H$ ,  $R\H$  is specified at the end of the command, a positioning completion reply is sent for each axis.

No, Send format

OI IC	<i>i</i> Oi	u ic	COI	IIIIIC	ai ia,
	1		2	3	3
Н		• •	/		
	0		1	,	R
	1		2		
	2		V		
	3				
	4				
	4	15 41-			

\* 1 If there is a number to be set to "None", pad the additional command with the specified number to the left.

No,	Item	Subcommand	Description					
		None	Depends on the parameter "ORG Mode Sel" * 2					
		0	ModeO					
1		1	Mode1					
* 1	Mode	2	Mode2					
		3	Mode3					
		4	Mode4					
		None	Depends on the parameter "AXIS Sel" * 3					
2	Axis	1	First axis					
* 1	AXIS	2	Second axis					
		W	Both axis					
2		None	Do not request a reply					
3 * 1	Reply request	,R	After positioning is complete, return "1" for the first axis and "2" for the second axis.					

<sup>\* 2</sup> Can be confirmed with the reply by Command "PRM\_GET: A30" or "PRM\_GET: A31".

34 Z: LOCAL REMOTE TEACH READY BUSY Send

Description Perform electrical home return. For details, see "9. Home Return". When the axis to be executed is non-excitation, the command error occurs. If  $\H$ ,  $\H$  is specified at the end of

the command, a positioning completion reply is sent for each axis.

No,

Send format



If there is a number to be set to "None", pad the additional command with the specified number to the left.

No,	Item	Subcommand	Description			
		None	Depends on the parameter "AXIS Sel" * 5			
1	Axis	1	First axis			
* 4		2	Second axis			
		W	Both axis			
2	Reply request	None	Do not request a reply			
2 * 4		,R	After positioning is complete,			
		,, ,	return "1" for the first axis and "2" for the second axis.			

<sup>\* 5</sup> It can be checked by "PRM\_GET: G01" command.



<sup>\* 3</sup> It can be checked by "PRM\_GET: G01" command.

Description Executes the electrical origin setting (zero set). For details, see "9.

Executes the electrical origin setting (zero set). For details, see "9. If the motor is in the demagnetized state, a command error will occur if the specified axis is demagnetized and executed.

No, 1
Send format R: 1
2
W

No,	Item	Subcommand	Description
		None	Depends on the parameter "AXIS Sel" * 1
1	Axis	1	First axis
'		2	Second axis
		W	Both axis

<sup>\* 1</sup> It can be checked by "PRM\_GET: G01" command.

36 LIMG: LOCAL REMOTE TEACH READY BUSY Send

Description Executes the stage stroke detection operation. This value can be obtained by "LIMR:"

command. When the axis to be executed is non-excitation, the command error occurs. When  $\H, R\H$  is added to the end of the send format, the axis number is replied when it is

completed.

No, 1 2
Send format L I M G : 1 , R
2 W

\*2 If there is a number to be set to "None", pad the additional command with the specified number to the left.

No,	Item	Subcommand	Description			
		None	Depends on the parameter "AXIS Sel" * 3			
1	Axis	1	First axis			
* 2		2	Second axis			
		W	Both axis			
2	Reply request	None	Do not request a reply			
2 * 2		,R	After positioning is complete,			
		,1 1	return "1" for the first axis and "2" for the second axis.			

<sup>\* 3</sup> It can be checked by "PRM\_GET: G01" command.

## (5) Commands related to the motion

37 L: LOCAL REMOTE TEACH READY BUSY Send

Description Execute stop and emergency stop.

No,

Send format



No,	ltem	Subcommand	Description
		None	Depends on the parameter "AXIS Sel" * 4
		1	First axis
1	Axis	2	Second axis
		W	Both axis
		Е	Emergency stop

<sup>\* 4</sup> It can be checked by "PRM\_GET: GO1" command.



<sup>\* 5</sup> Can be canceled with Command "BEC:".

1	2	3	4	7	6	7	8	9	10	11	12	13	14	15	16

#### 38 ACC: REMOTE TEACH READY Send

Description

Set the acceleration  $\/$  deceleration time. However, if the instrument is turned off, reset, restarted, or the GENERAL parameter is changed, the set contents are discarded. When without setting of this value, executing the "A:", "M:" or "JG:" commands, the setting value of the parameter "Acc Time" is applied for the acceleration and deceleration time. The setting value of the parameter "Acc Time" can be checked by the reply of "PRM\_GET: A18" and "PRM\_GET: A19" command.

No, Send format

/ (IO ai	I			<u>ا ۱</u>	.,	OOH	irriaria.
			1	2	3		4
AC	С	:	1		***	,	***
			2				
			W				

No,	Item	Subcommand	Description	Send format block No,	
		1	First axis		
1	1 Axis		Second axis	3	
		W	Both axis	3: First axis 4: Second axis	
2	Space	Space sign	Space	-	
3	Time ***		Set in milliseconds		
4	rime	***	(10~2000)	_	

#### 39 ACCR: LOCAL REMOTE TEACH READY BUSY Send/Reply

Description

Acquires the set acceleration / deceleration time.

No,

Send format ACCR

	No,	Item	Subcommand	Description	Reply format block No,
			None	Depends on the parameter "AXIS Sel" * 1	None
	1 Axis	1	First axis	4	
		2	Second axis	l	
			W	Both axis	1: First axis 2: Second axis

<sup>\* 1</sup> It can be checked by "PRM\_GET: G01" command.

No,

Reply format \*\*\* \*\*\*

No,	Item	Reply example	Description
1, 2	Time	100	In this case is 100 ms.



1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16

#### 40 LOCAL REMOTE TEACH READY D: Send

Description Set the movement speed of the stage. However, when the power is turned off, RESET,

restart or the GENERAL parameter is changed, this value will be lost. When without setting of this value, executing the "A:", "M:" or "JG:" commands, the setting value of the parameter "Max Speed" is applied for the movement speed. The setting value of the parameter "Max Speed" can be checked by the reply of "PRM\_GET: GO2" and "PRM\_GET:

No, Send format

G03" command.					
	1	2	3	4	5
D :	1	F	****	F	****
	2	Р		Р	
	W	Ν		N	
,		$\subset$		U	
		М		М	
		D		D	

No,	Item	Subcommand	Description	Send format block No,
		1	First axis	2
1	Axis	2	Second axis	3
		W	Both axis	3: First axis 4: Second axis
		F, P	No unit / secon	d * 1
		N	Nanometer / :	sec
2, 4	Unit	U	Micrometer / :	sec
		М	Millimeter / s	ec
		D	Degree / se	С
3, 5	Speed	****	The setting speed is from the mir to the parameter "Max	

<sup>\* 1</sup> The minimum digit of the operation speed is the digit of the minimum resolution.

### \*2 Example: Operating speed setting (when FC-114 is used and 12,3456 mm / sec)

Unit	Speed
No unit / sec	123456
nm/sec	12345600
μm/sec	12345.6
mm/sec	12,3456

Example: First axis side operation speed is 12.34565mm / sec, second axis side operation speed is  $500\,\mu$ m / sec, and when no unit is specified for both axes. (When using FC-414)

Axis	Command
First axis only	D:1F1234565
Second axis only	D:2F50000
Both axis	D:WF1234565F50000

Example: When the First axis side operation speed is 1.23456mm / sec and Millimeter is specified, and the second axis side operation speed is 500  $\mu$ m / sec and Micrometer is specified. (When using FC-514)

side oper	side operation speed is 500 $\mu$ m / sec and which othere is specified, (which dishing 1 $\circ$ 514)			
Axis	Command			
First axis only	D:1M1,23456			
Second axis only	D:2U500			
Both axis	D:WM1.23456U500			



# 41 DR: LOCAL REMOTE TEACH READY BUSY Send/Reply

Description

Acquires the set operation speed.

No, Send format

1 2 D R : / 1 F 2 P W N U M

\* 1 If there is a number to be set to "None", pad the additional command with the specified number to the left.

No,	Item	Subcommand	Description	Reply format block No,
		None	Depends on the parameter "AXIS Sel" * 2	None
1	A:-	1	First axis	1.2
* 1	Axis	2	Second axis	1, 2
		W	Both axis	1, 2: First axis 3, 4: Second axis
	Unit	None	Depends on the parameter	r "UNIT Sel" * 3
		F, P	No unit (minimum digit of operation speed	is digit of minimum resolution)
2		N	Nanometer	-
* 1		U	Micrometer	-
		M	Millimeter	
		D	Degree	

st 2 It can be confirmed by the reply contents by Command "PRM\_GET: GO1" .

No, Reply format

1	2		3	4
F	****	,	H	****
Р			Р	
Ν			Ν	
U			U	
М			М	
D			D	

No,	Item	Reply	Description
		F, P	No unit (minimum digit of operation speed is digit of minimum resolution)
		Ν	Nanometer
1, 3	Unit	U	Micrometer
		М	Millimeter
		D	Degree
2, 4	Speed	***	The setting speed is from the minimum resolution / sec to the parameter "Max Speed" . * 4

 $\ast$  4  $\,$  Example: Operating speed setting (when FC-114 is used and 12.3456 mm / sec)

Unit	Reply speed
No unit/sec	123456
nm/sec	12345600
um/sec	12345.6
mm/sec	12,3456

 $*\,4$   $\,$  Example: Operating speed setting (when FC-414 is used and 12.3455  $^{\circ}$   $\,$  / sec)

Unit	Reply speed
No unit/sec	123455
°/sec	1,23455

Example: When the First axis side operation speed is 1.23456mm / sec and Millimeter is specified, and the second axis side operation speed is 500  $\mu$ m / sec and Micrometer is specified. (When using FC-514)

side oper	side operation speed is 500 $\mu$ m / sec and Micrometer is specified. (When using FC-514)								
Send	Reply								
DR:	N123456,U500								
DR:M	M1 23456.M0.5								



st 3 It can be confirmed by the reply contents by Command "PRM\_GET: A03" and "PRM\_GET: A04" .

1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16

42 A: LOCAL REMOTE TEACH READY BUSY Send

Description

Set the absolute movement coordinate value. Set the coordinate value from the origin (position where the coordinate value is zero). This command alone will not work. To make it work, execute the command "G" or "GN:" . See commands "G" and "GN:" for details.

No, Send format

	JUU	io ii	IC COITII	Hai	iu v	a 01 C
1	2	3	4	5	6	7
1	+	Р	****	+	Ρ	****
2	-	Z		-	Ν	
W		$\supset$			$\supset$	
		М			М	
		D			D	
	1 1 2 W					1 2 3 4 5 6 1 + P **** + P 2 - N - N W U M D

No,	Item	Subcommand	Description	Send format block No,				
		1	First axis	2.2.4				
1	Axis	2	Second axis	2, 3, 4				
'	430	W	Both axis	2, 3, 4: First axis 5, 6, 7: Second axis				
2,5	Ciero	+	Plus					
2, 5	Sign	_	Minus					
		Р	No unit (minimum digit of coordinate value	is digit of minimum resolution)				
		N	Nanometer					
3, 6	Unit	U	Micrometer					
		М	Millimeter					
		D	Degree					
4, 7	Coordinate value	****	Coordinate value * 1 (Setting range dep	ends on the connected stage)				

### \* 1 Moving coordinate value setting example (when FC-114 is used and 12.3456mm)

Unit	Coordinate value
No unit	123456
nm	12345600
um	12345.6
mm	12,3456

### \* 1 Moving coordinate value setting example (when FC-414 is used and it is 1.23455 $^{\circ}$ )

Unit	Coordinate value
No unit	123455
0	1,23455

Example: When using FC-514 and moving the first axis (Unit: nm) from the origin (position where the coordinate value is zero) to -1.23456mm and the second axis (Unit: um) to + 0.5mm

IS ZEIO/ II	5 1,23 <del>4</del> 3011111 and	is zero/ to 1,25450mm and the second axis (or it: diff) to 10,5mm								
Send example	Send order	Command								
Evenente 1	1	A:W-N123456+U500								
Example1	2	G								
Example2	1	A:W-N123456+U500								
Example2	2	GN:W								



1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16

43 M: LOCAL REMOTE TEACH READY BUSY Send

Description Set the relative movement distance. This command alone will not work. To make it work, execute the command "G" or "GN:" . See commands "G" and "GN:" for details,

No, Send format

exe	Cuu		e 00	ווווכ	iariu C	a 0		11 10 . 36	;
		1	2	3	4	5	6	7	
М		1	+	Р	****	+	Р	****	
		2	_	Ν		_	Ν		
		W		U			$\Box$		
				М			М		
				D			D		

No,	Item	Subcommand	Description	Send format block No,				
		1	First axis	0.0.4				
1	Axis	2	Second axis	2, 3, 4				
	, , ,	W	Both axis	2, 3, 4: First axis				
				5, 6, 7: Second axis				
2,5	Sign	+	Plus					
2, 5	ושכ	_	Minus					
		Р	No unit (minimum digit of moving distance	is digit of minimum resolution)				
		N	Nanometer					
3, 6	Unit	U	Micrometer					
			M	Millimeter				
		D	Degree					
4, 7	Distance	****	Set the Movement distance * 1					
4, 7	Distance	<u>ተተተተ</u>	(The range depends on the connected stage)					

\* 1 Movement distance setting example (when FC-114 is used and 12.3456mm).

Unit	Movement distance
No unit	123456
nm	12345600
um	12345.6
mm	12,3456

\* 1 Movement distance setting example (when FC-414 is used and it is 1.23455  $^{\circ}$  ).

Unit	Movement distance 123455							
No unit	123455							
0	1,23455							

Example: When using FC-514 and moving from the current position to the first axis (Unit: nm) -12,34567mm and the second axis (Unit: um) to + 0,5mm,

30001 id a	XI3 (OF III. GITT) (O · O.	911111,
Send example	Send order	Command
Evenenda 1	1	M:W-N1234567+U500
Example1	2	G
Evenenda O	1	M:W-N1234567+U500
Example2	2	GN:W



#### 44 REMOTE TEACH READY G Send

The operation set by Command "A:" and "M:" starts. After execution, the values set with Description

the commands "A:" and "M:" are discarded. When ", R" is specified at the end, positioning

completion reply is sent for each axis.

No, Send format

45

No,	Item	Subcommand	Description
		None	Do not request a reply
1	Reply request	,R	After positioning is complete, return "1" for the first axis and "2" for the second axis.

#### LOCAL REMOTE TEACH READY GN: Send

The operation set by Command "A:" and "M:" starts. The values set by post-execution commands "A:" and "M:" are retained. When  $\tilde{\ }$ ,  $R^{\tilde{\ }}$  is specified at the end, positioning Description

completion reply is sent for each axis.

No, Send format GN:

No,	Item	Subcommand	Description				
		None	Depends on the parameter "AXIS Sel" * 1				
1	1 Axis	1	First axis				
'		2	Second axis				
		W	Both axis				
		None	Do not request a reply				
2	Axis 1 2 W	After positioning is complete, return $"1"$ for the first axis and $"2"$ for the second axis,					

 $<sup>*\,1</sup>$  It can be confirmed by the reply contents by Command "PRM\_GET: GO1" .

Example: When FC-114 is used and Command "A:" "M:" "G" "GN:" is used

(Before sending, confirm that the positioning status of the movement target axis is READY with the command "Q:", "SRQ:", "!:" And send it.)

No	Canal agreement			nt setting lue	Coordinate value		
INO	Send command	Description	First	Second	First	Second	
		Set the electrical origin for both axes	axis	axis	axis	axis	
1	R:W	(zero set)	None	None	Omm	Omm	
2	A:1-P100000	First axis: Absolute motion (A) -10mm	A -10mm	1	<b>↑</b>	<b>↑</b>	
3	GN:1	Move first axis	<b>↑</b>	<b>↑</b>	-10mm	<b>↑</b>	
4	GN:2	Command error	1	<b>↑</b>	<b>↑</b>	$\uparrow$	
5	A:2+P100000	Second axis: Absolute motion (A) +10mm	1	A +10mm	<b>↑</b>	$\uparrow$	
6	GN:2	Move second axis	1	1	<b>↑</b>	+10mm	
7	GN:W	It does not work because it is already moving to the destination.	1	1	<b>↑</b>	<b>↑</b>	
8	M:W-P10000-P10000	First axis: Relative motion (M) -1mm Second axis: Relative motion (M) -1mm	M -1mm	M -1mm	<b>↑</b>	<b>↑</b>	
9	GN:1	Move first axis	1	1	-11mm	<b>↑</b>	
10	GN:1	Move first axis	1	1	-12mm	<b>↑</b>	
11	GN:2	Move second axis	1	1	<b>↑</b>	+9mm	
12	GN:W	Move both axes	1	1	-13mm	+8mm	
13	A:1-P10000	First axis: Absolute motion (A) -1mm	A -1mm	1	<b>↑</b>	<b>↑</b>	
14	GN:1	Move first axis	1	1	-1mm	<b>↑</b>	
15	GN:W	Move both axes	1	1	<b>↑</b>	+7mm	
16	G	Move both axes	Discarded	Discarded	<b>↑</b>	+6mm	
17	GN:W	Command error	None	None	1	<b>↑</b>	



1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16

46 GC:

LOCAL REMOTE TEACH READY

Description

Deletes the value set by Command "A:" and "M:".

No, Send format G C : 1 2

No,	Item	Subcommand	Description
No,		None	Depends on the parameter "AXIS Sel" * 1
1	A:-	1	First axis
'	Axis	2	Second axis
		W	Both axis


<sup>\* 1</sup> It can be confirmed by the reply contents by Command "PRM\_GET: GO1".

47 GR: LOCAL REMOTE TEACH READY BUSY Send/Reply

Description

Get the coordinate value and movement distance of the set command "A:", "M:".

No, Send format




\*2 If there is a number to be set to "None", pad the additional command with the specified number to the left.

No,	Item	Subcommand	Description	Reply format block No,
		None	Depends on the parameter "AXIS Sel" $*$ 3	None
1		1	First axis	1, 2, 3, 4
*2	Axis	2	Second axis	1, 2, 3, 4
		W	Both axis	1, 2, 3, 4: First axis
		VV	DOLIT AXIS	5, 6, 7, 8: Second axis
		None	Depends on the parameter	"UNIT Sel" * 4
		Р	No unit (minimum digit of moving distance	is digit of minimum resolution)
2	l liait	P N	No unit (minimum digit of moving distance Nanometer	
2 * 2	Unit			
2 * 2	Unit		Nanometer	

<sup>\* 3</sup> It can be confirmed by the reply contents by Command "PRM\_GET: GO1".

No, Reply format



No,	Item	Reply	Description				
1, 5	Command	А	Absolute ( "A:" command)				
1,5	Command	A Absolute ("A:" command)  M Relative ("M:" command)  P No unit (minimum speed digit is minimum resolution di  N Nanometer  U Micrometer  M Millimeter  D Degree  + Move in the plus direction  Move in the minus direction	Relative ("M:" command)				
		Р	No unit (minimum speed digit is minimum resolution digit)				
	6 Unit	Ν	Nanometer				
2, 6		U	Micrometer				
					M	Millimeter	
		D	A Absolute ("A:" command)  M Relative ("M:" command)  P No unit (minimum speed digit is minimum resolution digit)  N Nanometer  U Micrometer  M Millimeter  D Degree  + Move in the plus direction  Move in the minus direction  ***** Coordinate value or moving distance				
3, 7	Sign	+	Move in the plus direction				
3, 7	B D Degree  + Move in the plus direction  - Move in the minus direction						
4, 8	Coordinate value	****	Coordinate value or moving distance				
	Common	*	When not set				



Send

<sup>\* 4</sup> It can be confirmed by the reply contents by Command "PRM\_GET: AO3" and "PRM\_GET: AO4".

4 14

48 JG: REMOTE TEACH READY Send

Description Executes movement by specifying the number of pulses (1 pulse = minimum resolution).

When ", R" is specified at the end, positioning completion reply is sent for each axis.

No, Send format JG \*\*\*


No,	Item	Subcommand	Description					
4	Axis	1	First axis					
'	<b>A S S</b>	2	Second axis					
2	Space	Space sign	Space					
3	Ciero	+	Move in the plus direction					
3	Sign -	Move in the minus direction						
4	Coordinate value	***	Set movement command value (1 $\sim$ 1000)					
	Reply	None	Do not request a reply					
5	request	,R	After positioning is complete, return "1" for the first axis and "2" for the second axis.					

49 JY: REMOTE TEACH READY BUSY Send

Description


Executes the operation without specifying the movement distance. The operation continues until a limit sensor input, operation stop command "L:" or speed stage "O" is transmitted. The operation speed can be selected from nine parameters "Jy Speed 1" to "Jy Speed 9" and "Jy Speed H". When operating at speed stages 1-9, the operating speed can be changed in stages 1-9. When operating at speed stage H, speed stages 1 to 9 cannot be selected. The control status during operation depends on the settings of the parameters "Jy Cont" and "Stage Cont Type". When ", R" is specified at the end, positioning completion reply is sent for each axis.

No, Send format



No,	Item	Subcommand	Description				
		1	First axis				
1	Axis	2	Second axis				
		W	Both axis				
2	2 Space Space sig		Space				
	Ciene	+	Move in the plus direction				
3	Sign	_	Move in the minus direction				
		0	Stop movement when executing JY command				
4	Speed stage	1~9	Set 1 to 9 (depends on parameter "Jy Speed 1 to 9")				
		Н	Set H (depends on parameter "Jy Speed H")				
	Reply	None	Do not request a reply				
5	request	,R	After positioning is complete, return "1" for the first axis and "2" for the second axis.				





	No,	Item	Item Subcommand Description					
	1	Number	***	Location registration designation number $(1 \sim 20)$				
		Reply request	None	Do not request a reply				
	2		,R	After positioning is complete, return "1" for the first axis and "2" for the second axis.				



1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16

## (7) General purpose I / O commands

4 I: LOCAL REMOTE TEACH READY

Description Get the status of the general-purpose input port.

Send format :

No, 1
Reply format \*

No,	ltem	Reply	Description
1	Status	*	Input state number

BUSY

Send/Reply

### Input status

Input state number * 1	Input1 (3pin)	Input2 (28pin)	Input3 (4pin)
0	OFF	OFF	OFF
1	ON	OFF	OFF
2	OFF	ON	OFF
3	ON	ON	OFF
4	OFF	OFF	ON
5	ON	OFF	ON
6	OFF	ON	ON
7	ON	ON	ON

<sup>\* 1</sup> See "(4) General-purpose I / O" for ON / OFF status and input circuit.

55 O: LOCAL REMOTE TEACH READY BUSY Send

Description Get the status of the general-purpose output port.

No, 1 Send format O: \*


No,	ltem	Subcommand	Description
1	Control	*	Output state number

### Output status

Output state number * 1	Output1 (1pin)	Output2 (26pin)	Output3 (2pin)
0	OFF	OFF	OFF
1	ON	OFF	OFF
2	OFF	ON	OFF
3	ON	ON	OFF
4	OFF	OFF	ON
5	ON	OFF	ON
6	OFF	ON	ON
7	ON	ON	ON

<sup>\* 1</sup> See "(4) General-purpose I / O" for ON / OFF status and input circuit.







1 2	2 3	4 5	7 8	9 10 11 12 13 14 15 16
62	TCR:		LOCA	AL REMOTE TEACH READY BUSY Send/Reply
	Description Send format		urrent teaching cha	annel.
	No, Reply forma	1 *		
	No,	Item	Reply	Description
	1	Channel	*	Current channel
63	TQ:		LOCA	AL REMOTE TEACH READY BUSY Send/Reply
	Description Send format		he status related t	to teaching.
	No,	1 2	3	4
	Reply forma			****
		M		
		P		
		0		
		+		
		R		
		E		
	No,	ltem	Reply	Description
			K M	Stopped * 1  During move * 1
			P	Paused
			0	During move (executed line by line) * 1
	1	Status	ı	General I/O operation status
			Т	Teaching edit mode
			R	loading teaching registration contents * 2
		Ol1	E *	Teaching command error
	3	Channel Line	***	Current channnel Current line number
	4	Command	****	Current command
			CH IF" is set to "I / O"	, the motion status can not be checked.
	* 2 All operat	ions are disabled (	until reading is comp	olete.
64	TG:		LOCA	AL REMOTE TEACH READY BUSY Send
	D: "	C+ · · · ·		alested decreed to the decree of the decree
	Description Send forma		execution of the s	selected channel. It also resumes execution when paused.
	Seria forma	ı <u>ı u</u> .	J	
65	TP:		1.00	AL REMOTE TEACH READY BUSY Send
00	11.		LOOP	LE LENGTE LE COLL TIEMET DOOT OFFICE
	Description	If it is sen	t during teaching,	it will pause. If you want to resume, execute the command
		″TG:″.	7	
	Send forma	T   P   :	J	
			LOCA	AL REMOTE TEACH READY BUSY Send
66	TO:		LOGA	AL HEMOTE TEACH THEADT BOOT Send
66		Eyecutes		
66	TO: Description		the contents line	by line in the paused state. If the stage is operating, this discarded until positioning is completed.



14 67 TEACH READY BUSY TL: Send Description Stops teaching and returns the line number to the first line. No. Send format No, Item Subcommand Description None Stop teaching 1 Axis Perform emergency stop \* Ε \* It can be canceled with Command "BEC:" . 68 TR: TEACH READY BUSY Send/Reply Description Check the registration status of the teaching channel. Send format TR No, Subcommand Item Description Set all channels None 1 Channel 1~5 Set each channel 3 5 No, 2 4 Reply format 0 0 0 Ο 0 1 1 1 No, Description Item Reply 1 Channel 1 or specified channel 2 Channel2 O: Unregistered 3 Status O, 1 Channel3 1: Registered 4 Channel4 5 Channel5 69 TFR: LOCAL REMOTE TEACH READY BUSY Send/Reply Description Check the number of loops currently being executed. It can be used only during teaching execution. When the reply content is O, it indicates that the target loop is not executed or unused. No, Send format TFR

I	No,	ltem	Subcommand	Description
Ī	1	1	None	Set all loops
	1	Loop	1~9	Set each loop



1	2	3	4	5		6		7	7		8		9		10		11		12		13	14	15	16
	No Re	o, eply forr	nat	1 **	,	2 **	,	3	,	4	,	5	,	6	,	7	,	8	,	9				

No,	Item	Reply	Description	1
1			Loop1 or Specified loop	
2			Loop2	
3			Loop3	
4			Loop4	Current loop count
5	Count of loops	**	Loop5	•
6			Loop6	(1~99999)
7			Loop7	
8			Loop8	
9			Loop9	

70 TM: LOCAL REMOTE TEACH READY Send

Description Set the teaching monitor. By turning this setting ON, the executed command is returned to

the interface set by the parameter "I / F Sel". Reply format is Reply format1 except "FE:", and "FE:" is Reply format2. The set value of the parameter "I / F Sel" can be confirmed by

the reply contents by the command "PRM\_GET: G24".

No, Send format TM О

No,	ltem	Subcommand	Description
1	Caustural	0	Monitor setting OFF
ı	Control	1	Monitor setting ON

No, Reply format1 \*\*\* \*\*\*\*

	No,	Item	Reply	Description
Ī	1	Line	***	Execution line number (001 $\sim$ 200)
[	2	Space	Space sign	Space
	3	Command	****	Execution command

Reply format2 \*\*\* \*\*\*\*

No,	Item	Reply	Description
1	Line	***	Execution line number (001 $\sim$ 200)
2	Space	Space sign	Space
3	Command	****	Execution command
4	Space	Space sign	Space
5	brackets	[	Use as separator
6	Count of loops	**	Current loop count (1 ~ 99999)
7	brackets	]	Use as separator

Example 1 When the line number is the second line and the execution command is "M:".

Auto reply
002,M:1+P10000

Example 2: When the line number is the second line, the execution command is "FE:", and the loop count is the third.

Auto reply 002 FE:1 [3]



14 71 TMR: TEACH READY REMOTE BUSY Send/Reply Description Gets the teaching monitor setting status. Send format T M R : No, 1 Reply format 0 1 No, ltem Reply Description Ο Monitor setting is OFF 1 Status 1 Monitor setting is ON 72 TNR: REMOTE TEACH READY BUSY Send/Reply Description Gets the teaching line number currently being executed or in standby. Send format T N R : No, 1 Reply format \*\*\* Item Reply No, Description Execution line number Line \*\*\* 1 73 TACR: LOCAL REMOTE TEACH READY BUSY Send/Reply Description Gets the teaching command currently being executed or in standby. Send format TACR: Reply format \*\*\*\*


No,	ltem	Reply	Description
1	Command	****	Execution command



(9) Commands related to teaching registration only LOCAL REMOTE TEACH READY Register Description Set the loop operation to repeat the specified range. The loop range is between the same loop numbers from this command to the command "FE:" . You can put another loop inside the loop. For a description example, see "(1) Transmission example 1". Be sure to use "FE:" when using this command. If "FE:" is not registered, unintended operation may occur. No, Send format \*\*\* FS No, Item Reply Description 1 \* Set loop  $(1 \sim 9)$ Loop 2 Space Space sign Space 3 Count of loops \*\*\* Loop count  $(1 \sim 99999)$ 75 REMOTE TEACH READY FE: Register Set the end of the loop range. The loop range is from Command "FS:" to the same loop Description number of this command. You can put another loop inside the loop. For a description example, see "(1) Transmission example 1". When using this command, be sure to use "FS:". If "FS:" is not registered, unintended operation may occur. No, \* Send format | F | E No, Item Subcommand Description Set loop  $(1 \sim 9)$ Loop 76 END: LOCAL REMOTE TEACH Register Description Set the last teaching line. Execution ends at this line. Send format E N D : 77 T: LOCAL REMOTE TEACH Description Set the wait time during teaching execution in increments of 0.1 seconds. No, Send format T \*\*\* No, Item Subcommand Description Time Set wait time  $(0.1 \sim 99.9)$ 78 GIS: REMOTE TEACH Register Description Wait until the specified general-purpose input port turns ON. Note that the ON hold time should be 10msec or more. For the ON status and input circuit, see "(4) General-purpose I / O". No. Send format G | S | : \*

No,	Item	Subcommand	Description
1	Control	*	Set general-purpose input $(0\sim7)$





No,	Item	Reply	Description
1	Contents	****	Parameter contents



# 7. Parameters

Describes the built-in setting parameters.

When changing parameters, please do after understanding the function fully.

7-1. List of parameters (1)Parameter type

> Perform settings related to stage operation and machine home position return AXIS Parameters:

direction.

GENERAL parameters: Perform settings related to stage control and communication.

Type	No,	Axis	Display	Description				
	01	First	AXIS Name	Setting the display axis name	Page 72			
	02	Second	ANOTAINC	Cotting the display axis harrie	12			
	03	First	UNIT Sel	Unit setting				
	04	Second	01411 001	OTHE GOLDING	72			
	05	First	Pos Dir	  Setting the coordinate count direction	72			
	06	Second	1 00 Bii	Social is a social late social train social i	' _			
	07	First	INPos Range	In-position range setting	72			
	80	Second	11 11 00 1 10.1130	in position range octaining				
	09	First	FB Speed	Setting the feedback speed	73			
	10	Second	. В бросс	Cottill is the recordant operation				
	11	First	ZERO Cont	Zero control setting	73			
	12	Second		25.5 25.14.5. 25.44.75				
	13	First	Stage Config	Stage configuration settings	73			
	14	Second						
	15		TEACH IF	TEACH operation interface settings	73			
	16	First	Acc Cont	Acceleration / deceleration control settings	74			
	17	Second	7 100 00111	7 tooder attorny according to the or octaining				
	18	First	Acc Time	Acceleration / deceleration time setting	74			
	19	Second	7 100 11110	7 toodict attorny according to the control of the c	<u>'</u>			
	20	First	Jog Speed 3	JOG speed 3 setting				
	21	Second	005 00000	oca opeca c setting	74			
	22	First	Jog Speed 2	JOG speed 2 setting	74			
	23	Second	008 00000 2					
	24	First	Jog Speed 1	JOG speed 1 setting	75			
AXIS	25	Second	008 Opcod 1					
	26	First	Jog Cont	Operation control settings when operating CCW and	75			
	27	Second	008 00111	CW buttons				
	28	First	ORG Dir	Setting the machine origin return direction				
	29	Second						
	30	First	ORG Mode Sel	Setting the machine origin return mode				
	31	Second		AA 1				
	32	First	ORG Mode3 Pos	Machine origin return Mode 3 specified				
	33	Second		position setting				
	34	First	ORG Speed H	Machine origin return speed H setting	76			
	35	Second		The state of the s	ļ			
	36	First	ORG Speed M	Machine origin return speed M setting				
	37	Second			76			
	38	First	ORG Speed L	Machine origin return speed L setting	77			
	39	Second		_				
	40	First	EORG Speed	Electric origin return speed setting	77			
	41	Second	•	_				
	42	First	Soft LMT Sel	Software limit function setting	77			
	43	Second						
	44	First	+ Soft LMT Pos	Setting the software limit position on the plus	78			
	45	Second		side				
	46	First	— Soft LMT Pos	Setting the software limit position on the	78			
	47	Second	L ( ) D C=1	minus side	70			
	48		L (-) R Sel	Jog controller left / right button operation axis setting	78			
	49	-	T <-> B Sel	Jog controller up / down button operation axis setting	78			



1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16

Туре	No,	Axis	Display	Description	Page
	50	-	Right Dir	Setting the count direction of the right button of the log controller	79
	51	-	Top Dir	Setting the count direction of the up button of the jog controller	79
	52 53	First Second	Jy Cont	Command "JY:" control settings	79
	54 55	First Second	Jy Speed H	JY speed H setting	79
	56 57	First Second	Jy Speed 9	JY speed 9 setting	80
	58 59	First Second	Jy Speed 8	JY speed 8 setting	80
AXIS	60 61	First Second	Jy Speed 7	JY speed 7 setting	81
	62 63	First Second	Jy Speed 6	JY speed 6 setting	81
	64 65	First Second	Jy Speed 5	JY speed 5 setting	81
	66 67	First Second	Jy Speed 4	JY speed 4 setting	81
	68 69	First Second	Jy Speed 3	JY speed 3 setting	82
	70 71	First Second	Jy Speed 2	JY speed 2 setting	82
	72 73	First Second	Jy Speed 1	JY speed 1 setting	82
	01	-	AXIS Sel	Control target axis setting	82
	02 03	First Second	Max Speed	Maximum speed setting	83
	04 05	First Second	Lin/Rot	Setting the control stage type	83
	06 07	First Second	Stop Sel	Stop control setting	83
	08 09	First Second	EMG Motor Excite	Setting of motor status at emergency stop	83
	10	_	EMG Connector	Enable / disable emergency stop function	84
	11 12	First Second	Motor Excite	Setting the motor status at startup	84
	13 14	First Second	Stage Cont Type	Setting the feedback stage control method	84
	15 16	First Second	Count Sel	Setting display contents of display unit counter	84
051501	17	First Second	CD Drive	Setting the current down drive	84
GENERAL	19 20	First Second	INP Dec	In-position judgment time setting	85
	21	First Second	FBT Sel	Setting the feedback start timing	85
	23		Ini Mode	Initial mode setting	85
	24		I/F Select	Interface settings	85
	25		USB Del	USB delimiter setting	85
	26		GP-IB Addr	GP-IB address setting	86
	27		GP-IB Del	GP-IB delimiter setting	86
	28		GP-IB EOI	GP-IB EOI setting	86
	29 30		GP-IB SRQ	GP-IB SRQ setting	86
	31		ETHER Del	Ethernet IP address setting	86 87
	32		IP Address	Ethernet default gateway setting	87
	32 33		Default Getway Subnet Mask	Ethernet default gateway settings Setting the Ethernet subnet mask	87
	34		ECHO BACK	Setting the Ethernet subhet mask Setting the command echo back function	87
	35		TEACH Monitor	Teaching monitor function setting	87
				General-purpose input port chattering check	
	36	-	GENERAL IN Chat	function setting	88




1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	

Туре	No,	Axis	Display	Description				
	37	-	TEACH IN Chat	TEACH input port chattering check function setting	88			
	38	-	Sleep Sel	Enable / disable sleep function	88			
	39	-	BEEP Sel	Enable / disable beep sound	88			
GENERAL	40	-	Disp bright	Setting the display brightness	88			
	41	First	lmt Stop Sel	Stop the limit sensor	89			
	42	Second	iiii Stop sei	Stop the IIIIIt sensor	09			
	43	-	Option type	Option type selection	88			

## 7-2. Description of display panel

(1) PARAMETER mode top screen



No	ltem		Description				
1	Туре	The type of param	eter.				
2	Cursor	Selection curso					
		Setting	Move to the parameter edit screen.				
		Save	Save the parameters. If the parameter has not been changed, it will not be saved.				
3	Menu	Return	Returns to the mode before entering PARAMETER mode. Coordinate values are maintained. This content is displayed except after changing the GENERAL parameter.				
	iviei iu	Reset To Start	Reboot with the same operation as the command "RESET:". Displayed when "2: Save" is executed with the GENERAL parameter changed.				
		Default	Default all parameters of the displayed type. * After that, if you execute "2: Save", it will be saved with default parameters.				
4	Range	The range of parameter No.					
5	Supported axes	The axis corresponding to the parameter. Parameters without axis display are common to both axes.					
6	Contents	The content of the	parameter.				

<sup>\*</sup> Default parameters may differ from factory parameters. See the attached "Parameter sheet at shipment".


Please do not operating by setting all parameters to default values. Be sure to set according to the performance of the connected stage.




Button	Display	Description
F1 / <b>⇒</b>	Тор	None
F1/ -	Edit	Move the cursor up 10 items.
F2 / <b></b>	Тор	None
FZ / •	Edit	Move the cursor down 10 items.
F3 / 🛧	Тор	Move the cursor up.
F3/ <b>T</b>	Edit	Move the cursor up 1 items.
F4 / <b>♣</b>	Тор	Move the cursor down.
Γ4 / 🔻	Edit	Move the cursor down 1 items.
F5 (Enter)	Тор	Determine the item of the cursor.
13 (Linter)	Edit	Move to the lower layer.
F6 (Cancel)	Тор	None
ro (Caricei)	Edit	Move to the upper layer,
F7	Top	None
1 7	Edit	I VOI IE
F8	Тор	None
1-0	Edit	
F7 + F8	Тор	Move to GENERAL parameter.
17110	Edit	None

## 7-3. Procedure for transition to GENERAL parameters Describes the procedure for moving to the GENERAL parameter.


• Procedure1 For example, when nothing is being executed in the LOCAL mode, pressing both the F7 (PLS / RTN) and F8 (MEU) buttons simultaneously displays the X SETUP mode menu.



• Procedure2 Press F2 button (PRM) Move to PARAMETER mode.



· Procedure3 same time.






F2

F5 Enter

F7

F8





· Procedure4 mode.

# PARAMETER mode [GENERAL]

GENERAL parameter editing [ PARAMETER - Menu] [GENERAL] Setting Save Return Default


#### 7-4. How to change parameter settings

The parameter settings can be changed using the front panel operation or commands. For the setting method, see the setting example below.


#### (1) Example 1

From the LOCAL mode, change the mechanical origin return method "ORG Mode Sel" on axis 1 by operating the front panel.

 Procedure1 For example, when nothing is being executed in the LOCAL mode, pressing both the F7 (PLS / RTN) and F8 (MEU) buttons simultaneously displays the X SETUP mode menu.




· Procedure2 Press F2 button (PRM) Move to PARAMETER mode.




· Procedure3 PARAMETER mode [AXIS]

cursor is at "1: Setting", press the F5 / ENTER button to move to the parameter edit mode.



 Procedure4 PARAMETER mode [AXIS] Setting mode


First axis "ORG ModeSel".



Procedure5

button,

Move the cursor down.





F5


F5 Enter Page

F6 Cance

F7

F8

F2





#### • Procedure11

the cursor to "3: Return". (Note: If you return without saving after making changes, the changes will not be Save reflected.)

PARAMETER mode [AXIS] Press F4 /  $\sqrt[]{}$  button to move [ PARAME TER-Menu] [AXIS]

Setting Return

Default

#### · Procedure12


to return from PARMETER mode to LOCAL mode, This completes the parameter change.

PARAMETER mode [AXIS]









#### (2) Example 2


From the LOCAL mode, change the JOG speed L "Jog Speed 1" on first axis by operating the front panel.

• Procedure  $1 \sim 3$ Same as example 1.

#### · Procedure4

F4 / ♥ button to display first axis "Jog Speed 1".







#### Procedure5

# PARAMETER mode [AXIS] Setting mode

Press the F5 (ENTER) [PARAMETER-Setting] button to move the cursor to the lower layer.

[AXIS 01-73] 1 Jog Speed 1 1.0000] mm/s[ 2 Jog Speed 1  $0\ 0\ 0\ 0\ ]$  m m / s



#### · Procedure6


# PARAMETER mode [AXIS] Setting mode

moves the Under bar from the first digit to the first decimal place.

Pressing the F1  $/ \Rightarrow$  button [ PARAMETER - Setting] [AXIS 01-73] 1 Jog Speed 2 4 ▶ [ <u>1</u> . 0 0 0 0 ] m m / s 2 Jog Speed 1 2 5 1 . 0 0 0 0 ] m m / s








the same as steps 9 to 12 of setting example 1.

After that, the procedure is [PARAMETER-Menu][AXIS]▶ Setting Save Return Default



#### 7-5. How to read parameter descriptions

The following explains how to read the parameter description page.





No	ltem	Contents					
1	Parameter No.	Parameter number. Used with the commands "PRM_GET:" and "PRM_SET:".					
2	Parameter display name	This is the name of the parameter displayed on the display unit.					
3	Parameter name	The name of the parameter.					
4	Parameter description	description of the parameter.					
5	Target axis	The target axis. If the target axis is not displayed, use the common setting for both axes or select the axis to be used.					
6	Choices or Setting value	Indicates the type to be set, Choices or settings,					
7	Details	Indicates a choices candidate or a settable range.					
8	Description	A description of the choice,					
9	Command setting, reply	The value to be sent or returned when using the commands "PRM_SET:" and "PRM_GET:".					
10	Default	Indicates the default of the parameter. * 1					
11	Whether to restart	Indicates whether or not a restart is performed automatically after saving parameters. Black text is executed, gray text is not executed. * 2					
12	Model	Indicates the model name of the feedback stage controller.					

<sup>\* 1</sup> The default parameters may be different from the factory parameters. See the attached "Parameter sheet at shipment".




<sup>\* 2</sup> Restarting is the same operation as the command "RESET:".
















1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16

32 ORG Mode3 Pos Mechanical origin return Mode3 specified position

33

Description Set the return position when mechanical origin return Mode3 is set. Please refer to the "About 9. Home Return for more information

ORG Mode First axis

Setting value m m R G Р 0 Mode О s Second axis Setting value m m

Model	Setting value	Default	Command setting, reply		
FC-114	$0.0001 \sim 999.9999$ mm	0.5000mm	5000		
FC-414	$0.00005 \sim 999.99995$ mm	0.5000mm	50000		
FC-514	$0.00001 \sim 999.99999$ mm	0.5000mm	50000		

o s

34 ORG Speed H Mechanical origin return speed H

35

Description Set the mechanical origin return speed High. Please refer to the "About 9. Home Return" for more information. Set it below the "Max Speed". If the setting value is large, it will return to the previous setting value.

ORG Spe е First axis Setting value m m Ο R G Н р е е d Second axis Setting value m m

Model	Setting value	Default *	Command setting, reply		
FC-114		10,0000mm/s	100000		
FC-414	"Max Speed" mm / s or less	10.0000mm/s	1000000		
FC-514		10.0000mm/s	1000000		

<sup>\*</sup> At the time of parameter default execution, if "Max Speed" is less than the default, it will be "Max Speed".

36 ORG Speed M Mechanical origin return speed M

S

S

ре

Setting value

ре е

500000

37

Description Set the mechanical origin return speed Middle, Please refer to the "About 9. Home Return" for more information. Set it below the "ORG Speed H". If the setting value is large, it will return to the previous setting value.

е

d

First axis

0 RG

0

R

Second axis

S		[ Setting value	] m m / s	
				_
	Model	Setting value	Default *	Command setting, reply
	FC-114		5.000mm/s	50000
	FC-414	"ORG Speed H" mm / s or less	5.0000mm/s	500000

Μ

M

m m



<sup>500000</sup>mm/s FC-514 \* At the time of parameter default execution, if "Max Speed" is less than the default, it will be "Max Speed".

1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
38 ORG Speed I				Mechanical origin return speed I									RESET		

39

Description Set the mechanical origin return speed Low, Please refer to the "About 9. Home Return" for more information, Set it below the "ORG Speed M", If the setting value is large, it will return to the previous setting value.

8 ORG Spe е First axis Setting value m m 9 O R G ре е Second axis Setting value Γ m m

ĺ	Model	Setting value	Default *	Command setting, reply
ſ	FC-114		1.0000mm/s	10000
	FC-414	$^{\prime\prime}$ ORG Speed M $^{\prime\prime}$ mm $/$ s or less	1.00000mm/s	100000
ĺ	FC-514		1.00000mm/s	100000

<sup>\*</sup> At the time of parameter default execution, if "Max Speed" is less than the default, it will be "Max Speed".

40 EORG Speed Electric origin return speed RESET 41

Description Set the electric origin return speed. Please refer to the "About 9. Home Return" for more information. Set it below the "Max Speed". If the setting value is large, it will return to the previous setting value.

First axis

4 0 1 E O R G S p e e d

Setting value ] m m / s

Second axis

4 1 2 E O R G S p e e d

Setting value ] m m / s

Model	Setting value	Default *	Command setting, reply		
FC-114		5.0000mm/s	50000		
FC-414	"Max Speed" mm / s or less	5.0000mm/s	500000		
FC-514		50000mm/s	500000		


<sup>\*</sup> If the "Max Speed" is less than the default of this parameter at the time of parameter default, it will be the value of "Max Speed".

42 Soft LMT Sel Software limit function setting
43 Description Set whether to use the limit function specified by the coordinate value.


S f L M Se 0 t First axis Choices 3 S М 0 е Second axis Choices

	Choices	Description	Command setting, reply
	OFF	Not use	0
	ON	Use	1
Default		0	









Model	Setting value	Default *	Command setting, reply		
FC-114	"Max Speed" mm / s or less	10,0000mm/s	100000		
FC-414		10,0000mm/s	1000000		
FC-514		10.00000mm/s	1000000		

\* At the time of parameter default execution, if "Max Speed" is less than the default, it will be "Max Speed".



1	2	3	4	5	6	7	00	9	10	11	12	13	14	15	16

56 Jy Speed 9 JY Speed 9 57 RESET

Description Set the operation speed 2 when pressing and holding the CW, CCW and Jog controller buttons, Set it below the "Jog Speed 3". If the setting value is large, it will return to the previous setting value.

5 6 J Spee d 9 First axis Setting value m m J S e d 9 р е У Second axis Γ Setting value m m

Model	Setting value	Default *	Command setting, reply		
FC-114		5.0000mm/s	50000		
FC-414	"Jy Speed H" mm $/$ s or less	5.0000mm/s	500000		
FC-514		5.0000mm/s	500000		

<sup>\*</sup> At the time of parameter default execution, if "Max Speed" is less than the default, it will be "Max Speed".

#### ■ About the minimum speed of Jy Speed1 to 8

The minimum operation speed of Jy Speed1 to 8 is determined by the setting of Jy Speed9. See the table below for details. For example, if Jy Speed 9 is set to 10 mm / s while using FC-114, the minimum operation speed is 0.0002 mm / s. The minimum speed of FC-114 can be set to 0.0001 mm / s, but in this example, it does not operate because the minimum speed is 0.0002 mm / s.

Model	Jy Speed 9 s	Minimum speed		
	0.0001mm/s	_	6.5535mm/s	0.0001mm/s
	6,5536mm/s	-	13.1070mm/s	0.0002mm/s
FC-114	13.1071mm/s	-	32.7675mm/s	0.0005mm/s
	32.7676mm/s	_	65,5350mm/s	0.0010mm/s
	65,5351mm/s	-	100,000mm/s	0.0020mm/s
	0.0005mm/s	-	3,27675mm/s	0.0005mm/s
	3,27680mm/s	-	6.55350mm/s	0.00010mm/s
FO 44.4	6,55355mm/s	-	16,38375mm/s	0.00025mm/s
FC-414	16,38380mm/s	-	32.76750mm/s	0.00050mm/s
	32,76755mm/s	-	65,53500mm/s	0.00100mm/s
	65,53550mm/s	-	100,000mm/s	0.00250mm/s
	0.00001mm/s	-	0,65535mm/s	0.00001mm/s
	0,65536mm/s	-	1.31070mm/s	0.0002mm/s
	1.31071mm/s	-	3,27675mm/s	0.0005mm/s
FC-514	3,27676mm/s	_	6,55350mm/s	0.00010mm/s
	6.55351mm/s	-	13.10700mm/s	0.00020mm/s
	13.10701mm/s	-	32.76750mm/s	0.00050mm/s
	32.76751mm/s	-	50,0000mm/s	0.00100mm/s

58 Jy Speed 8 JY Speed 8 59

Description Set the 8th operation speed of JY speed. Set it below the "Jy Speed 9". If the setting value is large, it will return to the previous setting value.

8 8 J S У ре е First axis Setting value m m 8 J e e d У р Second axis Setting value m m

Model	Setting value	Default *	Command setting, reply
FC-114		1.0000mm/s	10000
FC-414	"Jy Speed 9" mm / s or less	1,00000mm/s	100000
FC-514		1.00000mm/s	100000


<sup>\*</sup> At the time of parameter default execution, if "Max Speed" is less than the default, it will be "Max Speed".




14 60 Jy Speed 7 JY Speed 7 61 Description Set the 7th operation speed of JY speed. Set it below the "Jy Speed 8". If the setting value is large, it will return to the previous setting value ре е First axis Setting value m mр е е Second axis Setting value m m Model Setting value Default \* Command setting, reply FC-114 0.5000mm/s 5000 FC-414 0,50000mm/s 50000 "Jy Speed 8" mm / s or less 0,50000mm/s 50000 FC-514 \* At the time of parameter default execution, if "Max Speed" is less than the default, it will be "Max Speed". 62 Jy Speed 6 JY Speed 6 63 Description Set the 6th operation speed of JY speed. Set it below the "Jy Speed 7". If the setting value is large, it will return to the previous setting value. 2 J Speed 6 У First axis Setting value ] m m / s 6 3 J S e e d 6 У Second axis Setting value ] m m / s Model Setting value Default \* Command setting, reply FC-114 0.1000mm/s 1000 FC-414 "Jy Speed 7" mm / s or less 0.10000mm/s 10000 FC-514 0.10000mm/s 10000 \* At the time of parameter default execution, if "Max Speed" is less than the default, it will be "Max Speed". 64 Jy Speed 5 JY Speed 5 65 Description Set the 5th operation speed of JY speed. Set it below the "Jy Speed 6". If the setting value is large, it will return to the previous setting value. 4 S 6 е е р First axis Setting value m m S 5 5 6 J е У е Second axis Setting value m m Command setting, reply Model Setting value Default \* 500 FC-114 0.0500 mm/sFC-414 5000 "Jy Speed 6" mm / s or less 0.05000mm/s 5000 0.05000mm/s FC-514 \* At the time of parameter default execution, if "Max Speed" is less than the default, it will be "Max Speed". 66 JY Speed 4 Jy Speed 4 67 Description Set the 4th operation speed of JY speed. Set it below the "Jy Speed 5". If the setting value is large, it will return to the previous setting value. 6 J р е У First axis Setting value m m J е р е У Second axis Setting value m m / Default \* Model Setting value Command setting, reply FC-114 0.0100mm/s 100 FC-414 "Jy Speed 5" mm / s or less 0.01000mm/s 1000 0.01000mm/s FC-514 1000

\* At the time of parameter default execution, if "Max Speed" is less than the default, it will be "Max Speed".






























	2   3	4   5   6	8 9 10 11 12 13	3   14   15   16
36	GENER	RAL IN Chat Ger	neral-purpose input port chattering check	RESET
	Description	Set whether to ch	eck chattering of the general-purpose input port.	
		3 6 G	ENERAL IN Chat	
		[ ]	Choices ]	
	'			<del></del>
		Choices	 Description	Command setting, reply
		OFF	Disable	0
		ON	Enable	1
	Default		OFF	0
37	TEACI	LINI Ola a t	TEACL Lines to control attacks and also also	DECET
31	TEACE	H IN Chat	TEACH input port chattering check	RESET
	D	O - 4 do 4do 4 do		
	Description I		eck chattering of the teaching operation input port.	<del></del>
		3 7 T	EACH IN Chat	
	l	L	Choices ]	
	ı	01		0
		Choices	Description	Command setting, reply
		OFF	Disable	0
	5 ( )	ON	Enable	1
	Default		OFF	0
38	Sleep (	Sel Sle	ер	RESET
	Description	Set whether to us	e the sleep function.	
		3 8 S	leep Sel	
		[ ]	Choices ]	
	•			
		Choices	Description	Command setting, reply
		Choices OFF	Description  Not use	Command setting, reply O
			·	
	Default	OFF	Not use	0
	Default	OFF	Not use Use	O 1
39	l	OFF ON	Not use Use ON	0 1 1
39	Default BEEP	OFF ON	Not use Use ON	O 1
39	BEEP	OFF ON Sel Bee	Not use Use ON	0 1 1 RESET
39	BEEP	OFF ON  Sel Bee	Not use Use ON ep nit a beep sound. However, continuous sound cannot	0 1 1 RESET
39	BEEP	OFF ON  Sel Bee	Not use Use ON  ep nit a beep sound. However, continuous sound cannot E E P S e 1	0 1 1 RESET
39	BEEP	OFF ON  Sel Bee	Not use Use ON ep	0 1 1 RESET
39	BEEP	OFF ON  Sel Bee	Not use Use ON  Ep  nit a beep sound, However, continuous sound cannot  E E P S e 1 Choices  ]	O 1 1 RESET
39	BEEP	OFF ON  Sel Bee  Set whether to en  3 9 B [	Not use Use ON  Paper  Description	0 1 1 RESET be turned off.  Command setting, reply
39	BEEP	OFF ON  Sel Bee  Set whether to en  3 9 B [ Choices OFF	Not use Use ON  Paper Sep Sep Sep Sep Sep Sep Sep Sep Sep Sep	O 1 1 RESET be turned off.  Command setting, reply 0
39	BEEP :	OFF ON  Sel Bee  Set whether to en  3 9 B [	Not use Use ON  Paper  Description	0 1 1 RESET be turned off.  Command setting, reply
39	BEEP	OFF ON  Sel Bee  Set whether to en  3 9 B [ Choices OFF	Not use Use ON  Pep  mit a beep sound. However, continuous sound cannot  E E P S e l Choices ]  Description Do not beep Sound a beep	O 1 1 RESET be turned off.  Command setting, reply 0 1
	BEEP: Description  Default	OFF ON  Sel Bee  Set whether to en 3 9 B [ Choices OFF ON	Not use Use ON  The property of the property o	O 1 1 RESET be turned off.  Command setting, reply O 1 1
39	BEEP :	OFF ON  Sel Bee  Set whether to en 3 9 B [ Choices OFF ON	Not use Use ON  Pep  mit a beep sound. However, continuous sound cannot  E E P S e l Choices ]  Description Do not beep Sound a beep	O 1 1 RESET be turned off.  Command setting, reply 0 1
	BEEP: Description  Default	OFF ON  Sel Bee  Set whether to en  3 9 B [ Choices OFF ON  right Dis	Not use Use ON  Pep  mit a beep sound. However, continuous sound cannot  E E P S e 1 Choices  Description Do not beep Sound a beep ON  play brightness	O 1 1 RESET be turned off.  Command setting, reply O 1 1
	BEEP: Description  Default	OFF ON  Sel Bee  Set whether to en  3 9 B [ Choices OFF ON  right Dis  Set the display bri	Not use Use ON  ON  Pep  nit a beep sound, However, continuous sound cannot  E E P S e 1 Choices  Description Do not beep Sound a beep ON  play brightness ghtness.	O 1 1 RESET be turned off.  Command setting, reply O 1 1
	BEEP: Description  Default	OFF ON  Sel Bee  Set whether to en  3 9 B [ Choices OFF ON  right Dis	Not use Use ON  Pep  It a beep sound, However, continuous sound cannot E E P S e l Choices  Description Do not beep Sound a beep ON  play brightness ghtness. i s p b r i g h t	O 1 1 RESET be turned off.  Command setting, reply O 1 1
	BEEP: Description  Default	OFF ON  Sel Bee  Set whether to en  3 9 B [ Choices OFF ON  right Dis  Set the display bri	Not use Use ON  ON  Pep  nit a beep sound, However, continuous sound cannot  E E P S e 1 Choices  Description Do not beep Sound a beep ON  play brightness ghtness.	O 1 1 RESET be turned off.  Command setting, reply O 1 1
	BEEP: Description  Default	OFF ON  Sel Bee Set whether to en 3 9 B [ Choices OFF ON  right Dis Set the display bri 4 0 D [	Not use Use ON  Paper ON  Paper ON  Paper ON  Paper ON  Description Do not beep Sound a beep ON  Paper ON	O 1 1 RESET be turned off.  Command setting, reply O 1 1 RESET
	BEEP: Description  Default	OFF ON  Sel Bee  Set whether to en  3 9 B [ Choices OFF ON  right Dis  Set the display bri.  4 0 D [ Choices	Not use Use ON  Paper ON  Paper ON  Paper ON  Paper ON  Description Do not beep Sound a beep ON  Paper ON  Paper ON  Paper ON  Description Do not beep Sound a beep ON  Description Do not beep Sound a beep ON  Description Do not beep Sound a beep ON  Description  Description  Description	O 1 1 RESET be turned off.  Command setting, reply O 1 1 RESET  RESET
	BEEP: Description  Default	OFF ON  Sel Bee Set whether to en 3 9 B [ Choices OFF ON  right Dis Set the display bri 4 0 D [ Choices 25%	Not use Use ON  Pep  Init a beep sound. However, continuous sound cannot  E E P S e l Choices ]  Description Do not beep Sound a beep ON  play brightness ghtness. i s p b r i g h t Choices ]  Description Set to 25%	O 1 1 RESET be turned off.  Command setting, reply O 1 1 RESET  Command setting, reply O Command setting, reply O Command setting, reply O
	BEEP: Description  Default	OFF ON  Sel Bee Set whether to en 3 9 B [ Choices OFF ON  right Dis Set the display bri 4 0 D [ Choices 25% 50%	Not use Use ON  ON  Deposit a beep sound, However, continuous sound cannot  E E P S e l Choices ]  Description Do not beep Sound a beep ON  play brightness ghtness. i s p b r i g h t Choices ]  Description Set to 25% Set to 50%	O 1 1 RESET be turned off.  Command setting, reply O 1 RESET  Command setting, reply O 1 1 Command setting, reply O 1 1
	BEEP: Description  Default	OFF ON  Sel Bee  Set whether to en 3 9 B [  Choices OFF ON  right Dis  Set the display bri. 4 0 D [  Choices 25% 50% 75%	Not use Use ON  ON  Pep  It a beep sound, However, continuous sound cannot  E E P S e l Choices ]  Description Do not beep Sound a beep ON  play brightness ghtness, i s p b r i g h t Choices ]  Description Set to 25% Set to 50% Set to 75%	Command setting, reply  RESET  De turned off.  Command setting, reply  O  1  1  Command setting, reply  O  1  1  Command setting, reply  O  1  2
	BEEP: Description  Default	OFF ON  Sel Bee  Set whether to en  3 9 B [  Choices  OFF ON  right Dis  Set the display bri.  4 0 D [  Choices  25% 50% 75% 100%	Not use Use ON  ON  Pep  It a beep sound. However, continuous sound cannot  E E P S e l Choices  Description Do not beep Sound a beep ON  play brightness ghtness. i s p b r i g h t Choices  Description Set to 25% Set to 50% Set to 75% Set to 100%	Command setting, reply  O  1  RESET  be turned off.  Command setting, reply  O  1  1  Command setting, reply  O  1  2  3
	BEEP: Description  Default	OFF ON  Sel Bee  Set whether to en  3 9 B [  Choices OFF ON  right Dis  Set the display bri  4 0 D [  Choices 25% 50% 75% 100% 125%	Not use Use ON  ON  Deposit a beep sound. However, continuous sound cannot  E E P S e l Choices ]  Description Do not beep Sound a beep ON  play brightness ghtness. i s p b r i g h t Choices ]  Description Set to 25% Set to 50% Set to 75% Set to 100% Set to 125%	Command setting, reply  O 1 1 RESET  be turned off.  Command setting, reply  O 1 1 1  Command setting, reply  O 1 2 3 4
	BEEP: Description  Default	OFF ON  Sel Bee  Set whether to en  3 9 B [  Choices OFF ON  right Dis  Set the display bri  4 0 D [  Choices 25% 50% 75% 100% 125% 150%	Not use Use ON  Paper Description Description Do not beep Sound a beep ON  Play brightness  ghtness.  i s p b r i g h t Choices  Description Set to 25% Set to 50% Set to 75% Set to 100% Set to 125% Set to 150% Set to 150%	RESET  be turned off.  Command setting, reply  0 1 1 1  RESET  Command setting, reply  0 1 2 3 4 5
	BEEP: Description  Default	OFF ON  Sel Bee  Set whether to en  3 9 B [  Choices OFF ON  right Dis  Set the display bri  4 0 D [  Choices 25% 50% 75% 100% 125% 150% 175%	Not use Use ON  ON  Deposit a beep sound. However, continuous sound cannot  E E P S e l Choices ]  Description Do not beep Sound a beep ON  play brightness ghtness. i s p b r i g h t Choices ]  Description Set to 25% Set to 50% Set to 75% Set to 100% Set to 125% Set to 150% Set to 175% Set to 175%	RESET  be turned off.  Command setting, reply  0 1 1 1  RESET  Command setting, reply  0 1 2 3 4 5 6
	BEEP: Description  Default	OFF ON  Sel Bee  Set whether to en  3 9 B [  Choices OFF ON  right Dis  Set the display bri  4 0 D [  Choices 25% 50% 75% 100% 125% 150%	Not use Use ON  Paper Description Description Do not beep Sound a beep ON  Play brightness  ghtness.  i s p b r i g h t Choices  Description Set to 25% Set to 50% Set to 75% Set to 100% Set to 125% Set to 150% Set to 150%	RESET  be turned off.  Command setting, reply  0 1 1 1  RESET  Command setting, reply  0 1 2 3 4 5





Option type selection

Optiom

43

Option type

Default

4 3

Description Select the option type to connect.

 Choices
 Description
 Command setting, reply

 Type0
 JC-01, JC-01-3, JC-01-04
 0

 Type1
 MD-400
 1

TypeO

t у р е



RESET

0

1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16

### 8. Teaching function

The teaching function is a function to execute commands registered inside this equipment when in the TEACH mode. It can execute all lines automatically or one line at a time. The registered commands are executed sequentially from the first line, and execution ends when the command "END:" or the completion of the 200th line. If there is a line where no command is registered (blank line), skip that line. In the case of an error that can be cleared, executing the error clear returns to the first line. The error can be cleared by operating the "BEC" button of the target axis after switching to the command "BEC:" or LOCAL mode.

 $\triangle$ 

- · Before moving the stage, make sure that there is no effect on the surroundings.
- · Check the registration details before performing teaching.

#### 8-1. About channels

There are 1 to 5 channels, and the number of stored lines per channel is 200. The following table shows the image stored in the internal memory.

Line	Channel1	Channel2	Channel3	Channel4	Channel5
1	D:1M1	D:1M1	JG:1 +100	JG:1 -100	FS:1 100
2	M:1+M10	M:1-M10			FS:2 100
3	G	G	END:	END:	M:1+U100
4	END:	END:			G
5					FE:2
			}		
198					FE:1
199					END:
200					

#### 8-2. About registration and editing

Registration of contents is possible only by command. The commands used for registration and editing are shown below.

Item	No,	Commands	Function	Page
Teaching	56	T_ON:	Move to teaching edit mode	56
	57	T_OFF:	Save teaching content and return to TEACH mode	56
	58	T_DEL:	Delete the contents of the line registered for teaching	56
	59	T_SET:	Set teaching content for each line	56
	60	T_GET:	Get the contents of the line registered for teaching	56
	61	TC:	Select teaching channel	56
	62	TCR:	Get the current teaching channel	57
	63	TQ:	Get teaching status	57
	68	TR:	Get teaching channel registration status	58

#### (1) Example 1

Set the operation speed of the first axis to 1 mm / sec. After moving 10 times in increments of 1 mm from the current position, move the axis in the negative direction 10 mm and wait for 1 second. Repeat this 10 times.

Send command	Description
T_ON:	Move to TEACH edit mode
TC:1	Select channel 1
T_SET:1 D:1M1	First axis operation speed set to 1mm / sec
T_SET:2 H:1	Set First axis mechanical home return
T_SET:3 FS:1 10	Set 10 times for loop level 1
T_SET:4 M:1+M1	Set first axis relative movement set value to 1mm
T_SET:5 FS:2 10	Set 10 times for loop level 2
T_SET:6 GN:1	Start moving
T_SET:7 FE:2	End location of loop level 2
T_SET:8 M:1-M10	Set first axis relative movement set value to -10mm
T_SET:9 G	Start moving
T_SET:10 T:1.0	Wait 1 second
T_SET:11 FE:1	End location of loop level 1
T_SET:12 END:	Teaching end line setting
T_OFF:	Save teaching contents and return to normal mode



## (2) Example 2

Retrieve the first and fourth lines from Example 1.

Send command	Reply command	Description
T_ON:	-	Move to TEACH edit mode
TC:1	-	Select channel 1
T_GET:1	D:1M1	Get first line
T_GET:4	M:1+M1	Get line 4
T_GET:13	-	Not registered

#### (3) Example 3

Delete the first to fourth lines of Example 1.

Send command	Description
T_ON:	Move to TEACH edit mode
TC:1	Select channel 1
T_DEL:1	Delete first line
T_DEL:2	Delete line 2
T_DEL:3	Delete line 3
T_DEL:4	Delete line 4
T_OFF:	Save teaching contents and return to normal mode

# (4) Registration target command

This is a list of commands that can be registered for teaching.

Item	No,	Commands	Function	Page
Control	09	F:	Change the positioning control method(Closed or Open loop)	28
	11	C:	Change the excitation state of the motor	28
Origin	33	H:	Return to mechanical origin	44
	34	Z:	Return to electrical origin	44
	35	R:	Set the position coordinate values to zero	45
Motion	38	ACC:	Set the acceleration and deceleration time	46
	40	D:	Set the movement speed of the stage	47
	42	A:	Set the position to move in absolute motion	49
	43	M:	Set the amount of movement to move in relative motion	50
	44	G	Start moving (The setting values is lost after execution)	51
	45	GN:	Start moving (The setting values is not lost.)	51
	48	JG:	It moves according to the set number of pulses	53
General purpose I/O	55	O:	Set general-purpose Output status	55
Teaching	74	FS:	Loop setting	61
registration	75	FE:	Set loop end	61
only	76	END:	Set the execution end line	61
	77	T:	Set the wait time	61
	78	GIS:	Wait until the specified general-purpose input state is reached	61



																4
1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	ĺ


#### 8-3. Teaching operation

Perform teaching operations with the front panel, commands, general-purpose I / O, and jog controller. When "I / O" is selected in the parameter "TEACH IF", operation is possible only with general-purpose I / O. When "JOG / CMD" is selected, operation can be performed from other than general-purpose I / O. There is no need to unify the types of operations. For example, you can start execution on the front panel and stop it with a jog controller or command.

A Before moving the stage, make sure that there is no effect on the surroundings.

#### (1) Front panel operation

Operable when the parameter "TEACH IF" is "JOG / CMD".



No	ltem	Contents
1	Line number	Displays the current line number.
2	Mode	T (TEACH) is displayed.
3	Confirmation of registration	You can check the registration of teaching content.

Button	Display	Contents						
F1	C1 <b>I</b>	Select the teaching channel, It changes each time the F1 (C) button is pressed. If no command is registered, the right side of the channel number will be blank. When the display is blinking, the teaching contents are being registered, so the F1 to F8 buttons cannot be operated. The state on the left shows that channel 1 is selected and commands are registered.  C1 → C2 → C3 → C4 → C5 → Return to C1						
F2	GO	Perform teaching. The display flashes during execution.						
F3	STP	Stop the teaching execution, The blinking "GO", "PAU", and "ONE" will stop						
F4	SIP	blinking.						
F5	_	If the teaching execution has stopped, press the switch twice (double-click) to switch to the operation test mode.  * The operation test mode has nothing to do with teaching.						
F6	PAU	If pressed during teaching, pauses after the current execution is completed. The display flashes during pause. To restart, press F2 (GO) button.						
F7		Press when stopped or paused to execute one command line. The display						
F8 ONE		flashes during execution of one line. After execution of one line, the display stops blinking and pauses.						
F7 + F8	_	If both button are pressed at the same time switches to the SETUP mode, For details on the SETUP mode, see " $5-3$ , SETUP mode".						

Channel selection
 Select a channel with the F1 (C) button. It changes each time it is pressed.




1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16

#### · Start moving


Press F2 (GO) button to start execution, GO starts flashing when execution starts, Execution starts from the first line, and when the command "END:" or the 200th line is completed, execution and GO stop blinking, return to the state before execution, and wait. If an error occurs or an emergency stop occurs during execution, the operating stage stops immediately and stops at the line where it was executing. If a stop occurs due to a limit error, overflow error, or emergency stop, the line returns to the first line when released with the BEC button or command "BEC:" on the front panel. To clear any error other than the above, restart the power, restart, or send the command "RESET:".

When you press the F2 (GO) button, It will be executed. \* GO blinks




#### Pause

To pause during execution, press the F6 (PAU) button. Pauses after execution of the currently executed line is completed, and PAU blinks. To resume execution, press the F2 (GO) button.



#### One line execution


Press F7 (ONE) button or F8 (ONE) button to execute one line at a time. Can be executed while paused or stopped. When stopped, press F7 (ONE) button or F8 (ONE) button to execute one line and pause. When paused, one line is executed each time the button is pressed. However, if the stage is operating, no operation will be accepted. If an error occurs or an emergency stop occurs during execution, the operating stage stops immediately and stops at the line where it was executing. In the case of a stop due to a limit error, overflow error, or emergency stop, if it is released with the BEC button on the front panel or the command "BEC:", it returns to the first line and waits. To clear any error other than the above, restart the power, restart, or send the command "RESET:".








Press F2 (GO) button to execute from the next line. \* GO blinks.



003

#### Stop

If you press the F3 (STP) button or the F4 (STP) button, the stage will stop immediately if it is running, and will stop executing. After stopping, it returns to the first line.



#### (2) Command operation

Operable when the parameter  $\rm ^{''}TEACH~IF''$  is  $\rm ^{''}JOG~/~CMD''$  . The list of operation commands is shown below.

Item	No,	Commands	Description	Page
Teaching	61	TC:	Select a channel	56
	62	TCR:	Get current channel	57
	63	TQ:	Get teaching status	57
	64	TG:	Teaching start	57
	65	TP:	Pause	57
	66	TO:	Execute line by line	57
	67	TL:	Stop execution	58
	68	TR:	Get channel subscription status	58
	69	TFR:	Get the loop count	58
	70	TM:	Set the Teaching monitor function	59
	71	TMR:	Get the Teaching monitor function setting	60
	72	TNR:	Get current line number	60
	73	TACR:	Get current execution command	60

#### (3) Jog controller operation

Operable when the parameter  $\rm ^{''}TEACH~IF''$  is  $\rm ^{''}JOG~/~CMD''$  . See the jog controller instruction manual.



																4
1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	

(4) Operation with general-purpose I / O Operable when the parameter "TEACH IF" is "I / O". Teaching operation can be performed by switch input, For details, see "(4) General-purpose I / O" .

The input current waveform should be pulsed (rise and fall time <100  $\mu$ sec).

Terminal name	Terminal number	Description				
T_1	43					
T_2	19	Change channel				
T_3	44					
T_START	20	Start moving				
T_PAUSE	45	Pause				
T_ONE_STEP	21	One line execution				
T_STOP	46	Stop				

#### · Explanation of terms

l	ltem	Description
ſ	ON	ON means that a current flows through the input terminal of the photocoupler inside this
		equipment.
ſ	OFF	OFF means to cut off the current flowing to the input terminal of the photocoupler inside this
		equipment.

#### Explanation of Terminal

Item		Description									
T_1、T_2、T_3	T_STOP are av T_2, and T_3 a	$\Gamma_1$ , $\Gamma_2$ , and $\Gamma_3$ are pins for selecting channels. $\Gamma_2$ Tand $\Gamma_3$ and $\Gamma_4$ and $\Gamma_5$ and $\Gamma_6$ are available while you continue to select channels at $\Gamma_1$ , $\Gamma_2$ , and $\Gamma_3$ . If $\Gamma_1$ , $\Gamma_2$ , and $\Gamma_3$ are set to other than channels 1 to 5, the above terminal functions cannot be used. If you change the state of $\Gamma_1$ , $\Gamma_2$ , $\Gamma_3$ during execution, stop the execution.									
	Terminal name										
	T_1 T 2	ON OFF	OFF ON	ON ON	OFF OFF	ON OFF					
	T 3	OFF	OFF	OFF	ON	ON					
T_START	T_START is a terminal to start teaching execution. Turn ON for 10ms or more with pulse width.										
T_PAUSE	T_PAUSE is a to the next line.	erminal for tem	porarily stoppina	g teaching. Whil	e ON, pause wit	hout executing					
T_ONE_STEP	T_ONE_STEP is a terminal to execute the registered contents line by line. Executes each time a pulse is input. However, if the stage is operating, this command will not be permitted and will be discarded until positioning is completed. If a pulse is input to the T_ONE_STEP pin while the T_PAUSE pin is turned ON when T_START is not being executed, the registered contents can be executed line by line from the first line. The pulse width should be 10ms or more.										
T_STOP	T_STOP is a p stopping the te 10ms or more v	eaching execut	ion. After stop	-	_	operating and e. Turn ON for					



#### · Check execution status


The status can be checked with the following terminals.

Terminal name	Terminal number	Description
T_STATUS	17	Output HIGH during execution of registration line
/T_STATUS	42	Inverted output of T_STATUS

<sup>\*1</sup> The T\_STATUS signal is output regardless of the setting of the parameter "TEACH IF".

#### · Operation example and execution status of each terminal

Line number	Registered contents	Description
1	D:1M1	First axis operation speed set to 1mm / sec
2	H:1	Set First axis mechanical home return
3	M:1+M1	Set first axis relative movement set value to +1mm
4	G	Start moving
⑤	END:	Teaching end line setting



- \*2 ON of the T\_STATUS signal has the same meaning as T\_STATUS pin HIGH.
- \*3 The OFF time of the T\_STATUS signal during automatic execution is about 500  $\mu$  sec. If the OFF time of the T\_STATUS signal continues for 500  $\mu$  sec or more (for example, 1 msec) while the T\_PAUSE pin is not turned ON, determine that automatic execution has ended.
- \*4 After executing the last line in one line execution, it cannot be re-executed unless the T\_PAUSE pin is turned OFF for 10msec or more.
- \*5 If the T\_PAUSE pin remains ON after executing the last row in one row execution, it cannot be re-executed unless the T\_STOP pin is turned ON for at least 10 msec.



1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16

#### 9. Home return

There are two types of origin, mechanical origin and electric origin.

 $\Lambda$ 

Before moving the stage, make sure that there is no effect on the surroundings.

Туре	Contents
Machine origin	Position after Mode 0 to 4 operation based on the limit sensor mounted inside the stage
Electric origin	Zero position of coordinate value

#### 9-1. Mechanical origin return

#### (1) Mechanical origin return ModeO

Move to the position of 0.5mm from the limit sensor and set the coordinate value to zero.


i. Operates from any position to the limit sensor detection position.

Speed: Parameter "ORG Speed H" Direction: Parameter "ORG Dir"

ii. It operates until it passes through the limit sensor. Speed: 1/10 of parameter "ORG Speed L" Direction: Opposite direction to parameter "ORG Dir"

iii. It moves 0.5mm from the position where it passed the limit sensor. Speed: Parameter "ORG Speed L"

Direction: Opposite direction to parameter "ORG Dir"



#### (2) Mechanical origin return Mode1


Move to the position of 1mm from the limit sensor and set the coordinate value to zero.

i. Operates from any position to the limit sensor detection position.

Speed: Parameter "ORG Speed H" Direction: Parameter "ORG Dir"

ii. It operates until it passes through the limit sensor. Speed: 1/10 of parameter "ORG Speed L" Direction: Opposite direction to parameter "ORG Dir"

iii. Operates 1mm from the position where the limit sensor has passed.
 Speed: Parameter "ORG Speed L"
 Direction: Opposite direction to parameter "ORG Dir"



#### (3) Mechanical origin return Mode2

Move to the middle point between the CW and CCW limit sensors and set the coordinate value to zero.

i. Operates from any position to the limit sensor detection position.

Speed: Parameter "ORG Speed H" Direction: Parameter "ORG Dir"

ii. It operates until it passes through the limit sensor.
 Speed: 1/10 of parameter "ORG Speed L"
 Direction: Opposite direction to parameter "ORG Dir

iii. Operates to the limit sensor detection position on the opposite side,

Speed: Parameter "ORG Speed H"

Direction: Opposite direction to parameter "ORG Dir"


iv. It operates until it passes through the limit sensor.

Speed: 1/10 of parameter "ORG Speed L"

Direction: Parameter "ORG Dir"

v. Operates to the middle point between limit sensors.

Speed: Parameter "ORG Speed L' Direction: Parameter "ORG Dir"





#### (4) Mechanical origin return Mode3

Move from the limit sensor to the setting position of the parameter "ORG Mode3 Pos" and set the coordinate value to zero.

i. Operates from any position to the limit sensor detection position.


Speed: Parameter "ORG Speed H" Direction: Parameter "ORG Dir

ii. It operates until it passes through the limit sensor. Speed: 1/10 of parameter "ORG Speed L"

Direction: Opposite direction to parameter "ORG Dir"

iii. It operates from the position where the limit sensor has been passed to the position set in the parameter "ORG Mode3 Pos".

Speed: Parameter "ORG Speed M" Direction: Opposite direction to parameter "ORG Dir"



#### (5) Mechanical origin return Mode4

<u> /!\</u>

Do not use it on any stage other than the stage equipped with the ORG sensor. Unintended behavior.

Operates to the ORG sensor detection position and sets the coordinate value to zero.

i. Operates from any position to the limit sensor detection position.

Speed: Parameter "ORG Speed H" Direction: Parameter "ORG Dir

ii. It operates until it passes through the limit sensor. Speed: 1/10 of parameter "ORG Speed L

Direction: Opposite direction to parameter "ORG Dir

iii. Operates to the ORG sensor detection position.


Speed: Parameter "ORG Speed M"

Direction: Opposite direction to parameter "ORG Dir"

iv. It operates until it comes out of the ORG sensor.

Speed: Parameter "ORG Speed L"

Direction: Parameter "ORG Dir



#### 9-2. Electric home return

#### (1) Electric home return

Move to the position of zero (Omm) of the coordinate value from the current position,

Speed: Parameter "EORG Speed" Direction: Zero direction of coordinate value from current position

Current position Zero (Omm) position



1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16

#### 10. About Movement Test

You can check that the stage is connected correctly to this instrument. There are three types of operation, and operations are performed on the front panel. See "5-7. MOVEMENT TEST mode" for how to shift to the operation test mode.

 $\overline{\mathbb{A}}$ 

Before moving the stage, make sure that there is no effect on the surroundings.

#### 10-1. Movement Test

#### (1) Step Movement Test SMT (F1 button)

The travel distance obtained by dividing the stage stroke (nominal value) into 10 steps is one step, and the stage is reciprocated in 10 steps. After reciprocating, the stroke calculation operation is performed again. The stage stroke (nominal value) is calculated from the value by operating between the CW and CCW limit sensors.

i. Operates from any position to the limit sensor detection position.

Speed: Parameter "ORG Speed H" Direction: Parameter "ORG Dir"

ii. It operates until it passes through the limit sensor.Speed: 1/10 of parameter "ORG Speed L"Direction: Opposite direction to parameter "ORG Dir

iii. Operates to the limit sensor detection position on the opposite side.

Speed: Parameter "ORG Speed H"

Direction: Opposite direction to parameter "ORG Dir"

iv. It operates until it passes through the limit sensor. Speed: 1/10 of parameter "ORG Speed L"

Direction: Parameter "ORG Dir"

v. Operates 1/10 stroke steps. Speed: JOG speed 1-3

Direction: Parameter "ORG Dir"

vi. Operates 1/10 stroke steps. Speed: JOG speed 1-3

Direction: Opposite direction to parameter "ORG Dir"

vii. Return to i.

# Any position i Limit sensor ii vi

# (2) Movement test between limit sensors LMT (F2 button)

Reciprocates between the CW and CCW limit sensor detection positions.

i. Operates from any position to the limit sensor detection position.
 Speed: Parameter "ORG Speed H"

Direction: Opposite direction to parameter "ORG Dir"

ii. It operates until it passes through the limit sensor. Speed: 1/10 of parameter "ORG Speed L"

Direction: Parameter "ORG Dir

iii. Operates to the limit sensor detection position on the opposite side.

Speed: Parameter "ORG Speed H" Direction: Parameter "ORG Dir"

iv. It operates until it passes through the limit sensor.

Speed: Parameter "ORG Speed H"


Direction: Opposite direction to parameter "ORG Dir"

v. Operates to the limit sensor detection position on

the opposite side. Speed: Parameter "ORG Speed H"

Direction: Parameter "ORG Dir"

vi. Return to iv.





1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16

### (3) Milli Step Movement Test MMT (F5 button)

Calculate the stage stroke (nominal value) and reciprocate between strokes in 1mm steps. After reciprocating, the stroke calculation operation is performed again. The stage stroke (nominal value) is calculated from the value by operating between the CW and CCW limit sensors.

i. Operates from any position to the limit sensor detection position.

Speed: Parameter "ORG Speed H" Direction: Parameter "ORG Dir"

ii. It operates until it passes through the limit sensor. Speed: 1/10 of parameter "ORG Speed L" Direction: Opposite direction to parameter "ORG Dir"

iii. Operates to the limit sensor detection position on the opposite side.

Speed: Parameter "ORG Speed H"

Direction: Opposite direction to parameter "ORG Dir"

iv. It operates until it passes through the limit sensor. Speed: 1/10 of parameter "ORG Speed L" Direction: Parameter "ORG Dir"

v. Operates 1mm step.


Speed: JOG Speed 1 ~ 3 Direction: Parameter "ORG Dir"

vi. Operates 1mm step.

Speed: JOG Speed  $1 \sim 3$ 

Direction: Opposite direction to parameter "ORG Dir"

vii. Return to i.





### 11. About status

Status can be checked with the display and status command. For the display, see "5. Operation" . Refer to "6. Command" for the command.

#### 11-1. Operation related status

This status indicates the operation status of the connection stage and this equipment,

#### (1) Operation status list

No	Contents	Reply *
01	Normal stop (This state is READY)	K
02	During command move (This state is BUSY)	M
03	Out of the in-position range (After positioning is completed) (This state is FREADY)	F
04	Out of the in-position range (During fine adjustment) (This state is FBUSY)	G
05	During electrical origin return	Р
06	During mechanical origin return	0
07	CW side limit stop	С
80	CCW side limit stop	W
09	CW side software limit stop	В
10	CCW side software limit stop	V
11	CW side slowdown sensor area	А
12	CCW side slowdown sensor area	U
13	Error occurred	E
14	Motor is transitioning to excitation	Н
15	Motor is transitioning to non-excitation	
16	Disabled axis (Not set by parameter "AXIS Sel")	D

<sup>\*</sup> See the commands "Q:", "SRQ:" and "STS:".

#### (2) Operation status description

#### 01 Normal stop

Description	Positioning is completed within the in-position range.
Display	See "Positioning status" in "5.0perations".
Status command reply contents	K

# 02 During command move

Description	Command operation is being performed.
Display	See "Positioning status" in "5.0perations".
Status command reply contents	M

03 Out of the in-position range (After positioning is completed) (This state is READY)

Description	It is out of the in-position range after positioning is completed.
Display	See "Positioning status" in "5.0perations".
Status command reply contents	F

O4 Out of the in-position range (During fine adjustment) (This state is BUSY)

	This is the state where the positioning operation after the command operation is being executed.
Display	See "Positioning status" in "5.0perations".
Status command reply contents	G

05 During electrical origin return

Description	It is operating toward the electrical origin.
Display	-
Status command reply contents	P



2	3 4 5	6 7 8 9 10 11 12 13 14 15
06	During mechanical origin	return
	December	The incompation of the control of th
	Description Display	It is operating toward the machine origin.
	Status command reply	contents ()
	Status Communa Topy	
07	CW side limit stop	
	Description	It is a state stopping with CW limit sensor.
	Display	CW LMT (CW button flashing)
	Status command reply	
80	CCW side limit stop	
	Description	It is a state stopping with CCW limit sensor.
	Display	CCW LMT (CCW button flashing)
	Status command reply	
09	CW side software limit sto	qq
	Description	It is a state stopping with CW software limit.
	Display	CW SLMT (CW button flashing)
	Status command reply	
10	CCW side software limit s	stop
	Description	It is a state stopping with CCW software limit.
	Display	CCW SLMT (CCW button flashing)
11	Status command reply  CW side slowdown senso	·
11	Status command reply	·
11	Status command reply  CW side slowdown senso	or area
11	Status command reply  CW side slowdown senso  Description	It is in the CW side slowdown sensor.  CW S/D
	Status command reply  CW side slowdown senso  Description  Display  Status command reply	It is in the CW side slowdown sensor.  CW S/D  contents A
11	Status command reply  CW side slowdown senso  Description  Display	It is in the CW side slowdown sensor.  CW S/D  contents A
	Status command reply  CW side slowdown senso  Description  Display  Status command reply  CCW side slowdown senso	It is in the CW side slowdown sensor.  CW S/D  contents A  sor area
	Status command reply  CW side slowdown senso  Description  Display  Status command reply	It is in the CW side slowdown sensor.  CW S/D  contents A
	Status command reply  CW side slowdown senso  Description  Display  Status command reply  CCW side slowdown senso  Description	It is in the CW side slowdown sensor.  CW S/D  contents A  sor area  It is in the CCW side slowdown sensor.  CCW S/D
12	Status command reply  CW side slowdown senso  Description  Display  Status command reply  CCW side slowdown senso  Description  Display  Status command reply	It is in the CW side slowdown sensor.  CW S/D  contents A  sor area  It is in the CCW side slowdown sensor.  CCW S/D
	Status command reply  CW side slowdown senso  Description  Display  Status command reply  CCW side slowdown senso  Description  Display	It is in the CW side slowdown sensor.  CW S/D  contents A  sor area  It is in the CCW side slowdown sensor.  CCW S/D
12	Status command reply  CW side slowdown senso  Description  Display  Status command reply  CCW side slowdown senso  Description  Display  Status command reply	It is in the CW side slowdown sensor.  CW S/D  contents A  sor area  It is in the CCW side slowdown sensor.  CCW S/D
12	Status command reply  CW side slowdown sensor  Description  Display  Status command reply  CCW side slowdown sensor  Description  Display  Status command reply  Error occurred  Description  Display  Display	It is in the CW side slowdown sensor.  CW S/D  contents A  sor area  It is in the CCW side slowdown sensor.  CCW S/D  contents U  An error has occurred.  -
12	Status command reply  CW side slowdown sensor  Description  Display  Status command reply  CCW side slowdown sensor  Description  Display  Status command reply  Error occurred  Description	It is in the CW side slowdown sensor.  CW S/D  contents A  sor area  It is in the CCW side slowdown sensor.  CCW S/D  contents U  An error has occurred.  -
12	Status command reply  CW side slowdown sensor  Description  Display  Status command reply  CCW side slowdown sensor  Description  Display  Status command reply  Error occurred  Description  Display  Status command reply	It is in the CW side slowdown sensor.  CW S/D  contents A  sor area  It is in the CCW side slowdown sensor.  CCW S/D  contents U  An error has occurred.  - contents E
12	Status command reply  CW side slowdown sensor  Description  Display  Status command reply  CCW side slowdown sensor  Description  Display  Status command reply  Error occurred  Description  Display  Status command reply  Motor is transitioning to express the side of the sensor command reply  Motor is transitioning to express the side of the	It is in the CW side slowdown sensor.  CW S/D  contents A  sor area  It is in the CCW side slowdown sensor.  CCW S/D  contents U  An error has occurred.  - contents E
12	Status command reply  CW side slowdown sensor  Description  Display  Status command reply  CCW side slowdown sensor  Description  Display  Status command reply  Error occurred  Description  Display  Status command reply	It is in the CW side slowdown sensor.  CW S/D  contents A  sor area  It is in the CCW side slowdown sensor.  CCW S/D  contents U  An error has occurred.  - contents E  excitation  This is the state in which the motor is being shifted to the
12	Status command reply  CW side slowdown sensor  Description  Display  Status command reply  CCW side slowdown sensor  Description  Display  Status command reply  Error occurred  Description  Display  Status command reply  Motor is transitioning to express the side of the sensor command reply  Motor is transitioning to express the side of the	It is in the CW side slowdown sensor.  CW S/D  contents A  sor area  It is in the CCW side slowdown sensor.  CCW S/D  contents U  An error has occurred.  - contents E
12	Status command reply  CW side slowdown sensor  Description  Display  Status command reply  CCW side slowdown sensor  Description  Display  Status command reply  Error occurred  Description  Display  Status command reply  Motor is transitioning to expression  Description	It is in the CW side slowdown sensor.  CW S/D  contents A  Sor area  It is in the CCW side slowdown sensor.  CCW S/D  contents U  An error has occurred.  - contents E  excitation  This is the state in which the motor is being shifted to the excitation state.  -
12	Status command reply  CW side slowdown sensor  Description  Display  Status command reply  CCW side slowdown sensor  Description  Display  Status command reply  Error occurred  Description  Display  Status command reply  Motor is transitioning to expression  Display  Description  Display  Status command reply	It is in the CW side slowdown sensor.  CW S/D  contents A  sor area  It is in the CCW side slowdown sensor.  CCW S/D  contents U  An error has occurred.   contents E  excitation  This is the state in which the motor is being shifted to the excitation state.   contents H
12	Status command reply  CW side slowdown sensor  Description  Display  Status command reply  CCW side slowdown sensor  Description  Display  Status command reply  Error occurred  Description  Display  Status command reply  Motor is transitioning to expect the command reply  Description  Display  Status command reply  Motor is transitioning to expect the command reply  Motor is transitioning to expect the command reply  Motor is transitioning to expect the command reply	It is in the CW side slowdown sensor.  CW S/D  contents A  sor area  It is in the CCW side slowdown sensor.  CCW S/D  contents U  An error has occurred.   contents E  excitation  This is the state in which the motor is being shifted to the excitation state.   contents H
12	Status command reply  CW side slowdown sensor  Description  Display  Status command reply  CCW side slowdown sensor  Description  Display  Status command reply  Error occurred  Description  Display  Status command reply  Motor is transitioning to expression  Display  Status command reply  Motor is transitioning to expression  Display  Status command reply  Status command reply	It is in the CW side slowdown sensor.  CW S/D  contents A  sor area  It is in the CCW side slowdown sensor.  CCW S/D  contents U  An error has occurred.   contents E  excitation  This is the state in which the motor is being shifted to the excitation state.   contents H  This is the state in which the motor is being shifted to the excitation  This is the state in which the motor is being shifted to the excitation
12	Status command reply  CW side slowdown sensor  Description  Display  Status command reply  CCW side slowdown sensor  Description  Display  Status command reply  Error occurred  Description  Display  Status command reply  Motor is transitioning to expect the command reply  Description  Display  Status command reply  Motor is transitioning to expect the command reply  Motor is transitioning to expect the command reply  Motor is transitioning to expect the command reply	It is in the CW side slowdown sensor.  CW S/D  contents A  sor area  It is in the CCW side slowdown sensor.  CCW S/D  contents U  An error has occurred.   contents E  excitation  This is the state in which the motor is being shifted to the excitation state.   contents H



#### 16 Disabled axis (Not set by parameter "AXIS Sel")

	This is the state where all controls related to the axis are disabled. Set by the parameter "AXIS Sel". See "AXIS Sel" in "7. Parameterts"
Display content	-
Status command reply contents	D

#### 11-2.Error status

This status indicates the error status of the connection stage and this equipment.

#### (1) Error status list

No	Contents	Display	Reply *				
NO	Contents	Display	Format1	Format2 (13 bits)			
01	Normal (No error)	(Hidden)	K	1st bit is 1			
02	Command error	CMD ER	1	2nd bit is 1			
03	Scale error	SCALE ER	2	3rd bit is 1			
04	Limit stop	CW LMT • CCW LMT	3	4th bit is 1			
05	Over speed error	OS ER	4	5th bit is 1			
06	Overflow error	OF ER	5	6th bit is 1			
07	Emergency stop	EMERGENCY	6	7th bit is 1			
08	Interpolator error	IP ER	7	8th bit is 1			
09	Limit error	LIMIT ER	8	9th bit is 1			
10	System error	SYS ER	9	10th bit is 1			
11	Slowdown sensor area	CW S/D • CCW S/D	А	11th bit is 1			
12	Software Limit stop	CW SLMT • CCW SLMT	В	12th bit is 1			
13	TEACHING command error	TCMD ER	Т	13th bit is 1			

<sup>\*</sup> See the commands "Q:", "SRQ:" and "STS:".

#### Error status description

#### O1 Normal (No error)

Description		No error has occurred.
Display		(Hidden)
Status command	Reply format 1	K
reply contents	Reply format 2	1st bit is 1

#### 02 Command error CMD ER

Description		This error occurs when a command string that is not defined in this equipment or that does not match the status of this equipment is sent from the PC.
Display		CMD ER
Status command reply contents		1 2nd bit is 1

	1)	A command was sent that cannot be used in the current mode.
	2)	An ASCII code string other than the defined command was sent.
	3)	A code other than ASCII code was sent.
	4)	A normal delimiter was sent with a string that did not match the delimiter settings already sent to this equipment.
Occurrence case	5)	There is an unnecessary code in the communication buffer, and a command with this unnecessary code added was sent.
Case	6)	A command to specify an axis other than the controllable axis was sent.
	7)	The motor of the axis that issues the operation command is demagnetized.
	8)	In the limit sensor input state, a command to operate to the further limit sensor was transmitted.
	9)	In the state where the software limit is set, a command that operates to the further software limit is sent.
Release and recovery	whe	command can be sent even in the command error state, so it will be canceled on a command string defined in this equipment or in accordance with the status his equipment is sent.



	3 4	5	6	7	8	9	10	11	12	13	14	15	
03	Scale error SCALE ER												
	Description				Осси	rs whe	n no scale	e sign	al is in	nut			_
	Display					E ER	THE COUNT	0.61	<u> </u>	рац.			
	Status command Reply format 1												
	reply content	ts R	eply forr	nat 2	3rd b	it is 1							
		14)	Tlan and	ll-l-					I 4	+ -:		_	_
	Occurrence		The sca The sca				s disconn	iectec	rom	tnis eq	ulpmer	11.	_
	case		Encoder			(O) 1.							
	Release and	Turn	off the	oower o	of the ir	nstrume	ent, remov	ve the	cause	e, and t	hen tu	rn on t	he
	recovery	powe	r.										
													_
04	Limit stop			CW L	_MT • C	CW LN	lΤ						
	December			- t- t	+lc = 01	N as 0014	/ 1: :-	00:5 - :	r o t - · ·	0 0 m l= 1	tla =-:-	_	
	Description	F	or CW 1ir	mit			V or CCV W button			OLON	e or bo	ırı axes	<u>ة.</u>
	Display		or CCW 1				CCW but			)			_
	Status comn				3								_
	reply content		eply forr		4th b	it is 1							_
	D 1	lı.			1	1	1.						
	Release and recovery	It can direct		perate	d in the	limit se	ensor dire	ction	any m	ore. Op	erate ir	n the o	p
	1000 (01)	GII OO											_
	Description				Occurs when operating at a speed higher than the coucapability of the coordinate value counter.								
	Display Status comr	mand B	enly form	nat 1	OS ER 4								
	reply content		eply forr		5th bit is 1								
				dad wit	l_ +l								
			Somethi										
	Occurrence	2)	The stag	ge is un	dergoir	ng vibra	tion.	.1					_
	Occurrence case	2)	The stag Strong r	ge is un noise is	dergoir mixed i	ng vibra n the so		ıl.					_
		2) 3) 4) After	The stag Strong r A strong removin	ge is un noise is g flash l	dergoir mixed ii nit the s	ng vibra n the so stage.	tion.		ower s	upply o	r send	the cor	
	case	2) 3) 4)	The stag Strong r A strong removin	ge is un noise is g flash l	dergoir mixed ii nit the s	ng vibra n the so stage.	tion. cale signa		ower s	upply o	r send	the cor	าาเ
	case Release and recovery	2) 3) 4) After "RESI	The stag Strong r A strong removin	ge is un noise is g flash I	dergoin mixed ii nit the s ause, re	ng vibra n the so stage.	tion. cale signa		ower si	upply o	r send	the cor	
06	case Release and	2) 3) 4) After "RESI	The stag Strong r A strong removin	ge is un noise is g flash l	dergoin mixed ii nit the s ause, re	ng vibra n the so stage.	tion. cale signa		ower s	upply o	r send	the cor	 mi
<del>.</del> 26	Release and recovery  Overflow error	2) 3) 4) After "RESI	The stag Strong r A strong removin	ge is un noise is g flash I	dergoir mixed ii nit the s ause, re	ng vibra n the so stage. estart o	tion. cale signa r restart t	he po					
<del>.</del> 26	case Release and recovery	2) 3) 4) After "RESI	The stag Strong r A strong removin	ge is un noise is g flash I	dergoir mixed in nit the s ause, re	ng vibra n the so stage. estart o	tion. cale signa	he po	ce bet	ween t	the coc		
<b>D</b> 6	Case  Release and recovery  Overflow error  Description  Display	2) 3) 4) After "RESI	The stag Strong r A strong removin ET:".	ge is un noise is g flash I g the c	dergoir mixed in nit the sause, re	ng vibra n the so stage, estart o rs whe he spec	tion. cale signa r restart t	he po	ce bet	ween t	the coc		
<u></u>	Case  Release and recovery  Overflow error  Description  Display  Status comm	2) 3) 4) After "RESI	The stag Strong r A strong removin ET:".	ge is un noise is g flash I g the c	dergoir mixed in nit the s ause, re  R  Occu and t  OF El	ng vibra n the so stage. estart o rs whe he spec R	tion. cale signa r restart t	he po	ce bet	ween t	the coc		
06	Case  Release and recovery  Overflow error  Description  Display	2) 3) 4) After "RESI	The stag Strong r A strong removin ET:".	ge is un noise is g flash I g the c	dergoir mixed in nit the sause, re	ng vibra n the so stage. estart o rs whe he spec R	tion. cale signa r restart t	he po	ce bet	ween t	the coc		
<del>.</del> 26	Case  Release and recovery  Overflow error  Description  Display  Status comm	2) 3) 4) After "RESI	The stag Strong r A strong removin ET:".	ge is un noise is g flash I g the c OF E	dergoir mixed in nit the s ause, re  R  Occu and t  OF El 5 6th b	ng vibran the sestage.  Pestart of the special sestart of the special sessinates of the special sestart of the special sestart of the spe	tion. cale signa r restart t  n the dif- cified pos	he po	ce bet	ween t	the coc		
<del>-</del> 06	Case  Release and recovery  Overflow error  Description  Display  Status comm	2) 3) 4) After "RESI or	The stag Strong r A strong removin ET:".	ge is un noise is g flash I g the c OF E	dergoir mixed in nit the s ause, re  R  Occu and t OF El 5 6th b	ng vibran the sestage. estart of the special s	tion. cale signa r restart t  n the dif- cified pos	the po	ce bet 55mm	ween t	the coc re.	ordinate	e
206	Case  Release and recovery  Overflow error  Description  Display  Status comm	2) 3) 4) After "RESI or	The stag Strong r A strong removin ET:". eply forr eply forr The stag After the	ge is unnoise is g flash I g the cooperation of E	dergoin mixed in mixed in it the sause, resause, resause, resause, resause, results of the sause in its sause, results of the sause in its sause in	ng vibran the sestage.  rs whee he specential is 1  tep-out was core	tion. cale signa r restart t  n the dif- cified pos	ferencition is	ce bet 55mm	ween t	the coc re.	ordinate	e
<del>)</del>	Case  Release and recovery  Overflow error  Description  Display  Status common reply content	2) 3) 4) After "RESI or  mand R ts R	The stage Strong removing The stage After the by 5 mm When the motor care strong to the stage of	ge is unnoise is g flash I g the commat 1 mat 2 ge got ce position or mone stage able an	dergoir mixed in nit the sause, re R Occu and t OF El 5 6th b out of so oning v re. (Only	rs whee special specia	tion.  cale signa  r restart t  n the dif- cified posi	ferendition is	ce bets 5mm	ween to or mo	the coore. s turne	ordinate	e m
<del>)</del>	Case  Release and recovery  Overflow error  Description  Display  Status common reply content  Occurrence	2) 3) 4) After "RESI or  mand R ts R	The stage Strong removing The stage After the by 5 mm When the motor cathan 5 resistance.	ge is unnoise is g flash I bg the commat 1 mat 2 ge got ce position or mone stage able annm)	dergoir mixed in nit the sause, re  Occur and t OF El 5 6th b  out of so oning v re. (Only e moves d the s	rs whee special specia	tion. cale signa r restart t  n the dif- cified positive positive permise at permise	ferendition is	ce bets 5mm	ween to or mo	the coore. s turne	ordinate ed and on axis	e m



2	3 4	5 6 7	8 9 10 11 12 13 14 15			
07	Emergency st	top EM	MERGENCY			
	Description		Emergency stop,			
	Display		EMERGENCY			
	Status comn	nand Reply format 1				
	reply content	s Reply format 2	2   7th bit is 1			
	Release and	See "(5) Emergency	cy stop".			
	recovery					
00	1 . 1 .	ID F	<u></u>			
80	Interpolator e	error IP E	ER			
	Description		Occurs when the magnitude of the scale signal is out o			
		_	specified range.			
	Display Status comp	nand Reply format 1				
	reply content					
	Occurrence	1) The scale is d	· · · · · · · · · · · · · · · · · · ·			
	case	<ul><li>2) A strong flash</li><li>3) Encoder is bro</li></ul>	h hit the stage.			
	Release and		er of the instrument, remove the cause, and then turn on the			
	recovery	power.				
00	1		AT ED			
09	Limit error	LIIV	MIT ER			
	Description		Occurs when the CW and CCW limits are on at the sai			
	Display		time.			
		nand Reply format 1				
	reply content					
			able connector has disconnected from this equipment.			
			and CCW limit sensor wires of the motor cable are broken. Two or CCW limit sensor wire of the motor cable has l			
	Occurrence	- 5	d, and the limit has been entered on the opposite side from			
	case	disconnected	·			
		4) Dirt or foreig	n matter has entered both the CW and CCW limit senso ne stage			
			sor mounted on the stage has failed.			
		If you want to kee	ep the coordinate values after removing the cause, perfo			
	Release and		It it is not incompany to maintain the coordinate values roots			
	Release and recovery		, and the second second second second second second second second second second second second second second se			
			supply, or send the command "RESET:".			
10		restart the power s	· · · · · · · · · · · · · · · · · · ·			
10	System error	restart the power s	supply, or send the command "RESET:".  ('S ER			
10	System error  Description	restart the power s	supply, or send the command "RESET:".  'S ER  Occurs when this quipment system is out of order.			
10	System error  Description Display	restart the power's	supply, or send the command "RESET:".  /S ER  Occurs when this quipment system is out of order.  SYS ER			
10	System error  Description Display	restart the power's SYS	Supply, or send the command "RESET:".  /S ER  Occurs when this quipment system is out of order.  SYS ER  1 9			
10	System error  Description Display Status common reply content	restart the power's  SYS  mand Reply format 1 Reply format 2	Supply, or send the command "RESET:".  (SER  Occurs when this quipment system is out of order.  SYS ER  1 9 2 10th bit is 1			
10	System error  Description Display Status common reply content  Release and	restart the power's  SYS  mand Reply format 1 Reply format 2	supply, or send the command "RESET:".  (SER  Occurs when this quipment system is out of order.  SYS ER  1 9 2 10th bit is 1  beled or restored. Unplug the power cable from the outlet			
10	System error  Description Display Status common reply content	restart the power's  SYS  mand Reply format 1 Reply format 2	Supply, or send the command "RESET:".  (SER  Occurs when this quipment system is out of order.  SYS ER  1 9 2 10th bit is 1			
10	System error  Description Display Status common reply content  Release and	restart the power's  SYS  nand Reply format 1 Reply format 2  It cannot be cance contact our compa	Supply, or send the command "RESET:".  (SER  Occurs when this quipment system is out of order.  SYSER  1 9 2 10th bit is 1  celed or restored. Unplug the power cable from the outlet			
	System error  Description Display Status comm reply content  Release and recovery	restart the power's  SYS  nand Reply format 1 Reply format 2  It cannot be cance contact our compa	Supply, or send the command "RESET:".  (SER  Occurs when this quipment system is out of order.  SYS ER  1 9 2 10th bit is 1  celed or restored. Unplug the power cable from the outlet any or our distributor.			
	System error  Description Display Status comm reply content  Release and recovery	restart the power's  SYS  nand Reply format 1 Reply format 2  It cannot be cance contact our compa	Supply, or send the command "RESET:".  Occurs when this quipment system is out of order.  SYS ER  1 9 2 10th bit is 1  celed or restored. Unplug the power cable from the outlet any or our distributor.  W S/D • CCW S/D  Either one-axis or both-axis CW and			
	System error  Description Display Status commreply content  Release and recovery  Slowdown se	restart the power's  SYS  mand Reply format 1 s Reply format 2  It cannot be cance contact our compa	Supply, or send the command "RESET:".  Occurs when this quipment system is out of order.  SYS ER  1 9 2 10th bit is 1  celed or restored. Unplug the power cable from the outlet any or our distributor.  W S/D • CCW S/D  Either one-axis or both-axis CW and slowdown sensors are included.			
	System error  Description Display Status commreply content  Release and recovery  Slowdown se	restart the power's  SYS  mand Reply format 1  Reply format 2  It cannot be cance contact our comparensor input  CW  For CW slowdo	Supply, or send the command "RESET:".  Occurs when this quipment system is out of order.  SYS ER  1 9 2 10th bit is 1  Deled or restored. Unplug the power cable from the outlet any or our distributor.  W S/D • CCW S/D  Either one-axis or both-axis CW and slowdown sensors are included.  Own sensor input CW S/D			
	System error  Description Display Status common reply content  Release and recovery  Slowdown see  Description Display	restart the power's  SYS  mand Reply format 1  Reply format 2  It cannot be cance contact our comparensor input  CW  For CW slowdo	Occurs when this quipment system is out of order.  SYS ER  1 9 2 10th bit is 1  Deled or restored, Unplug the power cable from the outlet any or our distributor.  W S/D • CCW S/D  Either one-axis or both-axis CW and slowdown sensors are included.  Own sensor input CW S/D  rdown sensor input CCW S/D			



1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16

### 12 Software limit stop

$\bigcirc \bigvee$	Ō	Λ/T	• CCW	CI M	Ŧ
$\bigcirc$ v	0	IVII			

Description		One or both axes CW or CCW software limit is included.
II )ienlav		CW SLMT (CW button flashing) CCW SLMT (CCW button flashing)
Status command		B
reply contents	Reply format 2	12th bit is 1

Release and It cannot be operated in the limit sensor direction any more. Operate in the opposite direction.

## 13 Teaching command error TCMD ER

Description		The content of the teaching registration line to be executed does not match the status of this equipment.
Display		TCMD ER
Status command	Reply format 1	T
reply contents	Reply format 2	13th bit is 1

	In the limit sensor input state, it was the contents of the registration line that operates on the limit sensor side that is further included.
Occurrence case	Set the operation speed (command "D:") in the registration line, and then 2) set the parameter "Max Speed" to a value less than the registered operation speed (command "D:") and execute the registration line.
Release and recovery	If you want to keep the coordinate values after removing the cause, perform a busy error cancel. If it is not necessary to maintain the coordinate values, restart or restart the power supply, or send the command "RESET:".



# 12. Specification

These are the specifications for this equipment.

#### 12-1. Basic performance

	Model	Minimum command unit	In-position range * 1	Maximum speed setting	The maximum amount of movement setting range
	FC-114	100nm	$\pm$ 100, $\pm$ 300, $\pm$ 700nm	100mm/sec	$-13421.7728 \sim +13421.7727$ mm
	FC-414	50nm	$\pm$ 50, $\pm$ 150, $\pm$ 350nm	100mm/sec	$-6710.88640 \sim +6710.88635$ mm
ſ	FC-514 10nm $\pm$ 10, $\pm$ 30, $\pm$ 70nn		$\pm$ 10, $\pm$ 30, $\pm$ 70nm	50mm/sec	$-1342.17728 \sim +1342.17727$ mm

 $<sup>\</sup>ast$  1 Set the in-position range with the parameter "INPos Range" .

	Item	Contents	
Stage control axes		2	
		Command error, Scale error, Limit stop	
		Overspeed error, Overflow error, Emergency stop	
Error detection, etc.		Interpolator error, Limit error, System error	
		Software limit stop, Teaching command error	
Number of control	Out of the state	Jog controller, Emergency stop input	
interface ports	One for each	GP-IB, USB, Ethernet, General purpose I / O	

#### 12-2. General specifications

ltem	Contents
Power source	AC100V ~ 240V, 50/60Hz
Allowable variations of voltage	AC90V ~ 264V
Power consumption	110VA max
Fuse	250V, 2.5A, Time lag, 2 used
External dimensions	W220 × H88 × D290mm
Weight	5.2kg
Operating temperature	0°C~ 40°C
Operating ambient humidity	$20\% \sim 80\%$ RH (No condensation)
Storage temperature	-10℃~55℃
Ambient storage humidity	$20\% \sim 80\%$ RH (No condensation)
Place of use	Indoor
Storage altitude	up to 2000m
Operating altitude	up to 2000m
Maximum operating time	Continuous operation possible

#### 12-3. Safety and electromagnetic compatibility

Item	Contents			
Target model	FC-514			
Safety	EN61010-1:2010 compliant, Overvoltage category II, Pollution degree 2			
	EN61326-1:2013 compliant, EN61000-3-2:2014 compliant EN61000-3-3:2013 compliant			

Item	Conditions					
Item						
	Cable (2m or less)	Motor	Attach a ferrite core to one end (body side) of the cable. *3			
		Scale	Attach ferrite cores to both ends of the cable. *3			
		Jog controller	Attach femite cores to both ends of the cable. *3			
	Cable (3m or less)	Emergency stop	Wrap one end (body side) of the prepared cable			
Electromagnetic		(Twisted pair, no shield)	twice around the ferrite core. *4			
compatibility		GP-IB (shield) *2				
		USB (shield) *2	Attach a ferrite core that matches the external			
		Ethernet (shield) *2	shape of the prepared cable to one end (body side) of the cable. *4			
		汎用 I/O(shield)	Side/ Of the capie, '-			
	General	Connector hood	EMI measures			
	purpose I / O	Housing	General-purpose I / O board stored in metal box			

<sup>\*2</sup> Connect only the interface that communicates



<sup>\*3</sup> Ferrite core model: ZCAT 1730-0730A: made by TDK

<sup>\*4</sup> Ferrite core model: ZCAT series: made by TDK

# 12-4. Interface specifications (1) GP-IB

24 13

#### Connection

This instrument is equipped with a GP-IB connector and can be connected to a personal computer. When using, prepare a GP-IB cable.

#### Communication

In order for the PC to communicate with the instrument, the PC must have GP-IB communication-capable software (such as our sample software). Also, it is necessary to set the mode of this equipment to REMOTE and set the parameter "I / F Sel" to GP-IB. The GP-IB address must not be the same as other devices. Check the connection between the PC and this equipment by sending and receiving operation information commands (such as "Q:").

#### To disconnect

Close GP-IB from the PC software. When disconnecting the cable, there is no problem even if this equipment is turned on. However, it is recommended that this equipment be turned off for safety. For information on personal computers, see the instruction manual for your personal computer.

#### · Specification

ltem		Contents
	SH1	Source handshake all functions
	AH1	Acceptor handshake all functions
	T6	Basic talker function, Serial poll function, Talker cancellation by MLA
	TEO	No extended talker function
	L4	Basic listener function, Release of listener by MTA (no listen only)
Function	LEO	No function
Lanction	SR1	Service request all functions
	DC2	No SDC function
	RL2	No local lockout function
	DTO	No device trigger function
	PPO	No parallel pole function
	CO	Without controller function
Address		1~30
Delimiter		CR+LF, EOI, CR, LF
Service request		Enabled or disabled
Flow control		None (fixed)
Connector used		Manufacturer: DDK Corporation Model: 57LE-20240 (57LE Series)
Applicable plug, cable		IEEE-488 compliant product, GP-IB compliant product



• Pin assignment and circuit diagram

Pin	Signal		NAT9914		Vcc				
Number	Name		INAISSIA	Vcc		91	I75160		IEEE.488
1	DIO1	Vcc	/ACCRQ	TR	_	PE OI	GND GND	m	Compliant connector
2	DI02		/ACCGR	DIO1		DIO1	DIO1		1 O DIO1
3	DIO3	-	/CE	DIO2		DIO2	DIO2		2 O DIO2
4	DI04			DIO3		DIO3	DIO3		3 O DIO3
5	EOI	• -	/WE	DIO4		DIO4	DIO4		12 O DIO4
6	DAV		DBIN	DIO5		DIO5	DIO5		13 O DIO5
7	NRFD			DIO6		DIO6	DIO6		14 O DIO6
8	NDAC			DIO7		DIO7	DIO7		15 O DIO7
9	IFC			DIO8		DIO8	DIO8		—10 DIO8
10	SRQ				1/0/	Vcc	TE		18
11	ATN				Vc.	<u> </u>			18 O GND(6)
12	SHILD	-	RS0						O GND(7)
13	DI05	_	RS1			SN	l75161		$\int_{-\infty}^{20} \bigcirc GND(8)$
14	DIO6	-	RS2	/CONT		DC O	GND	m	10 0 0 0 0
15	DI07			SRQ		SRQ	SRQ		USRQ
16	DI08	_	/INT	ATN		ATN	ATN		CAIN
17	REN	_	D0(MSB)	EOI		EOI	EOI		TOEUI
18	GND	_	D1	DAV		DAV	DAV		
19	GND	_	D2	NRFD		NRFD	NRFD		NRFD 8 ONDAG
20	GND	_	D3	NDAC		NDAC	NDAC		ONDAC
21	GND	_	D4 D5	IFC REN		IFC REN	IFC		O IFC
22	GND	_	_	KEN		NC NC	REN TE		TO REN
23	GND	_	D6 D7			Vcc	SC		21 O GND(9)
24	GND		/CLK		Vc		30		22 O GND(10)
24	GIND		/RESET			_			23 O GND(11)
		_	/KESEI					/	A OBIND(11)
		m	Vss	TE					GND LOGIC  12  SHIELD

(2) USB



#### Connection

This instrument is equipped with a USB typeB connector (Function) and can be connected to a USB typeA connector (Host) of a personal computer. When using, please prepare USB1.0, 1.1, 2.0, 3.0, 3.1 compatible cable \*1.

\*1 Standard-A (male) -Standard-B (male) Connector straight cable (USB1.0, USB1.1, USB2.0, USB3.0, USB3.1)

#### Driver installation (for Windows 10)

When connecting for the first time, connect the USB cable, turn on the PC, and then turn on the power of the instrument. The driver is automatically installed on the PC and recognized as a COM port.

#### Driver installation (for Windows 7, 8, 8.1)

The setting information file "stage\_controller\_usb\_cdc\_drv\_w \*\*\*. Inf" is required. This setting information file can be obtained from the download page of our homepage (note that the setting information file differs depending on the version of Windows OS). Start the installation by manually specifying the setting information file when installing the driver. After the driver is installed, "Stage Controller Usb Cdc Port (COM \*)" will be created in the port of the PC

The COM number of "Stage Controller Usb Cdc Port (COM \*)" is changed for each USB port. Once the driver is installed, the newly added USB port will be automatically installed.

The port is recognized as a virtual COM port, In order for the PC to communicate with the DX, software that enables serial communication with the PC is required. Also, it is necessary to set the mode of this equipment to REMOTE and set the parameter "I / F Sel" to USB. Check the COM port number with Device Manager. The baud rate (communication speed) of the COM port can be set to any number (example: 9600 bps). Check the connection between the PC and this equipment by sending and receiving operation information commands (such as "Q:").

#### To disconnect

When disconnecting, be sure to close the COM port with the PC software before disconnecting. When disconnecting the cable, there is no problem even if this equipment is turned on. However, it is recommended that this equipment be turned off for safety. For information on personal computers, see the instruction manual for your personal computer. If the power to the instrument is restarted or restarted without closing the COM port from the PC software, it may be necessary to reconnect the cable after closing the COM port from the software to relink. (Some software automatically closes the COM port and automatically connects after restarting.)

#### Specification


Item		Contents		
Function		Used as a virtual COM port,		
Number of port		1port		
Transfer speed		Supports full-speed transfer (12Mbps)		
Delimiter		CR+LF, CR, LF		
0	Manufacturer	OMRON Corporation		
Connector used	Model	XM7B-0442		
useu	Туре	Standard-B		
Applicable plug, cable		「Standard-A (male)」connector - 「Standard-B (male)」connector Straight cable (USB1.0, USB1.1, USB2.0, USB3.0, USB3.1)		

#### Pin assignment

Pin number	Signal name	
1	VBUS	
2	D-	
3	D+	
4	GND	



#### (3) Ethernet



#### Connection

This device is equipped with an Ethernet connector, and can be connected to a personal computer directly or via a hub (router). LAN cable can be either cross or straight. Use a cable category that supports 10Mbps and 100Mbps transfers. When connecting, you need to set the IP address, default gateway, and subnet mask. Set the parameters of this equipment according to the information of the IP address, default gateway, and subnet mask of the personal computer used. For the IP address setting, if the IP address of the PC used as an example is "192.168.015.188", set "188" to another value. Set the default gateway and subnet mask to the same values as on the PC. If you use the same LAN area, you do not need to set the default gateway.

#### Communication

In order for the PC to communicate with this instrument, the PC must have software capable of TCP / IP communication (sample software provided by NEC). Also, it is necessary to set the mode of this equipment to REMOTE and set the parameter  $\H$ I / F Sel  $\H$ to Ethernet. Check the connection between the PC and this equipment by sending and receiving operation information commands (such as  $\H$ Q: $\H$ ). When communicating with a communication application, it is necessary to set the port number in addition to the IP address settings.

#### Specification

ltem	Contents
Number of port	1
Product specifications	IEEE802.3x Flow control compliant
Transfer speed	Supports 10Mbps and 100Mbps transfer
Delimiter	CR+LF, CR, LF
IP address	000.000.000.000 ~ 255.255.255.255
Default gateway	000.000.000.000 ~ 255.255.255.255
Subnet mask	000,000,000,000 ~ 255,255,255,255
Port number	60000 (fixed)
Connector used	RJ45

#### · Pin assignment


Pin number	Signal name
1	TX+
2	TX-
3	RX+
4	NC
5	NC
6	RX-
7	NC
8	NC

#### · LED

LED color	Contents			
Green lighting	Lights when the communication speed is 100Mbps, and turns off when the communication speed is 10Mbps or when no cable is connected.			
Yellow lighting	Illuminates when a valid link is detected and flashes when data transmission or reception is detected.			



#### (4) General purpose I / O



This section describes the functions of each pin of the general-purpose I / O.

#### Specification

l1	tem	Contents
		General purpose input $ imes$ 3 port
	Input	Teaching operation $\times$ 1
		Busy error cancel X 1
Function		General purpose output $ imes$ 3 port
Function		Scale division pulse signal × 2 axes
	Output	Alarm signal $\times$ 2 axes
		In-position signal $ imes$ 2 axes
		Teaching state $\times$ 1
Connector	Manufacturer	3M Japan Co., Ltd.
used	Model	50 pin half pitch connector (MDR) 10150-5202PL
Applicable	Manufacturer	3M Japan Co., Ltd.
Applicable plug	Model	50 pin half pitch connector (MDR) 10150-3000PE
pius	iviodei	50 pin half pitch connector (MDR) 10150-6000**

#### · Pin assignment and function explanation

i. Input terminal

The internal circuit and specifications are shown in "Internal Circuit and Specifications""Input Terminal".

Terminal Number	Terminal Name	Contents
3	General purpose input 1	
28	General purpose input 2	Select each general purpose input. The status can be checked with the command "I:".
4	General purpose input 3	Checked with the Command 1.
29	General purpose input COMMON	Common terminal for general purpose input.
43	T_1	Select a teaching number. See the table below for
19	T_2	patterns. Patterns other than those in the table below
44	T_3	are invalid.
20	T_START	Start teaching.
45	T_PAUSE	Pauses teaching execution.
21	T_ONE_STEP	Execute teaching one line at a time for each input.
46	T_STOP	Stop the teaching execution.
22	T_COMMON	Common terminal for teaching operation input.
24	BE_CANCEL	Executes busy error cancellation. See "BEC" or the command "BEC:" for details.
49	BE_COMMON	Common terminal for BE CANCEL input.

#### Teaching number pattern table

Terminal			Channel		
l errillial	1	2	3	4	5
T_1	ON	OFF	ON	OFF	ON
T_2	OFF	ON	ON	OFF	OFF
T_3	OFF	OFF	OFF	ON	ON



ii. Output terminal
The internal circuit and specifications are shown in "Internal Circuit and Specifications""Output Terminal".

Terminal Number	Terminal Name	Contents	Output circuit
26	General purpose output1 General purpose output2	Outputs the status selected by the "O:" command.	Open
2	General purpose output3	It is off at startup.	collector
27	General purpose output COMMON	Common terminal for general purpose output.	-
8	1PA	A-phase output of scale-divided pulse (2-phase square wave) on First axis side	
33	1/PA	Inverted output of 1PA	Line
9	1PB	B-phase output of scale-divided pulse (2-phase square wave) on First axis side	driver
34	1/PB	Inverted output of 1PB	
10	1AL	Alarm output for First axis. Outputs a HIGH level when any of an interpolator error, scale error, limit error, or overflow error occurs. See "11. Status" for details of each error. In the case of only a limit error, if a normal stage is connected to the motor cable, it returns to the LOW level.	
35	1/AL	Inverted output of 1AL	
11	1INP	In-position output for First axis. Outputs HIGH level upon completion of positioning and detection by limit sensor, and outputs LOW level when an operation command is input. However, if an operation command is input during detection of the limit sensor, a LOW level is output when the signal passes through the limit sensor.	
36	1/INP	Inverted output of 1INP	
37	2PA	A-phase output of scale-divided pulse (2-phase square wave) on Second axis side	
13	2/PA	2PA inverted output.	
38	2PB	B-phase output of scale-divided pulse (2-phase square wave) on Second axis side	Line
14	2/PB	Inverted output of 2PB	driver
	2AL	Alarm output for Second axis. Outputs a HIGH level when any of an interpolator error, scale error, limit error, or overflow error occurs. See "11. Status" for details of each error. In the case of only a limit error, if a normal stage is connected to the motor cable, it returns to the LOW level.	
15	2/AL	Inverted output of 2AL	
40	2INP	In-position output for Second axis. Outputs HIGH level upon completion of positioning and detection by limit sensor, and outputs LOW level when an operation command is input. However, if an operation command is input during detection of the limit sensor, a LOW level is output when the signal passes through the limit sensor.	
16	2/INP	2INP inverted output.	
17	T_STATUS	Execution status output of teaching registration contents	
42	/T_STATUS	Inverted output of T_STATUS	
5, 6, 25, 30, 47	+5V	5V power supply terminal.  The maximum output current of each terminal is 0,2A. Do not exceed 1A in total.	-
7, 12, 18, 23, 31, 32, 41, 48, 50	GND	GND terminal for 5V power supply.	-



- · Internal circuit and specifications
  - i. Input terminal


Input terminals include general-purpose input, teaching operation, and busy error cancel. The input is a photocoupler LED. The phototransistor is turned on by passing current through the LED, and the phototransistor is turned off by interrupting the current. In the case of general-purpose input, the status of this phototransistor can be checked with the command "I.". See "I:" in "6. Commands" for details. For teaching operation and busy error cancellation, the function "Pin assignment and function explanation" can be executed by turning on the phototransistor.


 $\bigwedge$ CAUTION The input current must not exceed 20mA. Exceeding this may cause a failure.

- Use LED current within the range of 5 to 20mA.
- The input current waveform should be pulsed (rise and fall time <100  $\mu$ sec).

• The time width of ON and OFF of general-purpose input should be at least longer than the transmission cycle of command "I.".

 When connecting the COMMON terminal to the GND of this instrument, use the power supply of the input terminal at 5V of this instrument,







#### ii. Output terminal


Ո


#### a. Open collector

The output is a phototransistor of a photocoupler, which is open collector. The photo transistor can be turned ON and OFF with the command "O:". See "O:" in "6. Commands" for details.


CAUTION Vceo should not exceed 30V and Ic should not exceed 150mA. Exceeding this may cause a failure.


When connecting the COMMON terminal to the instrument's GND, use the instrument's 5V power supply.





#### b. Line driver







#### (5) Emergency



The contact type of the input contact is "B contact". The contact type cannot be changed. After changing the parameter "EMG Connector" to "Enable" and releasing the short circuit between the input and GND, the operation of the stage connected to the instrument is immediately stopped and "EMERGENCY" is displayed on the display. The state of the motor at this time depends on the setting of the parameter "EMG Motor Excite". After that, operations related to the stage operation cannot be performed. To recover, short-circuit the input and GND, and then cancel the busy error to maintain the coordinate value. If it is not necessary to maintain the coordinate values, restart the power, restart, or send the command "RESET:" . See "BEC" or the command "BEC:" for busy error cancellation. See "EMG Motor Excite" for the parameter "EMG Connector" and the parameter "EMG Motor Excite" for setting the motor excitation and demagnetization during an emergency stop.

When using, be sure to change the parameter "EMG Connector" to Enable and perform a test operation to confirm that it can be used.

↑CAUTION Do not connect anything that outputs power, such as an AC adapter.

#### Specification


	Item	Contents
Connector	Manufacturer	Hosiden Corporation
used	Model	HEC3800-01-010 DC power jack (JEITA RC-5320A compliant)
Applicable p	plug, cable	JEITA RC-5320A TYPE4 Twisted pair

#### (6) Motor driver setting switch

The input terminal motor driver is set for each axis. When the right side of the stage controller is viewed from the front, the left side is the setting switch for First axis and the right side is the setting switch for Second axis. The settings to be made are the stop current, drive current, and number of divisions. At the time of shipment, it is set according to the stage purchased at the same time. If you want to change the settings, please contact us or our distributor.

Depending on your environment, you may be required to change the settings, but do not change anything else. If you change it intentionally, the stage may move unintentionally.

#### • Settings (FC-114, FC-414)



#### i, Stop current (STOP)

Set the stop current value as a percentage (%) of the drive current.

SW No.	0	1	2	თ	4	15	6	7	8	Ø	Α	В	O		Ш	F
%	25	30	35	41	45	50	55	59	63	67	71	75	79	83	87	91

#### ii. Drive current (RUN)

Set the current value when driving the motor.

SW No.	0	1	2	З	4	5	6	7	8	9	Α	В	С	D	Е	F
Current value(A)	0.35	0.44	0.52	0.59	0.67	0.75	0.83	0.9	0.98	1.05	1.12	1.19	1.27	1.34	1.4	1.48


# iii. Setting the number of divisions (M) Set the number of divisions, \*

SW No.	0	1	2	3	4	5	6	7	8	9	Α	В	С	D	Ε	F
Division number	1	2	4	5	8	10	20	40	80	16	25	50	100	125	200	250

<sup>\*</sup> The step angle for the number of divisions is "step angle = basic step angle (0.72 or 0.36  $^{\circ}$  ) / number of divisions"



#### • Settings (FC-514)



# iv. Setting the number of divisions (Step)

Set the number of divisions. \*

SW No.	0	1	2	თ	4	5	6	7	ω	9	Α	В	O		Ш	F
Division number	1	2	4	10	20	40	100	200	400	800	ı	ı	-	ı	-	-

<sup>\*</sup> The step angle for the number of divisions is "step angle = basic step angle (0.72 or 0.36  $^{\circ}$  ) / number of divisions"

#### v. Drive current (Drive)

Set the current value when driving the motor.

SW No.	0	1	2	3	4	5	6	7	8	9	Α	В	O	D	Ш	F
Current value(A)	0.3	0.325	0.35	0.375	0.4	0.425		0.475	0.5	0.525	0.55	0.575	0,6	0.65	0.7	0.75


vi. Stop current (Hold) Set the stop current value as a percentage (%) of the drive current.

	SW No.	0	1	2	თ	4	5	6	7	8	9
ſ	%	10	20	30	40	50	60	70	80	90	100




### 13. Dimensions

#### 13-1. Front panel





### 13-2.Rear panel




#### 13-3. Right side panel

(1) FC-114, FC-414










## 13-4.Left side panel



# 14. Trouble shooting

If a problem occurs, check the following. If this does not solve the problem, unplug the power cable from the outlet and contact our company or our distributor.

Contents	Possible cause	Workaround	Page
Can not turn on.	Power cable is not connected.	Connect the power cable.	6
The power turned off during use.	The fuse has blown.	Check the fuse and replace it if it is blown. If it cuts off frequently, the instrument may be malfunctioning. Unplug the power cable from the outlet and contact our company or our distributor.	7
	The MAIN POWER switch on the rear panel is not turned on.	Turn on the MAIN POWER switch.	8
	power supply is short-circuited	step 1 Turn off the MAIN POWER switch, disconnect all the cables of the connected peripheral devices, and wait at least 10 seconds.  step 2 Turn ON the MAIN POWER switch with only the power cable connected. If the beep continues to sound, proceed to step 4.  step 3 Check that the Off lamp is lit, then press and release the POWER button for one second.  step 4 If the power cannot be turned on, the beep continues to beep even if it can be turned on, or if the Off lamp goes off while the POWER button is being pressed, the instrument may be malfunctioning. Unplug the power cable from the outlet and contact our company or our distributor. If the power can be turned on, one of the peripheral devices may be faulty or a device with a different pin assignment may be connected, Check the peripheral devices.	8,
The display screen is off.	The POWER button on the front panel is not ON.	Press the POWER button to turn it on	3
	The MAIN POWER switch on the rear panel is not turned on	Turn on the MAIN POWER switch.	8
	The computer is sleeping.	Wake up from sleep.	5
<ul><li>The menu is off.</li><li>Certain buttons cannot be operated.</li></ul>	· ·	Switch to LOCAL or TEACH mode.	13
GP-IB communication is not possible.	Communication interface setting is other than GP-IB.	Select GP-IB in parameter "I / F Sel".	85
	Communication settings do not match.	Check the communication conditions and set the parameters.	86
USB communication is not possible.	Communication interface setting is other than USB.	Select USB in parameter "I / F Sel".	85
	Communication settings do not match.	Check the communication conditions and set the parameters.	85
Ethernet communication is not possible.	Communication interface setting is other than Ethernet	Select Ethernet in parameter "I / F Sel".	85
	Communication settings do not match.	Check the communication conditions and set the parameters.	86, 87
cannot be performed with the front panel, command, or jog controller.	"TEACH IF"	Change the parameter "TEACH IF" to "JOG / CMD".	73
	The circuit connected to the general-purpose I / O connector does not meet the specifications.		112
	"JOG / CMD" is selected in parameter "TEACH IF"	Change the parameter "TEACH IF" to "I / O".	73



Contents	Possible cause	Workaround	Page
Teaching operation cannot be performed from the optional jog controller.	Jog controller is disconnected.	Check the connection of the jog controller.  * See the jog controller instruction manual for details.	1, 5
<ul> <li>Emergency stop is not possible.</li> </ul>	You have not changed any parameters.	Change the parameter "EMG Connector" to Enable.	84
cannot be released.	The connector connected to the EMG connector is disconnected		116
<ul> <li>Make an unintended emergency stop.</li> </ul>	The cable connected to the EMG connector is broken	Check the connected cable.	116
	The EMERGENCY switch of the optional jog controller is pressed	Set the EMERGENCY switch of the jog controller to Off.  * See the jog controller instruction manual for details.	1
· Sleep cannot be performed.	The stage is operating.	Wait for positioning to be completed.	12, 43
	An error has occurred.	Please clear the error.	14, 29, 101
	Teaching is running.	Wait for the execution to be completed or stop the execution.	90
Stage does not work.	Motor not energized.	Check the lighting status of the front panel Motor lamp.	4
	selection parameter does not match	·	82
The stage moves on its own.	The axes of the motor cable and scale cable are not aligned.	Check the connection.	<b>®</b> 9
completed.	fixes the stage is vibrating	Make sure that vibration is not transmitted to the stage	-
Coordinate values are not stable during positioning	Strong light is on the stage	Avoid strong light on the stage	-
operation.	Receives strong magnetic fields, electric fields, and noise	Please do not receive them.	-
CMD ER is displayed.	A command string that is not defined in this equipment or that does not match the status of this equipment is sent from the PC.		103
SCALE ER is displayed.	No scale signal is input.	See "11, Status".	104
	Operating at a speed higher than the count capability of the coordinate value counter.		104
OF ER is displayed.	The difference between the coordinate value and the specified position is 5mm or more.		104
LIMIT ER is displayed.	CW and CCW limits are entered at the same time	See "11. Status" .	105
IP ER is displayed.	The magnitude of the scale signal is out of the specified range		105
EMERGENCY is displayed.		See "11. Status".	105
SYS ER is displayed.	The instrument system is broken	See "11. Status" .	105
TCMD ER is displayed.	The content of the teaching registration line to be executed does not match the status of this equipment.		106



1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16

# 15. Update history

Edition	Document control number	Revision date	Supported FV	Supported EN	Contents	
1	MF-1024-01.01	ı	01.030~	01.06	-	
1	MF-1024-01.02	Aug, 18, 2020	01.030~	01.06	<ul><li>Warranty change</li><li>Change of power cable specifications</li></ul>	
1	MF-1024-01.03	Aug, 20, 2020	01.030~	01.06	Vendor name change	
1	MF-1024-01.04	Sep, 2, 2020	01.030~	01.06	Typo correction	
1	MF-1024-01.05	Sep, 9, 2020	01.030~	01.06	lmage change	
2	MF-1024-02.01	Aug, 31, 2021	01.032~	01.06	Change contact URL	

Memo



# 16.Index

Α			_		
В	About commands AC100V AC adapter AC inlet Alarm ASCII AXIS parameter	21 6, 105 9, 114 6, 7 9, 110, 111 21, 102 72	F	F1 / ⇒ button F2 / ⇒ button F3 / ⇒ button F4 / ⇒ button F5 / ⇒ button F6 / ⇒ button F7 / ⇒ button	2 2, 2 2 2 2 2 2 2 2
C	BEC Beep BUSY Busy error cancel	14, 21, 88, 91 11, 65, 87, 116 12, 23, 73 14, 29, 103		F8 / button Feedback stage Feedback stage controller Firmware version Frame connection line Frame terminal	i, 1 i, 1 17, 21, 30, 31 1, 8 iv, 8
	Cancel button CCW CCW LMT CCW S/D CCW SLMT CCW Slowdown sensor input CCW software limit stop CCW limit sensor stop Cleaning Close Close1, 2 lamp Closed loop CLS CM-52 CMD ER Command error Contact information CW CW LMT CW S/D CW SLMT CW slowdown sensor input CW software limit stop	2, 66 13, 14 100, 101, 102 100, 101, 104 100, 101, 104 99, 100 99, 100 V 4, 14, 21, 83, 89 4 4, 14, 28, 83 14 1 101, 102 101, 102 i 13, 14 100, 101, 102 100, 101, 104 100, 101, 104 99, 100 99, 100	G	Front panel Fuse Fuse holder F/V	2 iii, 1, 7, 105, 116 7 17 105 110
E	Del Default gateway Display unit DGW  Earth ground Echo back Electric fan Electric origin setting Electric home return Elongated holes EMERGENCY EMG Connector Enter button Emergency stop EOG EOI Equipment Number Error Ethernet IP address Ethernet interface Ethernet cable	99, 100  18, 19, 64, 84 20, 64, 86, 109 2 20  iii, 7 30, 64 10 14, 22, 45 14, 44, 96 iv, 10 101, 103 9, 117 2, 66 101, 103 14 18, 85, 106, 107 17 101 20, 64, 85 20 9, 20, 84, 85, 86 9	J	Interpolator error  In-position range  Installation IP IP ER IP addoress  JC Connector JOF JON Jog controller  LAN LAN cable LIMIT LIMIT ER Limit error	38, 42, 63, 64, 72, 73, 83, 99, 105, 110, 111 35, 36, 37, 38, 41, 101, 103, 105, 111 12, 35, 36, 37, 38, 42, 63, 72, 73, 99, 105 ii, iv, 3, 8 20 101, 103, 117 20, 64, 85, 109  5 13 13 1, 5  9, 109 109 101, 103, 117 101, 103, 117 101, 103, 117
	Ethernet subnet mask Ethernet Connector Ethernet default gateway Ethernet delimiter Ethernet port number Ethernet lamp Exterior Dimensions E/N	20 9 20 20 20 4 115 17		LMT LOCAL mode LOCAL lamp LSTN lamp	16 13 4 4

1 2 3 4 5	6 7 8	9 10 11 12 13	14   15   16
MAC MAC address MAIN POWER switch Mechanical origin return (Mode0) Mechanical origin return (Mode1) Mechanical origin return (Mode2) Mechanical origin return (Mode3) Mechanical origin return (Mode4) MEU MMT MOD Mode Model Motor1, 2 lamp Model information Model name Motor cable connector Motor driver Movement Test MTR  O  OF ER Off lamp ONE On lamp Open Open loop Option ORG ORG sensor OS ER Over speed error Overflow error P  Page button Parameters PARAMETER mode PAU PLS PORT POWER button Power cable PRM Protective ground wire Precautions for use PST R  READY REMOTE mode Resturn button Return button RST RTN	20 20 3, 6, 7, 8, 116 95 95 96 96 13, 14, 66, 67 16, 98 13 12 17 4 17 17, 21, 30 8, 103, 111, 117 8 ii, 10, 114 16, 90, 97 14 101, 103, 117 3, 116 90 3, 5 4, 14, 21, 83, 89 4, 14, 28, 83 1, 116, 117 14 96 101, 102, 117 101, 102 101, 103 2 63 16 90 13, 14, 66, 67 20 3, 8, 11, 116 1 13 1, 7 ii 13, 14	SMT SRQ lamp SYS ER System error Stage control axes Status STP Subnet mask S/N  T  TALK lamp TCMD ER Teaching Teaching command error TEACH TEACH mode Transmission format Transmission / Reply format Transmission / Reply format Trouble shooting  U  USB cable USB connector USB delimiter USB interface USB lamp  V  Vendor Vendor render W  Wait time Warranty  Z  ZEO Zero control Zero set	14   15   16  16, 97 4 101, 104, 117 101, 104 105 99 14, 16, 90, 92 20, 64, 86, 109 17  4 101, 104, 117 88 101, 104 88 15 23 24 116  9, 108 9 19, 84 9, 18, 19, 84 4  17 17, 21, 30  22, 61, 89 i, ii  14 63, 73 14, 45
Repair Restart Rear panel S	i 3, 26 6		
Safety SCALE ER Scale error Scale cable Scale cable Scale cable connector Serial number Service request Side panel Sleep SLEEP SLEEP SLEEP SNM	ii 101, 102, 117 101, 102 8, 9, 102, 117 9 17, 21, 30, 31 4, 106 10, 115 5 5, 21, 26, 27 5 13 20		