



Highly accurate motorized stage applicable to slideglasses and schales.

- Applicable upright microscopes: Olympus BX series Nikon 80i·50i
- Convert your existing manual stages into precisely-positioning motorized stages.
- Thin Type and Thin Long Travel Type are available, considered space for illumination condenser lens.



▲Thin standard type



▲Thin long travel type

Thin Standard type (Cor	npatible microscopes: (Nikon 80i, 50i)			
Part Number	Resolution [µm]	Travel [mm]	Stage Controller	Joystick Controller	Feedback function
BIOS-205S	0.1	30 (±15)	FC-101G	BJ-02G	0
BIOS-105S	0.1		SC-101G		×

Accessories: Adapter for microscope, exclusive spring clip

Thin long travel type (Compatible microscopes: Olympus BX series Nikon 80i, 50i)							
Part Number	Resolution [µm]	Travel [mm]	Stage Controller	Joystick Controller	Feedback function		
BIOS-235S	0.1	X-axis 60 Y-axis 30	FC-101G	BJ-02G	0		
BIOS-135S	0.1		SC-101G		×		

Accessories: Adapter for microscope, exclusive spring clip

Microscope Unit

Index

Application Systems

Optics & Optical Coatings

Opto-Mechanics

Bases

Manual Stages

Actuators & Adjusters

Motoeized Stages

Light Sources & Laser Safety

Alignment

Interferometers

Inspection/ Observation

Bio-photonics

Laser Processing

■Option code

BIOS-

 $oxedsymbol{oxed}^{\mathsf{OL}}$ NIB

OL: Olympus BX series NIB: Nikon 80i, 50i

Guide

▶ Please inform us of the maker and part number (model) of your existing microscope before purchase.

Sample Software for Simple Stage Control (BIOSControl)





Device Setting



Prese

Guide

Sample programs are available for download from our website.

BIOSCONTROL (2000) FO 1010 (2000) FO 1010 (2000)

(SHOT-702/302GS, FC-101G/501G, SC-101G) Micro-Manager (FC-101G/SC-101G/SHOT-702)